


IPS® **PressVEST** Speed

Instructions for Use
Verarbeitungsanleitung
Mode d'emploi
Istruzioni d'uso
Instrucciones de uso
Instruções de Uso
Bruksanvisning
Brugsanvisning
Käyttöohjeet
Bruksanvisning
Productinformatie
Οδηγίες Χρήσεως

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


1

Direct the sprues towards the wax pattern (imagined continuation of wax pattern).

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


2

The total length of the sprue and pressed object should be max. 15 to 16 mm. Observe an angle of 45 to 60°.

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3


Place sprues in the direction of the flow of the ceramic material.

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
4

If the crown is viewed from the proximal, the longer side of the object (usually the buccal surface) points outwards. Additionally, the flow of the ceramic material must be observed.

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
5

The attachment points of the sprues must be rounded out. Observe an angle of 45 to 60°.

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
6

When divesting the object, blast from the direction indicated above.

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The distance between the objects to be pressed and the silicone ring/ring gauge must be at least 10 mm. This applies to all above schematics.

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Instructions for Use

Indication	<p>The IPS PressVEST Speed Investment material can be used for the following Ivoclar Vivadent pressable ceramics to press in the furnaces EP 500, EP 600, EP 600 Combi and Programat EP 5000:</p> <ul style="list-style-type: none"> – IPS e.max Press – IPS e.max ZirPress – IPS Empress 2 – IPS Empress Cosmo – IPS InLine PoM (Press-on-Metal Ceramic) 	
Contraindication	<ul style="list-style-type: none"> – Metal casting – IPS Empress layering technique (1180 °C/2156 °F) – All uses not explicitly listed as indications by the manufacturer 	
Delivery Form	<ul style="list-style-type: none"> – IPS PressVEST Speed 2,5 kg 25x bags, 100 g each IPS PressVEST Speed powder 0.5 l IPS PressVEST Speed Liquid 1 measuring cup – IPS PressVEST Speed 5 kg 50x bags, 100 g each, IPS PressVEST Speed powder 1 l IPS PressVEST Speed Liquid 1 measuring cup 	
Storage	<ul style="list-style-type: none"> – Ideal storage temperature is 18–23 °C / 64–73 °F – Store powder in a dry place – Do not store the liquid below +5 °C / 41 °C (sensitive to freezing) – Do not use liquid if it has been frozen. <p>If the liquid freezes accidentally, it can no longer be used (crystallization).</p>	
Working temperature	<ul style="list-style-type: none"> – Room temperature 18 °C to max. 23 °C / 64 °F to max. 73 °F) – Working temperature other than the one stipulated, considerably influences the setting behaviour. – A higher working temperature accelerates the chemical process. – A lower working temperature slows down the chemical process. 	
Investment ring system and Alox punger	<ul style="list-style-type: none"> – Please use only the IPS e.max investment system (100 g and 200 g) for the IPS e.max and IPS InLine PoM ingots. 	<ul style="list-style-type: none"> – Please use only the IPS Empress investment system (100 g and 200 g) for the IPS Empress ingots.
		
	<ul style="list-style-type: none"> – IPS e.max and IPS InLine PoM ingots must be used in conjunction with the IPS e.max Alox punger. The diameter of the IPS e.max Alox punger is 13 mm and therefore not compatible for the IPS Empress investment system. – Use the IPS Silicone Ring for investing. The IPS Silicone Ring is used for the IPS e.max as well as the IPS Empress investment system. 	
Powder/Liquid ratio	<ul style="list-style-type: none"> – 100 g powder : 27 ml diluted liquid – 200 g powder : 54 ml diluted liquid 	

Mixing ratio:

IPS e.max Press			
Indication	Concentration	100 g ring Liquid : distilled H ₂ O	200 g ring Liquid : distilled H ₂ O
Single crowns, veneers, partial crowns	60%	16 ml : 11 ml	32 ml : 22 ml
3-unit bridges	50%	16 ml : 11 ml	27 ml : 27 ml
Press on Galvano	60%	16 ml : 11 ml	32 ml : 22 ml

IPS e.max ZirPress			
Indication	Concentration	100 g ring Liquid : distilled H ₂ O	200 g ring Liquid : distilled H ₂ O
Press on ZrO ₂ single crowns	60%	16 ml : 11 ml	32 ml : 22 ml
Press on ZrO ₂ bridges	50%	16 ml : 11 ml	27 ml : 27 ml
Press on ZrO ₂ inlay-bridges	50%	16 ml : 11 ml	27 ml : 27 ml

IPS Empress 2			
Indication	Concentration	100 g Muffel Liquid : distilled H ₂ O	200 g Muffel Liquid : distilled H ₂ O
Single crowns	60%	16 ml : 11 ml	32 ml : 22 ml
3-unit bridges	50%	16 ml : 11 ml	27 ml : 27 ml

IPS Empress Cosmo			
Indication	Konzentration	100 g ring Liquid : distilled H ₂ O	200 g ring Liquid : distilled H ₂ O
Preprosthetic build-up in combination with CosmoPost	40–50%	11 ml : 16 ml 13.5 ml : 13.5 ml	22 ml : 32 ml 27 ml : 27 ml

IPS InLine PoM – Press-on-Metal ceramic			
Indication	Concentration	100 g ring Liquid : distilled H ₂ O	200 g ring Liquid : distilled H ₂ O
Pressing restorations on opaquerized metal frameworks	60%	16 ml : 11 ml	32 ml : 22 ml

The percentages are approximate and can be adjusted if necessary. Do not dilute concentrate to less than 40 %.

Vacuum mixing time

- 2.5 minutes (approx. 350 rpm at room temperature)
- If you use a 'High-Speed' mixing device, shorten the vacuum mixing time accordingly.

Investing

- Slightly vibrate the ring on a dental vibrator while pouring the material into the ring. Avoid vigorous shaking.

Working time

- Approx. 6–7 minutes (mixed at room temperature and 350 rpm)

Setting time

- 30 minutes (min. 30 to max. 45 minutes)
- Do not disturb the ring during the setting stage

- Preheating temperature**
- 850°C/1562 °F
 - Always place the set investment ring in the burnout furnace at 850 °C / 1562 °F.

IPS e.max Press / IPS e.max ZirPress / IPS InLine PoM		
	Ingot	Alox plunger
Preheat	do not preheat	do not preheat

IPS Empress 2 / IPS Empress Cosmo		
	Ingot	Alox plunger
Preheat	do not preheat	preheat

- Holding time, after the preheating temperature of 850 °C/1562 °F has been achieved again**
- 100 g ring – min. 45 minutes.
 - 200 g ring – min. 60 minutes.
 - If several rings (more than 2 rings) are preheated at the same time, the preheating time has to be increased to 15 min./ring

Important notes

- Preparation of the restorations, sprueing, investing, pressing, divesting and cleaning of the pressed objects are carried out according to the corresponding Instructions for Use.
- Completely remove excess separating liquid before investing, as these materials prevent the investment material from setting properly.
- Do not use surfactants (debubblizers), as reactions with the investment material may occur.
- Consistent results can only be achieved by following the Instructions for Use.
- The stipulated working time must be observed. The liquid can be stored in the refrigerator at 5–10 °C / 41–50 °F.
- Dilute the investment mixing liquid only with distilled or deionised water.
- Liquid concentration: Increasing the percentage of distilled water that is added to the investment mixing liquid will lower the amount of setting expansion and compromise the stability of the investment material. Do not use a liquid concentration that is less than 40 % investment mixing liquid.
- Always use clean and dry instruments with the investment material. Do not use mixing beakers that have been used with gypsum.
- The minimum time specifications (setting time/preheating time at 850 °C (1562 °F)) must be observed.
- **Working time**
The working time depends on the material's temperature, the quantity of material to be mixed, the mixing time and the mixing speed of your mixing device.
 - Higher material temperatures shorten the working time.
 - Longer mixing times shorten the working time.
 - The larger the quantities to be mixed, the shorter the working times.
- **Certain factors that influence the setting expansion of the investment material:**
 - Temperature of the investment material and liquid
 - Residual water in the mixing cup
 - Atmospheric humidity
- Generally, wax sprues of max. 3 to 8 mm in length should be used for sprueing of the waxed-up objects. To sprue oblong objects (e.g. delicate MODs), choose a short sprue rather than a long one so that the pressing time is not unnecessarily prolonged.
- Do not mix more than 400 g investment material at the same time.
- Strictly observe mixing ratios.
- Check the function of the vacuum mixing device on a regular basis.
- Do not place the invested ring in a pressure vessel during the setting phase. This procedure changes the physical properties of the investment material.
- Do not leave the investment over night or over the weekend (risk of cracks and flash).
- Always remove pressed rings from the press furnace immediately after the end of the program. Place the ring on a grid and let it cool.
- To achieve sufficient heating of the rings in the preheating furnace, the following instructions have to be observed:
 - Always place the set invested ring in the furnace preheated to 850 °C (1562 °F).
 - Use max. 50% of the available space in the preheating furnace.
 - Always place the rings in the rear part of the firing chamber.
 - Quickly place the invested ring in the burnout furnace. Make sure that the furnace temperature does not drop substantially.

- Always place set invested rings in the burnout furnace with the opening facing downward, preferably at a 45° angle.
- If several investments are carried out at once, the investment rings should be placed into the preheating furnace in 20-minutes' intervals.
- The investment rings must not touch each other. Permit air circulation.
- When divesting the object, make sure that you remove only excess investment material during preblasting at 4 bar pressure. Be careful not to blast the ceramic object at this pressure.

Safety Warnings

- Do not open the burnout furnace during the first 30 minutes of preheating since the wax vapour could ignite when exposed to air and cause a burn!
- Warning: The investment material contains quartz powder. Inhalation of dust should be avoided. (Prolonged exposure to dust may cause lung damage and silicosis).

Physical data according to**ISO 9694 – 1998**

- | | |
|----------------------------|---|
| – Compressive strength | 12–18 N/mm ² – liquid concentration 60 % |
| – Setting time | 9–14 minutes – liquid concentration 70 % |
| – Linear thermal expansion | 0.8–1.3 % – liquid concentration 70 % |

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damage resulting from failure to observe the Instructions for the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent AB
Ivoclar Vivadent UK Limited
Ivoclar Vivadent, Inc.