# bluephase® 20i

Licence to cure

new

LED for every use


CE 0123

For dental use only!  
Caution: U.S. Federal law restricts this device to  
sale by or on the order of a licensed dentist  
Made in Austria  
Ivoclar Vivadent AG, FL-9494 Schaan /  
Liechtenstein

ivoclar  
vivadent®  
clinical

# bluephase® 20i

LED for every use

## Technical data


Wavelength range:	385 – 515 nm
Light intensity:	2,000 – 2,200 mW/cm <sup>2</sup>
Curing modes:	Turbo, High Power, Soft Start, Low Power
Exposure time:	for composites: 10 sec. Tetric EvoCeram / IPS Empress Direct: 5 sec.
Light probe:	10>8 mm black
Power supply:	Lithium polymer battery, approx. 45 min. capacity (curing time); approx. 2 h charging time
Operating voltage:	100-240 VAC / 50-60 Hz
Warranty	3 years (battery 1 year)


## Overview

<b>E</b>	Operating Instructions	Page 6
<b>D</b>	Bedienungsanleitung	Seite 20
<b>F</b>	Mode d'emploi	Page 34
<b>I</b>	Istruzioni d'uso	Pagina 48
<b>Es</b>	Instrucciones de uso	Pagina 62
<b>P</b>	Instruções Operacionais	Página 76
	Appendix	Page 90

English	Deutsch	Français	Italiano	Español	Português
<b>List of parts</b>	<b>Teile- verzeichnis</b>	<b>Liste des pièces détachées</b>	<b>Elenco componenti</b>	<b>Despiece</b>	<b>Lista de partes</b>
① Handpiece	Handstück	Pièce à main	Manipolo	Pieza de mano	Peça de mão
② Start button	Starttaster	Bouton de mise en marche	Tasto di accensione	Interruptor de puesta en marcha	Botão de ligar
③ Display	Display	Affichage	Display	Pantalla	Display
④ Program selection button	Programmwahltaster	Bouton de sélection des programmes	Tasto di selezione programma	Tecla de selección de programa	Tecla de seleção do programa
⑤ Time selection button	Zeitwahltaster	Bouton de sélection du temps	Tasto di selezione tempo	Tecla de selección de tiempo	Tecla de seleção do tempo
⑥ Volume button	Lautstärketaster	Bouton de réglage du volume	Tasto volume	Tecla de volumen	Tecla de volume
⑦ Power indicator	Power-Anzeige	Indicateur de puissance	Indicatore alimentazione	Indicador de potencia	Indicador de potência
⑧ Charging base	Ladestation	Base de chargement	Base di ricarica	Base de carga	Base de carga
⑨ Power cord	Netzkabel	Cordon d'alimentation	Cavo di alimentazione	Cable eléctrico	Cabo de força
⑩ Power pack	Netzgerät	Unité d'alimentation	Alimentatore	Unidad de alimentación	Unidade de alimentação
⑪ Battery	Akku	Batterie	Batteria	Batería	Bateria
⑫ Light probe 10>8 mm	Lichtleiter 10>8 mm	Conducteur de lumière 10>8 mm	Conduttore ottico 10>8 mm	Conducto de luz 10>8 mm	Condutor de luz 10>8 mm
⑬ Anti-glare shield	Blendschutzschild	Ecran de protection	Schermo antiriflesso	Placa antideslumbrante	Escudo anti-ofuscante
⑭ Anti-glare cones	Blendschutzkegel	Embout de protection oculaire	Coni antiriflesso	Conos antideslumbrantes	Cone anti-ofuscante

English	Deutsch	Français	Italiano	Español	Português
<b>Indicators on the charging base</b>	<b>Anzeigen der Ladestation</b>	<b>Indicateurs sur la base de chargement</b>	<b>Indicatori sulla base di ricarica</b>	<b>Indicadores en la base de carga</b>	<b>Indicações da base de carga</b>
 Indicator lights up in blue – the charging base is connected to the power supply	Anzeige leuchtet blau – Ladestation ist an die Stromversorgung ange- schlossen	Les voyants s'allument en bleu – la base de chargement est reliée au courant électrique	Indicatore blu illuminato – la base di ricarica è collegata alla rete elettrica	Indicador se ilumina en azul – la base de carga está conectada a la red eléctrica	Indicador mostra uma luz azul – A base de carga está ligada com a rede elétrica
<b>Indicators on the handpiece</b>	<b>Anzeigen des Handstücks</b>	<b>Indicateurs sur la pièce à main</b>	<b>Indicatori sul manipolo</b>	<b>Indicadores en la pieza de mano</b>	<b>Indicações da peça de mão</b>
 Curing program and op- erating status	Belichtungsprogramm bzw. Betriebszustand	Programme de polymé- risation, temps de poly- mérisation, état de marche	Programma di polimerizzazione e stato operativo	Programa de polimerización y estado de funcionamiento	Programa de cura e es- tado de operação
 Selected indications and curing time	Angewählte Indikation bzw. Belichtungsdauer	Indications spécifiques du programme de poly- mérisation sélectionné	Indicazione selezionata e tempo di polimerizza- zione	Indicaciones seleccionadas y tiempo de polimerización	Indicação selecionada e tempo de cura
 Charging status of the battery	Ladezustand Akku	Témoïn de charge de la batterie	Stato di carica della batteria	Estado de carga de la batería	Status de carga da bате- ria
 Display is in stand-by mode	Display ist im Stand-by Modus	L'écran est en mode stand-by	Display in modalità Stand-by	Display in modo stand-by	Display está no modo "stand-by"


Signs and symbols	Zeichen- erklärung	Signes et symboles	Segni e simboli	Señales y símbolos	Sinais e símbolos
 Non-permissible use	Nicht zulässige Anwendung	Ne pas utiliser	Uso non previsto	Uso no permitido	Uso não autorizado
 Double insulation (Apparatus complies with safety class II)	Doppelt isoliert (Gerät der Schutzklasse II)	Double protection (l'appareil répond aux normes de sécurité de classe II)	Doppio isolamento (ap- parecchio conforme alla Classe di sicurezza II)	Doble aislamiento (El dispositivo cumple con la Clase de Seguri- dad II)	Isolamento duplo (o aparelho cumpre a norma de segurança da classe II)
 Protection against electrical shock (BF type apparatus)	Schutz gegen elektrischen Schlag (Gerätetyp BF)	Protection contre les décharges électriques (type d'appareil BF)	Protezione da scossa elettrica (apparecchio di tipo BF)	Protección frente al cho- que eléctrico (aparato tipo BF)	Proteção contra choque elétrico (aparelho tipo BF)
 Observe Operating Instructions	Bedienungsanleitung beachten	Respecter le mode d'emploi	Attenersi alle istruzioni d'uso	Observar las instruccio- nes de uso	Observar as Instruções de Uso
 The apparatus must not be disposed of in the normal domestic waste	Das Gerät darf nicht über den normalen Hausmüll entsorgt werden	Ne pas jeter l'appareil avec les ordures ménagères	L'apparecchio non può essere smaltito come normale rifiuto urbano	El aparato no se puede eliminar con la basura doméstica normal	O aparelho não pode ser descartado junto com o lixo doméstico normal

Light intensity	2,000 – 2,200 mW/cm <sup>2</sup>
Exposure time for composites	10 s
Exposure time for Tetric EvoCeram / IPS Empress Direct	5 s

### TURBO Program

Consistently high light intensity for the polymerization of restorative and cementation materials for direct and indirect restorations.

Filling materials	Curing time
<b>Composite</b>	
• 2 mm <sup>1)</sup>	
IPS Empress Direct / IPS Empress Direct Flow / Tetric EvoCeram / Tetric EvoFlow / Tetric / Tetric Basic White / Tetric N-Ceram / Tetric N-Flow	5 s
Heliomolar / Heliomolar HB / Heliomolar Flow / Others	2x 5 s
• 4 mm <sup>2)</sup>	
Tetric EvoCeram Bulk Fill	5 s
Tetric Basic White	2x 5 s
<b>Compomer</b> <sup>3)</sup>	
Compoglass F / Compoglass Flow	2x 5 s
<b>Indirect restorations / Luting materials</b>	
Variolink II Base <sup>4)</sup> / Variolink Veneer / Variolink N Base <sup>4)</sup> / Variolink N Clear Veneer / Variolink Ultra <sup>5)</sup> / Dual Cement <sup>3)</sup> / Variolink II <sup>3)</sup>	per mm ceramic: 1x 5 s per segment
<b>Miscellaneous</b>	
Heliosit Orthodontic	2x 5 s
Telio Add-On Flow	2x 5 s
Telio Stains	5 s
IPS Empress Direct Color	5 s
IPS Empress Direct Opaque	2x 5 s

### HIGH POWER Program

Consistently high light intensity for the polymerization of restorative and cementation materials for direct and indirect restorations.


Filling materials	Curing time
<b>Composite</b>	
• 2 mm <sup>1)</sup>	
IPS Empress Direct / IPS Empress Direct Flow / Tetric EvoCeram / Tetric EvoFlow / Tetric / Tetric Basic White / Tetric N-Ceram / Tetric N-Flow	10 s
Heliomolar / Heliomolar HB / Heliomolar Flow / Others	15 s
• 4 mm <sup>2)</sup>	
Tetric EvoCeram Bulk Fill / Tetric Basic White	10 s
<b>Compomer</b> <sup>3)</sup>	
Compoglass F / Compoglass Flow	20 s
<b>Indirect restorations / Luting materials</b>	
Variolink II Base <sup>4)</sup> / Variolink Veneer / Variolink N Base <sup>4)</sup> / Variolink N Clear Veneer / Variolink Ultra <sup>5)</sup> / Dual Cement <sup>3)</sup> / Variolink II <sup>3)</sup>	per mm ceramic: 10 s per segment
<b>Miscellaneous</b>	
Helioseal / Helioseal F / Helioseal Clear	10 s
Monopaque	20 s
MultiCore Flow / Multicore HB	20 s
Heliosit Orthodontic	10 s
Telio Add-On Flow	15 s
Telio Stains	10 s
IPS Empress Direct Color	10 s
IPS Empress Direct Opaque	20 s

### LOW POWER Program

Reduced light intensity with reduced heat development for the polymerization of adhesives, liners, and restorative materials in areas near the pulp when restoring Class V cavities.


Adhesive	Curing time
AdheSE / AdheSE One F	
ExcITE / ExcITE DSC	10 s
Tetric N-Bond / Tetric N-Bond Self-Etch	
Heliobond (Syntac)	
<b>Temporary materials</b>	
Telio CS Inlay/Onlay	10 s
System.inlay/onlay	10 s
Fermit / Fermit N	10 s
Telio CS Link / System.link	20 s per segment
<b>Miscellaneous</b>	
Heliosit Orthodontic	20 s
Vivaglass Liner	20 s

### SOFT START Program

Step-by-step increase of the light intensity with reduced shrinkage stress and reduced heat development for the polymerization of restorative materials.

Filling materials	Curing time
<b>Composite</b>	
• 2 mm <sup>1)</sup>	
IPS Empress Direct / IPS Empress Direct Flow / Tetric EvoCeram / Tetric EvoFlow / Tetric / Tetric Basic White / Tetric N-Ceram / Tetric N-Flow	15 s
Heliomolar / Heliomolar HB / Heliomolar Flow / Others	20 s
• 4 mm <sup>2)</sup>	
Tetric EvoCeram Bulk Fill / Tetric Basic White	15 s
<b>Compomer</b> <sup>3)</sup>	
Compoglass F / Compoglass Flow	20 s

1) Applies to a maximum layer thickness of 2 mm and provided that the Instructions for Use of the respective material do not state any other recommendation (might be the case, e.g. with dentin shades)  
 2) Applies to a maximum layer thickness of 4 mm and provided that the Instructions for Use of the respective material do not state any other recommendation (might be the case, e.g. with dentin shades)  
 3) Applies to a maximum layer thickness of 3 mm 4) Applies to light-curing (use of the base past only) 5) Applies to dual-curing

## Dear Customer

Optimum polymerization is an important requirement for all light cured materials in order to consistently produce high quality restorations. The polymerization light selected also plays a decisive role in this respect. Therefore, we would like to thank you for having purchased bluephase® 20i. bluephase 20i is a high-quality medical device which has been designed according to the latest standard of science and technology in compliance with the relevant industry standards.

These Operating Instructions will help you to safely start up the device, make full use of its capabilities and ensure a long service life.

Should you have any further questions, please do not hesitate to contact us (see addresses on the back cover).

Your

Ivoclar Vivadent Team

# Table of Contents

<b>Safety</b>	
Intended use	8
Indications	8
Safety notes	8
<b>Start-up</b>	
Charging base	10
Handpiece	11
Battery	11
Click & Cure corded operation	12
<b>Operation</b>	
Selecting the curing program and the curing time	13
Cure Memory function	16
<b>Start</b>	
Acoustic signals	16
Light intensity	16
Measuring the light intensity	16
<b>Maintenance and cleaning</b>	
Housing	17
Light probes	17
Battery contacts	17
Disposal	17
<b>What if ...</b>	<b>18</b>
<b>Warranty / Procedure in case of repair</b>	<b>19</b>
<b>Product specifications</b>	
Delivery form	19
Accessories	19
Technical data	19
Operating conditions	19
Transportation and storage conditions	19

# Safety

## Intended use

bluephase 20i is an LED polymerization light that produces energy-rich blue light. It is used for the polymerization of light-curing dental materials immediately at the dental unit. The intended use also includes the observation of the notes and regulations in these Operating Instructions.

## Indication

With its "polywave" broadband spectrum, bluephase is suitable for the polymerization of all light-curing dental materials curing in the wavelength range of 385–515 nm. These materials include restoratives, bonding agents/adhesives, bases, liners, fissure sealants, temporary materials as well as luting materials for brackets and dental-lab materials such as those used in the fabrication of ceramic inlays.

## Contraindications

Materials, the polymerization of which is activated outside the wavelength range of 385–515 nm (no materials known to date). If you are not sure about certain products, please ask the manufacturer of the corresponding material.

Do not charge or use the appliance near flammable or combustible substances.

Portable and mobile high-frequency communication devices may interfere with medical equipment. The use of mobile phones during operation is not allowed.


Caution – The use of controls or adjustment devices or the performance of procedures other than those specified herein may result in hazardous radiation exposure.

## Safety notes

The bluephase 20i is a medical electronic device which is subject to IEC 601-1 (EN 60601-1) and EMC regulations. It complies with the relevant EU regulations and is classified as LED Class 2.

The apparatus has been shipped from the manufacturer in a safe and technically sound condition. In order to maintain this condition and to ensure risk-free operation, the notes and regulations in these Operating Instructions have to be observed. To prevent damage to equipment and risks for patients, users and third parties, the following safety instructions have to be observed.

## Usage and liability

- bluephase 20i must only be employed for the intended use. Any other uses are contraindicated. Liability cannot be accepted for damage resulting from misuse or failure to observe the Operating Instructions.
- The user is responsible for testing bluephase 20i for its use and suitability for the intended purposes. This is particularly important if other equipment is used in the immediate vicinity of bluephase at the same time.
- Only original spare parts and accessories from Ivoclar Vivadent (see Accessories) should be used. The manufacturer does not accept any liability for damage resulting from the use of other spare parts or accessories.

## Operating voltage

Before switching on, make sure that

- a) the voltage indicated on the rating plate complies with the local power supply and
- b) the unit has acquired the ambient temperature.

If the battery or power pack are used separately, e.g. during start-up or Click & Cure corded operation, contact with patients or third parties must be prevented. Do not touch the exposed contacts of the battery or connection plug (power pack).

## Assumption of impaired safety

If it is assumed that safe operation is no longer possible, the power must be disconnected and the battery removed to avoid accidental operation. This may be the case, for example, if the apparatus is visibly damaged or no longer works correctly. Complete disconnection from the power supply is only ensured when the power cord is disconnected from the power source.

## Eye protection

Direct or indirect exposure of the eyes must be prevented. Prolonged exposure to the light is unpleasant to the eyes and may result in injury. Therefore, using the supplied anti-glare cones is recommended. Individuals who are generally sensitive to light, who take photosensitizing drugs, have undergone eye surgery, or people who work with the apparatus or in its vicinity for long periods of time should not be exposed to the light of this device. Therefore, they should wear protective goggles (orange) that absorb light below 515 nm.

## Battery

Caution: Use only original spare parts, particularly Ivoclar Vivadent batteries and charging bases. Do not short circuit battery. Do not store at temperatures above 40 °C / 104 °F (or 60 °C / 140 °F for a short period). Always store batteries charged. The storage period must not exceed 6 months. May explode if disposed of in fire.

Please note that the lithium polymer battery may react with explosion, fire, and smoke development in case of inappropriate handling and mechanical damage. Damaged lithium polymer batteries must no longer be used.

The electrolytes and electrolyte fumes released during explosion, fire, and smoke development are toxic and corrosive. In case of contact with the eyes and skin, immediately rinse with copious amounts of water. Avoid the inhalation of the fumes. See a physician immediately in case of indisposition.

## Heat development (burn hazard)

As it is the case with all high-performance lights, the high light intensity results in a certain heat development. Prolonged exposure of areas near the pulp and soft tissues may result in irreversible or reversible damage. Therefore, this high-performance curing light must only be operated by trained professionals.

The recommended curing times have to be observed. Uninterrupted curing times of more than 5 seconds in the Turbo-Program on the same tooth surface, as well as direct contact with the gingiva, oral mucous membrane or skin have to be prevented. If required, polymerize in intermittent intervals of at most 5 seconds each. Prevent heat accumulation when polymerizing with a rubber dam in place. If irradiation of soft tissue cannot be avoided, operate with reduced light intensity (use Low Power program). Polymerize indirect restorations at intermittent intervals of 10 seconds in the High-Power program and 1x 5 seconds each in the Turbo program or use external cooling with an air stream.

The instructions regarding curing programs and curing times must be observed (see Selecting the curing program and the curing time). In addition, the light emission window must be placed directly over the material to be cured at all times (e.g. by holding it in place using a finger).

It is not recommended to use the High Power/Turbo program for sensitive patients. In this case, operate with reduced light intensity (use Low Power program).

## Start-up

### Delivery form

Check the delivery for completeness and any possible transportation damage (see Delivery form). If components are missing or damaged, immediately contact your dealer or your Service Center.

### Charging base

Before you switch on the device, make sure that the voltage mentioned on the rating plate complies with your local power supply. The rating plate is attached to the underside of the charging station.

Place the charging base on a suitable, flat table top. Remove the protective foil from the connection plug of the power pack. Slide the connection plug of the power pack into the socket on the underside of the charging base. Tilt it slightly and apply light pressure until you hear and feel it snap into place.


Connect the power cord with the power supply and the power pack.

The **Power ON** indicator on the left side of the housing lights up in blue (see Indicators on the charging base).

## Handpiece

Before you insert the light probe, remove the protective foil on the corresponding opening of the handpiece.

Slightly turn the light probe while attaching it to the handpiece.

After that, mount the anti-glare shield on the light probe.

## Battery

We recommend fully charging the battery before the first use.

Slide the battery straight into the handpiece until you hear and feel it click into place.

Gently place the handpiece in the corresponding cradle in the charging base. If a sleeve is used, please remove it before you charge the battery. If the battery is fully charged, it features a curing capacity of approximately 45 minutes.

## Tip

If possible, use the light always with a fully charged battery. This will prolong its service life. It is therefore recommended to place the handpiece into the charging base after each appointment. The charging time for empty batteries is 2 hours.

## Start-up

### Click & Cure corded operation

bluephase 20i can be used in corded operation at any time, but particularly when the battery is completely empty.

For this purpose, remove the battery from the handpiece. Then remove the power pack from the underside of the charging base. Do not pull on the power cord.

Insert the connection plug straight into the handpiece until you hear and feel it click into place.

During corded operation, the charging base cannot charge the battery, since it is not connected to a power source.

The apparatus is only completely disconnected from the power supply if the power cord has been pulled out of the plug-in socket.

## Charging status of the battery

With the handpiece switched on, the current charging status is shown on the display as follows:


Battery fully charged  
(curing capacity of approximately 45 minutes).

Battery half full


Reserve (The last bar in the battery indicator is red in the reserve mode. The battery has to be charged as soon as possible.)

When the battery is being charged, "Charging Battery" briefly appears in the display before the display returns to the stand-by mode.

If the battery is completely empty, the handpiece automatically switches to the stand-by mode. The light can no longer be called up and the curing program and curing time can no longer be set. However, the handpiece can be used in the Click & Cure corded operation.

Since the battery is a consumable, it has to be replaced after its typical life cycle has expired after 2½ years. See battery label for the age of the battery.

S000001 04/04/09  
dd/mm/yy

# Operation

Disinfect or autoclave contaminated surfaces of the curing light as well as light probes and anti-glare cones before each use.

Furthermore, make sure that the stipulated light intensity permits adequate polymerization. For this purpose, check the light probe for contamination and damage as well as the light intensity at regular intervals (e.g. with the Ivoclar Vivadent radiometer "bluephase meter").

## Selecting the curing program and the curing time

The curing program and the curing time can be individually set. bluephase 20i is equipped with the following 4 curing programs for the different indications. The desired curing program is selected with the program selection buttons. The display changes accordingly (see Indicators on the hand-piece). The device comes equipped with the following pre-set programs:

The intended curing time is selected using the time selection buttons. Users may choose between 5, 10, 15, 20 and 30 seconds.

Observe the Instructions for Use of the material applied when selecting the curing time.


The curing recommendations for composite materials apply to all shades and, if not mentioned otherwise in the Instructions for Use, to a maximum layer thickness of 2 mm. Generally, these recommendations apply to situations where the emission window of the light probe is placed directly over the material to be polymerized. Increasing the distance between the light source and the material will require the curing time to be extended accordingly. For instance, if the distance to the material is 8 mm, the effective light output is reduced by approx. 50%. In this case, the recommended curing time has to be doubled.

- 1) The information provided herein applies to the light probe 10>8 mm supplied in the delivery form.
- 2) The information regarding heat development and burn hazard must be taken into consideration (see Safety notes).

## Factory settings

	<b>TURBO</b>	<b>5 seconds</b>
	<b>HIGH (High Power)</b>	<b>10 seconds</b>
	<b>LOW (Low Power)</b>	<b>10 seconds</b>
	<b>SOFT (Soft start)</b>	<b>15 seconds</b>

## Curing program



### TURBO Program

Consistently high light intensity for the polymerization of restorative and cementation materials for direct and indirect restorations.

Light intensity	2,000 – 2,200 mW/cm <sup>2</sup>
Exposure time for composites	10 s
Exposure time for Tetric EvoCeram / IPS Empress Direct	5 s

Filling materials	Curing time
<b>Composite</b>	
• 2 mm <sup>1)</sup>	
IPS Empress Direct / IPS Empress Direct Flow / Tetric EvoCeram / Tetric EvoFlow / Tetric / Tetric Basic White / Tetric N-Ceram / Tetric N-Flow	5 s
Heliomolar / Heliomolar HB / Heliomolar Flow / Others	2x 5 s
• 4 mm <sup>2)</sup>	
Tetric EvoCeram Bulk Fill	5 s
Tetric Basic White	2x 5 s
<b>Compomer</b> <sup>3)</sup>	
Compoglass F / Compoglass Flow	2x 5 s
<b>Indirect restorations / Luting materials</b>	
Variolink II Base <sup>4)</sup> / Variolink Veneer / Variolink N Base <sup>4)</sup> / Variolink N Clear Veneer / Variolink Ultra <sup>4)</sup> / Dual Cement <sup>5)</sup> / Variolink II <sup>5)</sup>	per mm ceramic: 1x 5 s per segment
<b>Miscellaneous</b>	
Heliosit Orthodontic	2x 5 s
Telio Add-On Flow	2x 5 s
Telio Stains	5 s
IPS Empress Direct Color	5 s
IPS Empress Direct Opaque	2x 5 s

### HIGH POWER Program

Consistently high light intensity for the polymerization of restorative and cementation materials for direct and indirect restorations.


Filling materials	Curing time
<b>Composite</b>	
• 2 mm <sup>1)</sup>	
IPS Empress Direct / IPS Empress Direct Flow / Tetric EvoCeram / Tetric EvoFlow / Tetric / Tetric Basic White / Tetric N-Ceram / Tetric N-Flow	10 s
Heliomolar / Heliomolar HB / Heliomolar Flow / Others	15 s
• 4 mm <sup>2)</sup>	
Tetric EvoCeram Bulk Fill / Tetric Basic White	10 s
<b>Compomer</b> <sup>3)</sup>	
Compoglass F / Compoglass Flow	20 s
<b>Indirect restorations / Luting materials</b>	
Variolink II Base <sup>4)</sup> / Variolink Veneer / Variolink N Base <sup>4)</sup> / Variolink N Clear Veneer / Variolink Ultra <sup>4)</sup> / Dual Cement <sup>5)</sup> / Variolink II <sup>5)</sup>	per mm ceramic: 10 s per segment
<b>Miscellaneous</b>	
Helioclear / Helioclear F / Helioclear Clear	10 s
Monopaque	20 s
MultiCore Flow / Multicore HB	20 s
Heliosit Orthodontic	10 s
Telio Add-On Flow	15 s
Telio Stains	10 s
IPS Empress Direct Color	10 s
IPS Empress Direct Opaque	20 s

### LOW POWER Program

Reduced light intensity with reduced heat development for the polymerization of adhesives, liners, and restorative materials in areas near the pulp when restoring Class V cavities.

Adhesive	Curing time
AdheSE / AdheSE One F	
ExcTE / ExcTE DSC	10 s
Tetric N-Bond / Tetric N-Bond Self-Etch	
Heliobond (Syntac)	
<b>Temporary materials</b>	
Telio CS Inlay/Onlay	10 s
System.inlay/onlay	10 s
Fermit / Fermit N	10 s
Telio CS Link / System.link	20 s per segment
<b>Miscellaneous</b>	
Heliosit Orthodontic	20 s
Vivaglass Liner	20 s

### SOFT START Program

Step-by-step increase of the light intensity with reduced shrinkage stress and reduced heat development for the polymerization of restorative materials.

Filling materials	Curing time
<b>Composite</b>	
• 2 mm <sup>1)</sup>	
IPS Empress Direct / IPS Empress Direct Flow / Tetric EvoCeram / Tetric EvoFlow / Tetric / Tetric Basic White / Tetric N-Ceram / Tetric N-Flow	15 s
Heliomolar / Heliomolar HB / Heliomolar Flow / Others	20 s
• 4 mm <sup>2)</sup>	
Tetric EvoCeram Bulk Fill / Tetric Basic White	15 s
<b>Compomer</b> <sup>3)</sup>	
Compoglass F / Compoglass Flow	20 s

- 1) Applies to a maximum layer thickness of 2 mm and provided that the Instructions for Use of the respective material do not state any other recommendation (might be the case, e.g. with dentin shades)
- 2) Applies to a maximum layer thickness of 4 mm and provided that the Instructions for Use of the respective material do not state any other recommendation (might be the case, e.g. with dentin shades)
- 3) Applies to a maximum layer thickness of 3 mm
- 4) Applies to light-curing (use of the base past only)
- 5) Applies to dual-curing

# Operation

## Cure Memory function

The last settings used, together with the combination of curing program and curing time, are automatically saved.

## Start

The light is switched on by means of the start button. Once the selected curing time has elapsed, the curing program is automatically terminated. If desired, the light can be switched off before the set curing time has elapsed by pressing the start button again. The fan is activated simultaneously to the light. Once the curing time has elapsed, the fan continues to run for a certain time to cool the apparatus. The battery must not be removed as long as the fan is still running.

## Acoustic signals

Acoustic signals can be heard for the following functions:

- Start (Stop)
- Every 10 seconds
- Program change
- Curing time change
- Connecting the battery to the charging base
- Inserting battery
- Error message

If desired, the volume of the acoustic signals can be adjusted.

For this purpose, press the blue volume button (see Product overview) while the light is off and reduce the volume or turn the acoustic signals off with the left time selection button. If you want to turn the acoustic signals back on or increase the volume, press the blue volume button and then the right time selection button.

## Light intensity

The light intensity is maintained at a consistent level during operation. If the supplied 10>8 mm light probe is used, the light intensity has been calibrated to 2,000 – 2,200 mW/cm<sup>2</sup>.

The use of a light probe other than the one provided has a direct influence on the light intensity emitted.

In parallel-walled light probes (10 mm), the diameter is equal at both ends. In focussing light probes (10>8 mm light probe, Pin-Point light probe 6>2 mm), the diameter of the rear end is larger than that of the light emission window. The incident blue light is thus bundled to a smaller surface, which increases the light intensity emitted.

Pin-Point light probes are suitable for the polymerization of confined areas, such as the attachment of veneers prior to excess removal. For thorough curing, it is necessary to change the light probe.

## Measuring the light intensity

The light intensity of the bluephase 20i and the supplied 10>8 mm light probe can be checked, for instance, by means of the bluephase meter.

If the measured value does not correspond with the expected light intensity, proceed as follows:

- Check the selected curing program
- Clean light sensor if it is contaminated
- Remove the light probe and clean the light emission window of the handpiece with a cotton swab dipped in alcohol
- Clean light probe if it is contaminated (see Maintenance and cleaning)
- Replace a damaged light probe with a new one

If these measures do not improve the results, please contact your dealer or your local Service Center.

# Maintenance and cleaning

For reasons of hygiene, it is recommended to use a disposable protective sleeve for each patient. Contaminated surfaces of the device and anti-glare cones must be disinfected (e.g. FD333, FD366/Dürr Dental, Incidin Liquid/Ecolab) and the light probe must be sterilized before each use, unless disposable protective sleeves are used. Make sure that no liquids or other foreign substances enter the handpiece, charging base and particularly the power pack during cleaning (risk of electrical shock). Disconnect the charging base from the power source when cleaning it.

## Battery contacts

For good conductivity at all times, keep the battery contacts free from possible contamination (e.g. composite residue). For this purpose, clean the affected contacts regularly with disinfectant wipes (after each patient).

## Housing

Wipe the handpiece and handpiece holder with a customary aldehyde-free disinfecting solution. Do not clean with highly aggressive disinfecting solutions (e.g. solutions based on orange oil or with an ethanol content of more than 40%), solvents (e.g. acetone), or pointed instruments, which may damage or scratch the plastic. Clean dirty plastic parts with a soap solution.

## Light probes



Clean the light probe with a disinfectant and a soft cloth. Carefully remove any possible material residue, such as composite material, from the light probe with a fingernail or a plastic spatula, for instance. Do not use sharp or pointed instruments, since they may scratch the surface of the light probe and thus reduce light transmission. Light probes and anti-glare cones may be autoclaved (e.g. 121 °C, 2 bar, 15 minutes). Check light probes for damage. If a light probe is held against a light source and certain segments appear black, then the glass fibres are broken. In this case, the light probe has to be replaced with a new one.

## Disposal

The apparatus must not be disposed of in normal domestic waste. Please visit your local Ivoclar Vivadent website for detailed information about the disposal of the unit. Dispose of unserviceable batteries and polymerization lights according to the corresponding legal requirements in your country. Batteries must not be incinerated.

# What if ...

Symbol	Causes	Error rectification
 <p><b>ERROR</b> 06</p>	Electronic defect in the handpiece or battery	Remove and reinsert the battery. If the error persists, replace the battery with the power pack (Click & Cure). If the error persists, please contact your dealer or your local Service Center.
 <p><b>COLD/ HOT</b></p>	Apparatus is overheated or under-cooled	Allow the apparatus to cool down (or assume room temperature if it is undercooled) and try again after a certain time. If the error persists, please contact your dealer or your local Service Center.
 <p><b>CABLE</b></p>	Electronic defect in the battery	Remove and reinsert the battery. If the error continues to be displayed, place the apparatus into the charging base. If the error persists nonetheless, replace the battery with the power pack (Click & Cure). Please contact your dealer or your local Service Center.
 <p><b>LOW BATTERY</b></p>	Battery empty	Place the handpiece into the charging base. The corresponding message („Charging“) will be displayed. If the battery is not charged, the contacts have to be cleaned. If the error persists, please contact your dealer or your local Service Center.
 <p><b>No display during charging</b></p>	Charging contacts contaminated Battery in safety mode (integrated protective circuit)	Check if the handpiece has been correctly placed into the charging base. Clean battery contacts. Information: "Charging Battery" briefly appears in the display before the display returns to the stand-by mode. If cleaning the battery contacts does not resolve the problem, remove the battery from the apparatus. Charge the battery separately from the handpiece in the charging base for approximately 10 minutes. If the error persists, please contact your dealer or your local Service Center.
 <p><b>LED of the charging base is not illuminated</b></p>	Power pack not connected or defective	Check the battery contacts and if the power pack is correctly positioned in the charging base or if the power pack is connected to the power supply by means of the power cord (display on the power pack lights up in green if it works correctly). If the error persists, please contact your dealer or your local Service Center.

## Warranty / Procedure in case of repair

The warranty period for bluephase 20i is 3 years from the date of purchase (battery 1 year). Malfunctions resulting from faulty material or manufacturing errors are repaired free of charge during the warranty period. The warranty does not provide the right to recover any material or non-material damage other than the ones mentioned. The apparatus must only be used for the intended purposes. Any other uses are contraindicated. The manufacturer does not accept any liability resulting from misuse and warranty claims cannot be accepted in such cases. This is particularly valid for:

- damage resulting from improper handling, especially incorrectly stored batteries (see Technical data: Transportation and storage conditions).
- damage to components resulting from wear under standard operating conditions (e.g. battery).
- damage resulting from external influences, e.g. blows, drop to the floor.
- damage resulting from incorrect set-up or installation.
- damage resulting from connecting the unit to a power supply, the voltage and frequency of which do not comply with the ones stated on the rating plate.
- damage resulting from improper repairs or modifications that have not been carried out by certified Service Centres.

In case of a claim under warranty, the complete apparatus (handpiece, charging base, power cord and power pack) must be returned, carriage paid, to the dealer or directly to Ivoclar Vivadent, together with the purchase document. Use the original packaging with the corresponding cardboard inserts for transportation.

Repair work may only be carried out by a certified Ivoclar Vivadent Service Center. In case of a defect that cannot be rectified, please contact your dealer or your local Service Center (see addresses on the reverse side). A clear description of the defect or the conditions under which the defect occurred will facilitate locating the problem. Please enclose this description when returning the apparatus.

## Product specifications

### Delivery form

1 Charging base, 1 Power cord, 1 Power pack, 1 Handpiece  
1 Battery (Li-Po), 1 10>8 mm light probe, 1 Anti-glare shield,  
3 Anti-glare cones, Sleeves, 1 Set of Operating Instructions

### Accessories

The following accessories are available for bluephase 20i:

REF	Description
627389	Light probe 10>8 mm, black (G2)
608538	Pin-point light probe 6/2 mm, black (G2)
551756	Anti-glare cones
592496	Anti-glare shield
627299	bluephase 20i charging base (G2)
627300	bluephase 20i battery (G2)
627298	bluephase 20i handpiece (G2)
607922	bluephase meter
608554	bluephase sleeves (G2)

### Technical data

Operating voltage charging base	5 VDC
Operating voltage handpiece	3.7 VDC with battery 5 VDC with power pack
Power Pack	100-240 VAC / 50-60 Hz / max. 0.4 A / Output 5 VDC TYPE 15.2630 Manufacturer: Friwo

### Operating conditions

Temperature	+10 °C to +40 °C (+50 to +104°F)
Relative humidity	30 % to 75%
Ambient pressure	700 hPa to 1060 hPa
Dimensions of the charging base	L=205mm; W=150mm; H=85mm
Weight of the charging base	250 g
Charging time	approx. 2 h (with the battery empty)
Power supply of the handpiece	Li-Po battery
Max. battery time	approx. 45 min. (with a new, fully charged battery)
Light source	polywave LED
Wavelength range	385–515 nm
LED-Class 2	LED LIGHT DO NOT STARE INTO BEAM CLASS 2 LED PRODUCT Max. 5.11 mW / Wavelength 385 - 515 nm
Light intensity	2,000 - 2,200 mW/cm <sup>2</sup>
Operation	6 min. on / 4 min. off (intermittently)
Dimensions of the handpiece	L=260mm; W=42mm; H=120mm
Weight of the handpiece	(incl. battery) 245 g

### Transportation and storage conditions:

Temperature	-20 °C to +60 °C (-4 to 140°F)
Relative humidity	10% to 75%
Ambient pressure	500 hPa to 1060 hPa

bluephase 20i has to be stored in closed, roofed rooms.

Protect the device from severe jarring.

### Battery

- Do not store at temperatures above 40 °C / 104 °F (or 60 °C / 140 °F for a short period). Recommended storage temperature 15-30 °C (59-86 °F).
- Always store the battery charged and
- not for longer than 6 months.

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# Appendix

## Guidance and manufacturer's declaration - electromagnetic emission

The bluephase 20i is intended for use in the electromagnetic environment specified below. The customer or the user of the bluephase 20i should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The bluephase 20i uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The bluephase 20i is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	$P < 75W$ (pass without test)
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	$P < 75W$ (pass without test)

## Guidance and manufacturer's declaration - electromagnetic immunity

The bluephase 20i is intended for use in the electromagnetic environment specified below. The customer or the user of the bluephase 20i should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	$\pm 6$ kV contact $\pm 8$ kV air	$\pm 6$ kV contact $\pm 8$ kV air	Floors should be concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	$\pm 2$ kV for power supply lines $\pm 1$ kV for input/output lines	$\pm 2$ kV for power supply lines $\pm 1$ kV for input/output lines	Mains power quality should be that of typical commercial or dental environment.
Surge IEC 61000-4-5	$\pm 1$ kV line - line $\pm 2$ kV line - earth	$\pm 1$ kV line - line no prot. earth	Mains power quality should be that of typical commercial or dental environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	$<5\% U_T$ ( $>95\%$ dip in $U_T$ ) for 0.5 cycle $40\% U_T$ ( $60\%$ dip in $U_T$ ) for 5 cycles $70\% U_T$ ( $30\%$ dip in $U_T$ ) for 25 cycles $<5\% U_T$ ( $>95\%$ dip in $U_T$ ) for 5 sec	$<5\% U_T$ ( $>95\%$ dip in $U_T$ ) for 0.5 cycle $40\% U_T$ ( $60\%$ dip in $U_T$ ) for 5 cycles $70\% U_T$ ( $30\%$ dip in $U_T$ ) for 25 cycles $<5\% U_T$ ( $>95\%$ dip in $U_T$ ) for 5 sec	Mains power quality should be that of typical commercial or dental environment. If the user of the bluephase 20i requires continued operation during power mains interruptions, it is recommended that the bluephase 20i be powered from an uninterruptible power supply or battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or dental environment.

NOTE:  $U_T$  is the a.c. mains voltage prior to application of the test level.

## Guidance and manufacturer's declaration - electromagnetic immunity

The bluephase 20i is intended for use in the electromagnetic environment specified below. The customer or the user of the bluephase 20i should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should not be used closer to any part of the bluephase 20i, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
			<b>Recommended separation distance</b>
Conducted RF IEC 61000-4-6	3 V <sub>rms</sub> 150 kHz to 80 MHz	10 V	$d = 0.35 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	10 V/m	$d = 0.35 \sqrt{P}$ 80 MHz to 800 MHz $d = 0.70 \sqrt{P}$ 800 MHz to 2.5 GHz
			where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).  Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup>  Interference may occur in the vicinity of equipment marked with the following symbol:  

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic

environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the bluephase 20i is used exceeds the applicable RF compliance level above, the bluephase 20i should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the bluephase 20i.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strength should be less than 10 V/m.

## Recommended separation distances between portable and mobile RF communications equipment and the bluephase 20i

The bluephase 20i is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled.

The customer or the user of the bluephase 20i can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the bluephase 20i as recommended below, according to the maximum output power of the communication equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 0.4 \sqrt{P}$	80 MHz to 800 MHz $d = 0.4 \sqrt{P}$	800 MHz to 2.5 GHz $d = 0.7 \sqrt{P}$
0.01	0.04	0.04	0.07
0.1	0.13	0.13	0.22
1	0.40	0.40	0.70
10	1.3	1.3	2.2
100	4.0	4.0	7.0

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 4: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

# Ivoclar Vivadent – worldwide

## Ivoclar Vivadent AG

Bendererstrasse 2  
FL-9494 Schaan  
Liechtenstein  
Tel. +423 235 35 35  
Fax +423 235 33 60  
www.ivoclarvivadent.com

## Ivoclar Vivadent Pty. Ltd.

1 – 5 Overseas Drive  
P.O. Box 367  
Noble Park, Vic. 3174  
Australia  
Tel. +61 3 979 595 99  
Fax +61 3 979 596 45  
www.ivoclarvivadent.com.au

## Ivoclar Vivadent Ltda.

Alameda Caiapós, 723  
Centro Empresarial Tamboré  
CEP 06460-110 Barueri – SP  
Brazil  
Tel. +55 11 2424 7400  
Fax +55 11 3466 0840  
www.ivoclarvivadent.com.br

## Ivoclar Vivadent Inc.

2785 Skymark Avenue, Unit 1  
Mississauga  
Ontario L4W 4Y3  
Canada  
Tel. +1 905 238 5700  
Fax +1 905 238 5711  
www.ivoclarvivadent.com

## Ivoclar Vivadent Marketing Ltd.

Rm 603 Kuen Yang  
International Business Plaza  
No. 798 Zhao Jia Bang Road  
Shanghai 200030  
China  
Tel. +86 21 5456 0776  
Fax +86 21 6445 1561  
www.ivoclarvivadent.com

## Ivoclar Vivadent Marketing Ltd.

Calle 134 No. 7-B-83, Of. 520  
Bogotá  
Colombia  
Tel. +57 1 627 33 99  
Fax +57 1 633 16 63  
www.ivoclarvivadent.com

## Ivoclar Vivadent SAS

B.P. 118  
F-74410 Saint-Jorioz  
France  
Tel. +33 450 88 64 00  
Fax +33 450 68 91 52  
www.ivoclarvivadent.fr

## Ivoclar Vivadent GmbH

Dr. Adolf-Schneider-Str. 2  
D-73479 Ellwangen, Jagst  
Germany  
Tel. +49 (0) 79 61 / 8 89-0  
Fax +49 (0) 79 61 / 63 26  
www.ivoclarvivadent.de

## Ivoclar Vivadent Marketing

(India) Pvt. Ltd.  
503/504 Raheja Plaza  
15 B Shah Industrial Estate  
Veera Desai Road, Andheri (West)  
Mumbai, 400 053  
India  
Tel. +91 (22) 2673 0302  
Fax +91 (22) 2673 0301  
www.ivoclar-vivadent.in

## Ivoclar Vivadent s.r.l.

Via Isonzo 67/69  
40033 Casalecchio di Reno (BO)  
Italy  
Tel. +39 051 611 35 55  
Fax +39 051 611 35 65  
www.ivoclarvivadent.it

## Ivoclar Vivadent K.K.

1-28-24-4F Hongo  
Bunkyo-ku  
Tokyo 113-0033  
Japan  
Tel. +81 3 6903 3535  
Fax +81 3 5844 3657  
www.ivoclarvivadent.jp

## Ivoclar Vivadent Ltd.

12F W-Tower, 1303-37  
Seocho-dong, Seocho-gu,  
Seoul 137-855  
Republic of Korea  
Tel. +82 (2) 536 0714  
Fax +82 (2) 596 0155  
www.ivoclarvivadent.co.kr

## Ivoclar Vivadent S.A. de C.V.

Av. Insurgentes Sur No. 863,  
Piso 14, Col. Napolés  
03810 México, D.F.  
México  
Tel. +52 (55) 50 62 10 00  
Fax +52 (55) 50 62 10 29  
www.ivoclarvivadent.com.mx

## Ivoclar Vivadent Ltd.

12 Omega St, Albany  
PO Box 5243 Wellesley St  
Auckland, New Zealand  
Tel. +64 9 914 9999  
Fax +64 9 814 9990  
www.ivoclarvivadent.co.nz

## Ivoclar Vivadent Polska Sp. z o.o.

Al. Jana Pawla II 78  
00-175 Warszawa  
Poland  
Tel. +48 22 635 54 96  
Fax +48 22 635 54 69  
www.ivoclarvivadent.pl

## Ivoclar Vivadent Marketing Ltd.

Derbenevskaja Nabereshnaya 11,  
Geb. W  
115114 Moscow  
Russia  
Tel. +7 495 913 66 19  
Fax +7 495 913 66 15  
www.ivoclarvivadent.ru

## Ivoclar Vivadent Pte. Ltd.

171 Chin Swee Road  
#02-01 San Centre  
Singapore 169877  
Tel. +65 6535 6775  
Fax +65 6535 4991  
www.ivoclarvivadent.com

## Ivoclar Vivadent S.L.U.

c/ Emilio Muñoz Nº 15  
Entrada c/ Albarracín  
E-28037 Madrid  
Spain  
Tel. + 34 91 375 78 20  
Fax + 34 91 375 78 38  
www.ivoclarvivadent.es

## Ivoclar Vivadent AB

Dalvägen 14  
S-169 56 Solna  
Sweden  
Tel. +46 (0) 8 514 93 930  
Fax +46 (0) 8 514 93 940  
www.ivoclarvivadent.se

## Ivoclar Vivadent Liaison Office

:Teskiviye Mahallesi  
Sakayik Sokak  
Nisantas' Plaza No:38/2  
Kat:5 Daire:24  
34021 Sisli – Istanbul  
Turkey  
Tel. +90 212 343 08 02  
Fax +90 212 343 08 42  
www.ivoclarvivadent.com

## Ivoclar Vivadent Limited

Ground Floor Compass Building  
Feldspar Close  
Warrens Business Park  
Enderby  
Leicester LE19 4SE  
United Kingdom  
Tel. +44 116 284 78 80  
Fax +44 116 284 78 81  
www.ivoclarvivadent.co.uk

## Ivoclar Vivadent, Inc.

175 Pineview Drive  
Amherst, N.Y. 14228  
USA  
Tel. +1 800 533 6825  
Fax +1 716 691 2285  
www.ivoclarvivadent.us