Bluephase® 20i

Licence to cure

LED for every use

EN Instructions for Use

Page 6

SV Bruksanvisning

Sidan 20

DA Brugsanvisning

Side 34

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Sivu 48

NO Bruksanvisning

Side 62

NL Productinformatie

Pagina 76

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CE 0123

For dental use only!

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed dentist


Made in Austria

Ivoclar Vivadent AG, FL-9494 Schaan / Liechtenstein

ivoclar
vivadent®
clinical


Bluephase® 20i

LED for every use


Technical data

Wavelength range:	385 – 515 nm
Light intensity:	2,000 – 2,200 mW/cm²
Curing modes:	Turbo, High Power, Soft Start, Low Power
Exposure time:	for composites: 10 sec. Tetric EvoCeram / IPS Empress Direct: 5 sec.
Light probe:	10>8 mm black
Power supply:	Lithium polymer battery, approx. 45 min. capacity (curing time); approx. 2 h charging time
Operating voltage:	100-240 VAC / 50-60 Hz
Warranty	3 years (battery 1 year)


English	Svenska	Dansk	Suomi	Norsk	Nederlands	Ελληνικά	Русский
List of parts	Lampans delar	Komponentoversigt	Osaluettelo	Deleliste	Onderdelenoverzicht	Κατάλογος μερών	Перечень комплектующих
① Handpiece	Handstycke	Håndstykke	Käsikappale	Håndstykke	Handstuk	χειρολαβή	Наконечник
② Start button	Startknapp	Starttast	Start-näppäin	Starttast	Startknop	Κουμπί έναρξης	Кнопка «Старт»
③ Display	Display	Display	Näyttö	Display	Display	Οθονήξ	Дисплей
④ Program selection button	Programvalsknapp	Programvælgertast	Ohjelman valintanäppäin	Programvalgtast	Programmakeuzeknop	Κουμπί επιλογής προγράμματος	Кнопка выбора программы
⑤ Time selection button	Tidvalsknapp	Tidvælgertast	Ajan valintanäppäin	Tidsvalgtast	Tijdstellingsknop	Κουμπί επιλογής χρόνου	Кнопка выбора времени
⑥ Volume button	Volymknapp	Volumenknop	Äänenvoimakkuuden valintanäppäin	Volumtast	Volumeknop	Κουμπί έντασης	Кнопка громкости
⑦ Power indicator	Strömindikator	Tændt-indikator (Power)	Virtanäyttö	Driftsindikator	Powerindicator	Δείκτης ισχύος	Индикатор питания
⑧ Charging base	Laddningsstation	Ladestation	Latausasema	Ladestasjon	Laadstation	Βάση φόρτισης	Зарядный блок
⑨ Power cord	Nätkabel	Netkabel	Verkkojohto	Strømledning	Netsnoer	Καλώδιο τροφοδοσίας ηλεκτρικού ρεύματος	Сетевой кабель
⑩ Power pack	Transformator	Netapparat	Voimaosa	Nettadapter	Netvoedingseenheid	Τροφοδοτικό	Сетевой блок
⑪ Battery	Batteri	Batteri	Akku	Batteri	Batterij	Μπαταρία	Акумулятор
⑫ Light probe 10>8 mm	Ljusledare 10>8 mm	Lysleder 10>8 mm	Valojohdin 10>8 mm	Lysleder 10>8 mm	Lichtgeleider 10>8 mm	Ρύγχος φωτός 10>8 χιλ.	Световод 10>8 mm
⑬ Anti-glare shield	Avbländningsskydd	Blændingskærm	Häikäisy suoja kilpi	1Blendingsbeskyttelse (-skjerm)	Antiverblindingschermpje	Αντιθαμβωτική ασπίδα	Защитный экран
⑭ Anti-glare cones	Antibländningskydd (koner)	Anti-blændingskegle	Häikäisy suoja (-kartio)	Blendingsbeskyttelse (-kjegle)	antiverblindingskapje	Αντιθαμβωτικοί κώνοι	Защитный козырек

	Indicators on the charging base	Indikatorer på laddningsstationen	Indikatorer på ladestationen	Latausaseman indikaattorit	Indikatorer på ladestasjonen	Indicatoren op het laadstation	Δείκτες στη βάση φόρτισης	Индикация на зарядном блоке
	Indicator lights up in blue – the charging base is connected to the power supply	Indikatorn lyser blått - laddningsstationen är kopplad till strömkällan.	Indikatoren lyser blått – ladestationen er forbundet til forsyningsnettet	Indikaattorissa palaa sininen valo - latausasema on liitetty virtalähteeseen	Indikatoren lyser blått – ladestasjonen er koblet til strømforsyningen	Het indicatielampje brandt blauw – het laadstation is aangesloten op de stroomvoorziening	ο δείκτης ανάβει σε μπλε χρώμα – η βάση φόρτισης είναι συνδεδεμένη με την τροφοδοσία ρεύματος	Индикатор горит синим светом – зарядный блок подключен к сети
	Indicators on the handpiece	Indikatorer på handstycket	Indikatorer på håndstykket	Käsikappaleen indikaattorit	Indikatorer på håndstykket	Indicatoren op het handstuk	Δείκτες στη χειρολαβή	Индикация на наконечнике
	Curing program and operating status	Ljushärdningsprogram och driftsstatus	Belysningsprogram og driftstilstand	Kovetusohjelma ja käyttötila	Herdeprogram og driftstatus	Belichtingsprogramma en status	Πρόγραμμα πολυμερισμού και κατάσταση λειτουργίας	Программа полимеризации и состояние готовности
	Selected indications and curing time	Vald indikation och ljushärdningstid	Valgte indikatorer og belysningstid	Valittuja indikaatioita ja kovetusaika	Valgte indikeringer og herdetid	Gekozen indicaties en belichtingsduur	Επιλεγμένοι δείκτες και χρόνος πολυμερισμού	Показания и время полимеризации
	Charging status of the battery	Batteriets laddningsstatus	Batteriets tilstand	Akun latausaika	Batteriets ladestatus	Laadtoestand van de batterij	Κατάσταση φόρτισης της μπαταρίας	Состояние зарядки аккумулятора
	Display is in stand-by mode	Display är i standby-läge	Displayet er i stand-by-modus	Näyttö on valmiustilassa.	Displayet er i standby-modus	Display staat in de stand-by modus	Η οθόνη βρίσκεται σε κατάσταση αναμονής	Дисплей в режиме Stand-by
	Signs and symbols	Tecken och symboler	Symbolforklaring	Merkit ja symbolit	Tegnforklaring	Verklaring van de symbolen	Σήματα και σύμβολα	Символы и обозначения
	Non-permissible use	Kontraindikation	Ikke-tilladt brug	Käyttö ei sallittu	Ikke tillatt bruk	Niet-toegestane toepassing	Μη επιτρεπτή χρήση Διπλή μόνωση	Недопустимое применение
	Double insulation (Apparatus complies with safety class II)	Dubbelisolering (Apparaten överensstämmer med säkerhetsklass II)	Dobbeltisoleret (apparat i beskyttelsesklasse II)	Kaksoiserite (laite vastaa turvaluokan II määräyksiä)	Dobbel isolering (apparatet svarer til sikkerhetsklasse II)	Dubbel geïsoleerd (beschermingsklasse II)	(Η συσκευή συμμορφώνεται με την κατηγορία ασφαλείας II)	Двойная изоляция (аппарат соответствует II классу безопасности)
	Protection against electrical shock (BF type apparatus)	Skydd mot elektrisk chock (BF apparat)	Beskyttelse mod elektriske stød (BF-luokan laite)	Suojaus sähköiskulta (BF-luokan laite)	Beskyttelse mot elektrisk stød (apparat av type BF)	Beveiligd tegen elektrische schokken (type BF)	Προστασία από ηλεκτροπληξία (συσκευή τύπου BF)	Защита от поражения электротоком (аппарат типа BF)
	Observe Instructions for Use	Observera bruksanvisningen	Bemærk brugsanvisningen	Noudata käyttöohjeita	Følg brugsanvisningen	Raadpleeg de gebruiksaanwijzing	Τηρείτε τις οδηγίες χειρισμού	Смотри инструкцию по эксплуатации
	Observe Instructions for Use	Observera bruksanvisningen	Bemærk brugsanvisningen	Noudata käyttöohjeita	Følg brugsanvisningen	Raadpleeg de gebruiksaanwijzing	Τηρείτε τις οδηγίες χειρισμού	Смотри инструкцию по эксплуатации
	Caution	Försiktigt!	Forsigtig	Varoitus	Forsiktig	Let op!	Προσοχή	Осторожно!
	The apparatus must not be disposed of in the normal domestic waste	Apparaten får inte kastas bland hushållssopor	Apparatet må ikke bortskaffes via det normale husholdningsaffald	Laitetta ei saa hävittää tavallisen talousjätteenä	Apparatet må ikke kastes i det normale husholdningsavfallet.	Het apparaat mag niet als normaal huishoudelijk afval worden afgevoerd	Η συσκευή δεν πρέπει να απορριπτεται με τα κανονικά οικιακά απορρίμματα	Аппарат не должен утилизироваться с обычными бытовыми отходами
	Recyclable	Kan materialåtervinnas	Genanvendelig	Kierrätettävä	Resirkulerbar	Recyclebaar	Ανακυκλώσιμο	Повторное использование после утилизации

TURBO Program
Consistently high light intensity for the polymerization of restorative and cementation materials for direct and indirect restorations.

Light intensity	2,000 – 2,200 mW/cm²
Exposure time for composites	10 s
Exposure time for Tetric EvoCeram / IPS Empress Direct	5 s
Filling materials	Curing time
Composite	
• 2 mm ¹⁾	
IPS Empress Direct / IPS Empress Direct Flow / Tetric EvoCeram / Tetric EvoFlow / Tetric / Tetric Basic White / Tetric N-Ceram / Tetric N-Flow	5 s
Heliomolar / Heliomolar HB / Heliomolar Flow / Others	2x 5 s
• 4 mm ²⁾	
Tetric EvoCeram Bulk Fill	5 s
Tetric Basic White	2x 5 s
Compomer ³⁾	
Compoglass F / Compoglass Flow	2x 5 s
Indirect restorations / Luting materials	
Variolink II Base ⁴⁾ / Variolink Veneer / Variolink N Base ⁴⁾ / Variolink N Clear Veneer / Variolink Ultra ⁵⁾ / Dual Cement ⁵⁾ / Variolink II ⁵⁾	per mm ceramic: 1x 5 s per segment
Miscellaneous	
Heliosit Orthodontic	2x 5 s
Telio Add-On Flow	2x 5 s
Telio Stains	5 s
IPS Empress Direct Color	5 s
IPS Empress Direct Opaque	2x 5 s

HIGH POWER Program
Consistently high light intensity for the polymerization of restorative and cementation materials for direct and indirect restorations.

Filling materials	Curing time
Composite	
• 2 mm ¹⁾	
IPS Empress Direct / IPS Empress Direct Flow / Tetric EvoCeram / Tetric EvoFlow / Tetric / Tetric Basic White / Tetric N-Ceram / Tetric N-Flow	10 s
Heliomolar / Heliomolar HB / Heliomolar Flow / Others	15 s
• 4 mm ²⁾	
Tetric EvoCeram Bulk Fill / Tetric Basic White	10 s
Compomer ³⁾	
Compoglass F / Compoglass Flow	20 s
Indirect restorations / Luting materials	
Variolink II Base ⁴⁾ / Variolink Veneer / Variolink N Base ⁴⁾ / Variolink N Clear Veneer / Variolink Ultra ⁵⁾ / Dual Cement ⁵⁾ / Variolink II ⁵⁾	per mm ceramic: 10 s per segment
Miscellaneous	
Helioseal / Helioseal F / Helioseal Clear	10 s
Monopaque	20 s
MultiCore Flow / Multicore HB	20 s
Heliosit Orthodontic	10 s
Telio Add-On Flow	15 s
Telio Stains	10 s
IPS Empress Direct Color	10 s
IPS Empress Direct Opaque	20 s

LOW POWER Program
Reduced light intensity with reduced heat development for the polymerization of adhesives, liners, and restorative materials in areas near the pulp when restoring Class V cavities.

Adhesive	Curing time
AdheSE / AdheSE One F	
ExGTE / ExGTE DSC	10 s
Tetric N-Bond / Tetric N-Bond Self-Etch	
Heliobond (Syntac)	
Temporary materials	
Telio CS Inlay/Onlay	10 s
Systemp.inlay/onlay	10 s
Fermit / Fermit N	10 s
Telio CS Link / Systemp.link	20 s per segment
Miscellaneous	
Heliosit Orthodontic	20 s
Vivaglass Liner	20 s

SOFT START Program
Step-by-step increase of the light intensity with reduced shrinkage stress and reduced heat development for the polymerization of restorative materials.

Filling materials	Curing time
Composite	
• 2 mm ¹⁾	
IPS Empress Direct / IPS Empress Direct Flow / Tetric EvoCeram / Tetric EvoFlow / Tetric / Tetric Basic White / Tetric N-Ceram / Tetric N-Flow	15 s
Heliomolar / Heliomolar HB / Heliomolar Flow / Others	20 s
• 4 mm ²⁾	
Tetric EvoCeram Bulk Fill / Tetric Basic White	15 s
Compomer ³⁾	
Compoglass F / Compoglass Flow	20 s

1) Applies to a maximum layer thickness of 2 mm and provided that the Instructions for Use of the respective material do not state any other recommendation (might be the case, e.g. with dentin shades)
2) Applies to a maximum layer thickness of 4 mm and provided that the Instructions for Use of the respective material do not state any other recommendation (might be the case, e.g. with dentin shades)
3) Applies to a maximum layer thickness of 3 mm 4) Applies to light-curing (use of the base past only) 5) Applies to dual-curing

Dear Customer

Optimum polymerization is an important requirement for all light cured materials in order to consistently produce high quality restorations. The polymerization light selected also plays a decisive role in this respect. Therefore, we would like to thank you for having purchased Bluephase® 20i. Bluephase 20i is a high-quality medical device which has been designed according to the latest standard of science and technology in compliance with the relevant industry standards.

These Instructions for Use will help you to safely start up the device, make full use of its capabilities and ensure a long service life.

Should you have any further questions, please do not hesitate to contact us (see addresses on the back cover).

Your
Ivoclar Vivadent Team

Table of Contents

Safety	
Intended use	8
Indications	8
Safety notes	8
Start-up	
Charging base	10
Handpiece	11
Battery	11
Click & Cure corded operation	12
Operation	
Selecting the curing program and the curing time	14
Cure Memory function	16
Start	
Acoustic signals	16
Light intensity	16
Measuring the light intensity	16
Maintenance and cleaning	
Housing	17
Light probes	17
Battery contacts	17
Disposal	17
What if ...	18
Warranty / Procedure in case of repair	19
Product specifications	19

Safety

Intended use



Bluephase 20i is an LED polymerization light that produces energy-rich blue light. It is used for the polymerization of light-curing dental materials immediately at the dental unit. The intended use also includes the observation of the notes and regulations in these Instructions for Use.

Indication

With its "Polywave®" broadband spectrum, Bluephase is suitable for the polymerization of all light-curing dental materials curing in the wavelength range of 385–515 nm. These materials include restoratives, bonding agents/adhesives, bases, liners, fissure sealants, temporary materials as well as luting materials for brackets and dental-lab materials such as those used in the fabrication of ceramic inlays.

Safety notes

Bluephase 20i is an electronic device and a medical product which is subject to IEC 60601-1 (EN 60601-1) and EMC directives, as well as the 93/42/EEC Medical Device Directive. The appliance complies with the relevant EU regulations and is classified as an LED Class 2 product.

CE 0123

The curing light has been shipped from the manufacturer in a safe and technically sound condition. In order to maintain this condition and to ensure risk-free operation, the notes and regulations in these Instructions for Use have to be observed. To prevent damage to equipment and risks for patients, users and third parties, the following safety instructions have to be observed.

Contraindications

Materials, the polymerization of which is activated outside the wavelength range of 385–515 nm (no materials known to date). If you are not sure about certain products, please ask the manufacturer of the corresponding material.

Do not charge or use the appliance near flammable or combustible substances.

Portable and mobile high-frequency communication devices may interfere with medical equipment. The use of mobile phones during operation is not allowed.

Caution – The use of controls or adjustment devices or the performance of procedures other than those specified herein may result in hazardous radiation exposure.

Usage and liability

- Bluephase 20i must only be employed for the intended use. Any other uses are contraindicated. Liability cannot be accepted for damage resulting from misuse or failure to observe the Operating Instructions.
- The user is responsible for testing Bluephase 20i for its use and suitability for the intended purposes. This is particularly important if other equipment is used in the immediate vicinity of Bluephase at the same time.
- Only original spare parts and accessories from Ivoclar Vivadent (see Accessories) should be used. The manufacturer does not accept any liability for damage resulting from the use of other spare parts or accessories.
- The light probe is an applied part and may warm up to a maximum of 45 °C (113 °F) at the interface to the hand-piece during operation.

Operating voltage

Before switching on, make sure that

- a) the voltage indicated on the rating plate complies with the local power supply and
- b) the unit has acquired the ambient temperature.

If the battery or power pack are used separately, e.g. during start-up or Click & Cure corded operation, contact with patients or third parties must be prevented. Do not touch the exposed contacts of the battery or connection plug (power pack).

Assumption of impaired safety

If it is assumed that safe operation is no longer possible, the power must be disconnected and the battery removed to avoid accidental operation. This may be the case, for example, if the apparatus is visibly damaged or no longer works correctly. Complete disconnection from the power supply is only ensured when the power cord is disconnected from the power source.

Eye protection

Direct or indirect exposure of the eyes must be prevented. Prolonged exposure to the light is unpleasant to the eyes and may result in injury. Therefore, using the supplied anti-glare cones is recommended. Individuals who are generally sensitive to light, who take photosensitizing drugs, have undergone eye surgery, or people who work with the apparatus or in its vicinity for long periods of time should not be exposed to the light of this device. Therefore, they should wear protective goggles (orange) that absorb light below 515 nm.

Battery

Caution: Use only original spare parts, particularly Ivoclar Vivadent batteries and charging bases. Do not short circuit battery. Do not store at temperatures above 40 °C / 104 °F (or 60 °C / 140 °F for a short period). Always store batteries charged. The storage period must not exceed 6 months. May explode if disposed of in fire.

Please note that the lithium polymer battery may react with explosion, fire, and smoke development in case of inappropriate handling and mechanical damage. Damaged lithium polymer batteries must no longer be used.

The electrolytes and electrolyte fumes released during explosion, fire, and smoke development are toxic and corrosive. In case of contact with the eyes and skin, immediately rinse with copious amounts of water. Avoid the inhalation of the fumes. See a physician immediately in case of indisposition.

Heat development (burn hazard)

As it is the case with all high-performance lights, the high light intensity results in a certain heat development. Prolonged exposure of areas near the pulp and soft tissues may result in irreversible or reversible damage. Therefore, this high-performance curing light must only be operated by trained professionals.

The recommended curing times have to be observed. Uninterrupted curing times of more than 5 seconds in the Turbo-Program on the same tooth surface, as well as direct contact with the gingiva, oral mucous membrane or skin have to be prevented. If required, polymerize in intermittent intervals of at most 5 seconds each. Prevent heat accumulation when polymerizing with a rubber dam in place. If irradiation of soft tissue cannot be avoided, operate with reduced light intensity (use Low Power program). Polymerize indirect restorations at intermittent intervals of 10 seconds in the High-Power program and 1x 5 seconds each in the Turbo program or use external cooling with an air stream.

The instructions regarding curing programs and curing times must be observed (see Selecting the curing program and the curing time). In addition, the light emission window must be placed directly over the material to be cured at all times (e.g. by holding it in place using a finger).

It is not recommended to use the High Power/Turbo program for sensitive patients. In this case, operate with reduced light intensity (use Low Power program).

Start-up

Delivery form

Check the delivery for completeness and any possible transportation damage (see Delivery form). If components are missing or damaged, immediately contact your dealer or your Service Center.

Charging base

Before you switch on the device, make sure that the voltage mentioned on the rating plate complies with your local power supply. The rating plate is attached to the underside of the charging station.

Place the charging base on a suitable, flat table top. Remove the protective foil from the connection plug of the power pack. Slide the connection plug of the power pack into the socket on the underside of the charging base. Tilt it slightly and apply light pressure until you hear and feel it snap into place.


Connect the power cord with the power supply and the power pack.

The **Power ON** indicator on the left side of the housing lights up in blue (see Indicators on the charging base).

Handpiece

Before you insert the light probe, remove the protective foil on the corresponding opening of the handpiece.

Slightly turn the light probe while attaching it to the handpiece.

After that, mount the anti-glare shield on the light probe.

Battery

We recommend fully charging the battery before the first use.

Slide the battery straight into the handpiece until you hear and feel it click into place.

Gently place the handpiece in the corresponding cradle in the charging base. If a sleeve is used, please remove it before you charge the battery. If the battery is fully charged, it features a curing capacity of approximately 45 minutes.

Tip

If possible, use the light always with a fully charged battery. This will prolong its service life. It is therefore recommended to place the handpiece into the charging base after each appointment. The charging time for empty batteries is 2 hours.

Click & Cure corded operation

Bluephase 20i can be used in corded operation at any time, but particularly when the battery is completely empty.

For this purpose, remove the battery from the handpiece. Then remove the power pack from the underside of the charging base. Do not pull on the power cord.

Insert the connection plug straight into the handpiece until you hear and feel it click into place.

During corded operation, the charging base cannot charge the battery, since it is not connected to a power source.

The apparatus is only completely disconnected from the power supply if the power cord has been pulled out of the plug-in socket.

Charging status of the battery

With the handpiece switched on, the current charging status is shown on the display as follows:

Battery fully charged

(curing capacity of approximately 45 minutes).

Battery half full

Reserve


(The last bar in the battery indicator is red in the reserve mode. The battery has to be charged as soon as possible.)


When the battery is being charged, „Charging Battery“ briefly appears in the display before the display returns to the stand-by mode.

If the battery is completely empty, the handpiece automatically switches to the stand-by mode. The light can no longer be called up and the curing program and curing time can no longer be set. However, the handpiece can be used in the Click & Cure corded operation.

Since the battery is a consumable, it has to be replaced after its typical life cycle has expired after 2½ years. See battery label for the age of the battery.

Operation

Disinfect or autoclave contaminated surfaces of the curing light as well as light probes and anti-glare cones before each use.

Furthermore, make sure that the stipulated light intensity permits adequate polymerization. For this purpose, check the light probe for contamination and damage as well as the light intensity at regular intervals (e.g. with the Ivoclar Vivadent radiometer "Bluephase Meter").

Selecting the curing program and the curing time

The curing program and the curing time can be individually set. Bluephase 20i is equipped with the following 4 curing programs for the different indications. The desired curing program is selected with the program selection buttons. The display changes accordingly (see Indicators on the handpiece). The device comes equipped with the following preset programs:

Factory settings

TURBO

5 seconds

**HIGH
(High Power)**

10 seconds

**LOW
(Low Power)**

10 seconds

**SOFT
(Soft start)**

15 seconds



The intended curing time is selected using the time selection buttons. Users may choose between 5, 10, 15, 20 and 30 seconds.

Observe the Instructions for Use of the material applied when selecting the curing time.

The curing recommendations for composite materials apply to all shades and, if not mentioned otherwise in the Instructions for Use, to a maximum layer thickness of 2 mm. Generally, these recommendations apply to situations where the emission window of the light probe is placed directly over the material to be polymerized. Increasing the distance between the light source and the material will require the curing time to be extended accordingly. For instance, if the distance to the material is 8 mm, the effective light output is reduced by approx. 50%. In this case, the recommended curing time has to be doubled.

- 1) The information provided herein applies to the light probe 10>8 mm supplied in the delivery form.
- 2) The information regarding heat development and burn hazard must be taken into consideration (see Safety notes).

Curing program



TURBO Program

Consistently high light intensity for the polymerization of restorative and cementation materials for direct and indirect restorations.

Light intensity	2,000 – 2,200 mW/cm ²
Exposure time for composites	10 s
Exposure time for Tetric EvoCeram / IPS Empress Direct	5 s



Filling materials	Curing time
Composite	
• 2 mm ¹⁾	
IPS Empress Direct / IPS Empress Direct Flow / Tetric EvoCeram / Tetric EvoFlow / Tetric / Tetric Basic White / Tetric N-Ceram / Tetric N-Flow	5 s
Heliomolar / Heliomolar HB / Heliomolar Flow / Others	2x 5 s
• 4 mm ²⁾	
Tetric EvoCeram Bulk Fill	5 s
Tetric Basic White	2x 5 s
Compomer ³⁾	
Compoglass F / Compoglass Flow	2x 5 s
Indirect restorations / Luting materials	
Variolink II Base ⁴⁾ / Variolink Veneer / Variolink N Base ⁴⁾ / Variolink N Clear Veneer / Variolink Ultra ⁴⁾ / Dual Cement ⁴⁾ / Variolink II ⁴⁾	per mm ceramic: 1x 5 s per segment
Miscellaneous	
Heliosit Orthodontic	2x 5 s
Telio Add-On Flow	2x 5 s
Telio Stains	5 s
IPS Empress Direct Color	5 s
IPS Empress Direct Opaque	2x 5 s

HIGH POWER Program

Consistently high light intensity for the polymerization of restorative and cementation materials for direct and indirect restorations.



Filling materials	Curing time
Composite	
• 2 mm ¹⁾	
IPS Empress Direct / IPS Empress Direct Flow / Tetric EvoCeram / Tetric EvoFlow / Tetric / Tetric Basic White / Tetric N-Ceram / Tetric N-Flow	10 s
Heliomolar / Heliomolar HB / Heliomolar Flow / Others	15 s
• 4 mm ²⁾	
Tetric EvoCeram Bulk Fill / Tetric Basic White	10 s
Compomer ³⁾	
Compoglass F / Compoglass Flow	20 s
Indirect restorations / Luting materials	
Variolink II Base ⁴⁾ / Variolink Veneer / Variolink N Base ⁴⁾ / Variolink N Clear Veneer / Variolink Ultra ⁴⁾ / Dual Cement ⁴⁾ / Variolink II ⁴⁾	per mm ceramic: 10 s per segment
Miscellaneous	
Helioseal / Helioseal F / Helioseal Clear	10 s
Monopaque	20 s
MultiCore Flow / Multicore HB	20 s
Heliosit Orthodontic	10 s
Telio Add-On Flow	15 s
Telio Stains	10 s
IPS Empress Direct Color	10 s
IPS Empress Direct Opaque	20 s

LOW POWER Program

Reduced light intensity with reduced heat development for the polymerization of adhesives, liners, and restorative materials in areas near the pulp when restoring Class V cavities.

Adhesive	Curing time
AdheSE / AdheSE One F	
ExoTE / ExoTE DSC	10 s
Tetric N-Bond / Tetric N-Bond Self-Etch	
Heliobond (Syntac)	
Temporary materials	
Telio CS Inlay/Onlay	10 s
System.inlay/onlay	10 s
Fermit / Fermit N	10 s
Telio CS Link / System.Link	20 s per segment
Miscellaneous	
Heliosit Orthodontic	20 s
Vivaglass Liner	20 s

SOFT START Program

Step-by-step increase of the light intensity with reduced shrinkage stress and reduced heat development for the polymerization of restorative materials.

Filling materials	Curing time
Composite	
• 2 mm ¹⁾	
IPS Empress Direct / IPS Empress Direct Flow / Tetric EvoCeram / Tetric EvoFlow / Tetric / Tetric Basic White / Tetric N-Ceram / Tetric N-Flow	15 s
Heliomolar / Heliomolar HB / Heliomolar Flow / Others	20 s
• 4 mm ²⁾	
Tetric EvoCeram Bulk Fill / Tetric Basic White	15 s
Compomer ³⁾	
Compoglass F / Compoglass Flow	20 s

1) Applies to a maximum layer thickness of 2 mm and provided that the Instructions for Use of the respective material do not state any other recommendation (might be the case, e.g. with dentin shades)

2) Applies to a maximum layer thickness of 4 mm and provided that the Instructions for Use of the respective material do not state any other recommendation (might be the case, e.g. with dentin shades)

3) Applies to a maximum layer thickness of 3 mm 4) Applies to light-curing (use of the base past only) 5) Applies to dual-curing

Cure Memory function

The last settings used, together with the combination of curing program and curing time, are automatically saved.

Start

The light is switched on by means of the start button. Once the selected curing time has elapsed, the curing program is automatically terminated. If desired, the light can be switched off before the set curing time has elapsed by pressing the start button again. The fan is activated simultaneously to the light. Once the curing time has elapsed, the fan continues to run for a certain time to cool the apparatus. The battery must not be removed as long as the fan is still running.

Acoustic signals

Acoustic signals can be heard for the following functions:

- Start (Stop)
- Every 10 seconds
- Program change
- Curing time change
- Connecting the battery to the charging base
- Inserting battery
- Error message

If desired, the volume of the acoustic signals can be adjusted.

For this purpose, press the blue volume button (see Product overview) while the light is off and reduce the volume or turn the acoustic signals off with the left time selection button. If you want to turn the acoustic signals back on or increase the volume, press the blue volume button and then the right time selection button.

Light intensity

The light intensity is maintained at a consistent level during operation. If the supplied 10>8 mm light probe is used, the light intensity has been calibrated to 2,000 - 2,200 mW/cm².

The use of a light probe other than the one provided has a direct influence on the light intensity emitted.

In parallel-walled light probes (10 mm), the diameter is equal at both ends. In focussing light probes (10>8 mm light probe, Pin-Point light probe 6>2 mm), the diameter of the rear end is larger than that of the light emission window. The incident blue light is thus bundled to a smaller surface, which increases the light intensity emitted.

Pin-Point light probes are suitable for the polymerization of confined areas, such as the attachment of veneers prior to excess removal. For thorough curing, it is necessary to change the light probe.

Measuring the light intensity

The light intensity of the Bluephase 20i and the supplied 10>8 mm light probe can be checked, for instance, by means of the Bluephase Meter.

If the measured value does not correspond with the expected light intensity, proceed as follows:

- Check the selected curing program
- Clean light sensor if it is contaminated
- Remove the light probe and clean the light emission window of the handpiece with a cotton swab dipped in alcohol
- Clean light probe if it is contaminated (see Maintenance and cleaning)
- Replace a damaged light probe with a new one

If these measures do not improve the results, please contact your dealer or your local Service Center.

Maintenance and cleaning

For reasons of hygiene, we recommend using a disposable protective sleeve for each patient. Make sure to fit the protective sleeve closely to the light probe. Disinfect contaminated surfaces of the device and anti-glare cones (FD 366/Dürr Dental, Incidin Liquid/Ecolab) and sterilize the light probe before each use if disposable protective sleeves are not used. Make sure that no liquids or other foreign substances enter the handpiece, charging base and particularly the power pack during cleaning (risk of electrical shock). Disconnect the charging base from the power source when cleaning it.

Housing

Wipe the handpiece and handpiece holder with a customary aldehyde-free disinfecting solution. Do not clean with highly aggressive disinfecting solutions (e.g. solutions based on orange oil or with an ethanol content of more than 40%), solvents (e.g. acetone), or pointed instruments, which may damage or scratch the plastic. Clean dirty plastic parts with a soapy solution.

Light probe

Before cleaning and/or disinfecting the light probe, pretreat it. This applies to both automated and manual cleaning and disinfection.

Pretreatment

- Remove gross contamination immediately after use or within 2 hours. For this purpose, thoroughly rinse the light probe under running water (for at least 10 seconds). Alternatively, use a suitable aldehyde-free disinfecting solution to prevent the fixation of blood.
- To remove contamination manually, use a soft brush or soft cloth. Partially polymerized composite can be removed with alcohol and a plastic spatula, if necessary. Do not use sharp or pointed objects, as they may scratch the surface.

Cleaning and disinfection

For cleaning, immerse the light probe in a cleaning solution and ensure that it is sufficiently covered with liquid (ultrasound or careful brushing with a soft brush may support the effect). A neutral-enzymatic cleaning agent is recommended.

When cleaning and disinfecting, please make sure that the agents used are free of:

- organic, mineral and oxidizing acids (the minimum admissible pH value is 5.5)
- alkaline solution (the maximum admissible pH value is 8.5)
- oxidizing agent (e.g. hydrogen peroxide)

Afterwards, remove the light probe from the solution and thoroughly rinse it under running water (for at least 10 seconds). Cleaning in a thermal disinfectant is an effective alternative.

Sterilization

Thorough cleaning and disinfecting is imperative to ensure that the subsequent sterilization is effective. Use only autoclave sterilization for this purpose. The sterilization time (exposure time at sterilization temperature) is 4 minutes at 134 °C (273 °F); pressure should be 2 bar (29 psi). Dry the sterilized light probe using either the special drying program of your steam autoclave or hot air. The light probe has been tested for up to 200 sterilization cycles.

After that, check the light probe for damage. Hold it against light. If individual segments appear black, glass fibres are broken. If this is the case, replace the light probe with a new one.

Battery contacts

To ensure reliable conductivity at all times, keep the battery contacts free from possible contamination (e.g. composite residue). For this purpose, clean the affected contacts regularly in the course of the usual wipe disinfection (after each patient).




Disposal

The curing light must not be disposed of as urban waste. Dispose unserviceable batteries and polymerization lights according to the corresponding legal requirements in your country. Batteries must not be incinerated.

What if ...

Symbol	Causes	Error rectification
	Electronic defect in the handpiece or battery	Remove and reinsert the battery. If the error persists, replace the battery with the power pack (Click & Cure). If the error persists, please contact your dealer or your local Service Center.
	Apparatus is overheated or under-cooled	Allow the apparatus to cool down (or assume room temperature if it is undercooled) and try again after a certain time. If the error persists, please contact your dealer or your local Service Center.
	Electronic defect in the battery	Remove and reinsert the battery. If the error continues to be displayed, place the apparatus into the charging base. If the error persists nonetheless, replace the battery with the power pack (Click & Cure). Please contact your dealer or your local Service Center.
	Battery empty	Place the handpiece into the charging base. The corresponding message („Charging“) will be displayed. If the battery is not charged, the contacts have to be cleaned. If the error persists, please contact your dealer or your local Service Center.
	Charging contacts contaminated Battery in safety mode (integrated protective circuit)	Check if the handpiece has been correctly placed into the charging base. Clean battery contacts. Information: „Charging Battery“ briefly appears in the display before the display returns to the stand-by mode. If cleaning the battery contacts does not resolve the problem, remove the battery from the apparatus. Charge the battery separately from the handpiece in the charging base for approximately 10 minutes. If the error persists, please contact your dealer or your local Service Center.
	Power pack not connected or defective	Check the battery contacts and if the power pack is correctly positioned in the charging base or if the power pack is connected to the power supply by means of the power cord (display on the power pack lights up in green if it works correctly). If the error persists, please contact your dealer or your local Service Center.

Warranty / Procedure in case of repair

The warranty period for Bluephase 20i is 3 years from the date of purchase (battery 1 year).

Malfunctions resulting from faulty material or manufacturing errors are repaired free of charge during the warranty period. The warranty does not provide the right to recover any material or non-material damage other than the ones mentioned. The apparatus must only be used for the intended purposes. Any other uses are contraindicated. The manufacturer does not accept any liability resulting from misuse and warranty claims cannot be accepted in such cases. This is particularly valid for:

- damage resulting from improper handling, especially incorrectly stored batteries (see Technical data: Transportation and storage conditions).
- damage to components resulting from wear under standard operating conditions (e.g. battery).
- damage resulting from external influences, e.g. blows, drop to the floor.
- damage resulting from incorrect set-up or installation.
- damage resulting from connecting the unit to a power supply, the voltage and frequency of which do not comply with the ones stated on the rating plate.
- damage resulting from improper repairs or modifications that have not been carried out by certified Service Centres.

In case of a claim under warranty, the complete apparatus (handpiece, charging base, power cord and power pack) must be returned, carriage paid, to the dealer or directly to Ivoclar Vivadent, together with the purchase document. Use the original packaging with the corresponding cardboard inserts for transportation.

Repair work may only be carried out by a certified Ivoclar Vivadent Service Center. In case of a defect that cannot be rectified, please contact your dealer or your local Service Center (see addresses on the reverse side). A clear description of the defect or the conditions under which the defect occurred will facilitate locating the problem. Please enclose this description when returning the apparatus.

Product specifications

Delivery form

1 Charging base, 1 Power cord, 1 Power pack, 1 Handpiece 1 Battery (Li-Po),
1 10>8 mm light probe, 1 Anti-glare shield, 3 Anti-glare cones, Sleeves,
1 Set of Instructions for Use

Accessories

The following accessories are available for Bluephase 20i:

REF	Description
627389	Light probe 10>8 mm, black (G2)
608538	Pin-point light probe 6/2 mm, black (G2)
551756	Anti-glare cones
592496	Anti-glare shield
627299	Bluephase 20i charging base (G2)
627300	Bluephase 20i battery (G2)
627298	Bluephase 20i handpiece (G2)
607922	Bluephase Meter
608554	Bluephase sleeves (G2)

Technical data

Operating voltage charging base	5 VDC
Operating voltage handpiece	3.7 VDC with battery 5 VDC with power pack
Power Pack	100-240 VAC / 50-60 Hz / max. 0.4 A / Output 5 VDC TYPE 15.2630 Manufacturer: Friwo

Operating conditions

Temperature	+10 °C to +35 °C (+50 to +95°F)
Relative humidity	30 % to 75%
Ambient pressure	700 hPa to 1060 hPa
Dimensions of the charging base	L=205mm; W=150mm; H=85mm
Weight of the charging base	250 g
Charging time	approx. 2 h (with the battery empty)
Power supply of the handpiece	Li-Po battery
Max. battery time	approx. 45 min. (with a new, fully charged battery)
Light source	Polywave® LED
Wavelength range	385–515 nm
LED-Class 2	LED LIGHT DO NOT STARE INTO BEAM CLASS 2 LED PRODUCT Max. 5.11 mW / Wavelength 385 - 515 nm
Light intensity	2,000 - 2,200 mW/cm²
Operation	5 min. on / 6 min. off (intermittently)
Dimensions of the handpiece	L=260mm; W=42mm; H=120mm
Weight of the handpiece	(incl. battery) 225 g

Transportation and storage conditions:

Temperature	-20 °C to +60 °C (-4 to 140°F)
Relative humidity	10% to 75%
Ambient pressure	500 hPa to 1060 hPa

Bluephase 20i has to be stored in closed, roofed rooms.

Protect the device from severe jarring.

Battery

- Do not store at temperatures above 40 °C / 104 °F (or 60 °C / 140 °F for a short period). Recommended storage temperature 15-30 °C (59-86 °F).
- Always store the battery charged and
- not for longer than 6 months.

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Appendix

Guidance and manufacturer’s declaration - electromagnetic emission

The Bluephase 20i is intended for use in the electromagnetic environment specified below. The customer or the user of the Bluephase 20i should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Bluephase 20i uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Bluephase 20i is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	P < 75W (pass without test)
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	P < 75W (pass without test)

Guidance and manufacturer’s declaration - electromagnetic immunity

The Bluephase 20i is intended for use in the electromagnetic environment specified below. The customer or the user of the Bluephase 20i should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of typical commercial or dental environment.
Surge IEC 61000-4-5	± 1 kV line - line ± 2 kV line - earth	± 1 kV line - line no prot. earth	Mains power quality should be that of typical commercial or dental environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U _i (>95 % dip in U _i) for 0.5 cycle 40 % U _i (60 % dip in U _i) for 5 cycles 70 % U _i (30 % dip in U _i) for 25 cycles <5 % U _i (>95 % dip in U _i) for 5 sec	<5 % U _i (>95 % dip in U _i) for 0.5 cycle 40 % U _i (60 % dip in U _i) for 5 cycles 70 % U _i (30 % dip in U _i) for 25 cycles <5 % U _i (>95 % dip in U _i) for 5 sec	Mains power quality should be that of typical commercial or dental environment. If the user of the Bluephase 20i requires continued operation during power mains interruptions, it is recommended that the Bluephase 20i be powered from an uninterruptible power supply or battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or dental environment.

NOTE: U_i is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer’s declaration - electromagnetic immunity

The Bluephase 20i is intended for use in the electromagnetic environment specified below. The customer or the user of the Bluephase 20i should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should not be used closer to any part of the Bluephase 20i, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
			Recommended separation distance
Conducted RF IEC 61000-4-6	3 V _{ms} 150 kHz to 80 MHz	10 V	$d = 0.35 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	10 V/m	$d = 0.35 \sqrt{P}$ 80 MHz to 800 MHz $d = 0.70 \sqrt{P}$ 800 MHz to 2.5 GHz
			where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic

environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Bluephase 20i is used exceeds the applicable RF compliance level above, the Bluephase 20i should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Bluephase 20i.

^b Over the frequency range 150 kHz to 80 MHz, field strength should be less than 10 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the Bluephase 20i

The Bluephase 20i is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Bluephase 20i can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Bluephase 20i as recommended below, according to the maximum output power of the communication equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 0.4 \sqrt{P}$	80 MHz to 800 MHz $d = 0.4 \sqrt{P}$	800 MHz to 2.5 GHz $d = 0.7 \sqrt{P}$
0.01	0.04	0.04	0.07
0.1	0.13	0.13	0.22
1	0.40	0.40	0.70
10	1.3	1.3	2.2
100	4.0	4.0	7.0

Ivoclar Vivadent – worldwide

Ivoclar Vivadent AG

Bendererstrasse 2
FL-9494 Schaan
Liechtenstein
Tel. +423 235 35 35
Fax +423 235 33 60
www.ivoclarvivadent.com

Ivoclar Vivadent Pty. Ltd.

1 – 5 Overseas Drive
P.O. Box 367
Noble Park, Vic. 3174
Australia
Tel. +61 3 979 595 99
Fax +61 3 979 596 45
www.ivoclarvivadent.com.au

Ivoclar Vivadent Ltda.

Alameda Caiapós, 723
Centro Empresarial Tamboré
CEP 06460-110 Barueri – SP
Brazil
Tel. +55 11 2424 7400
Fax +55 11 3466 0840
www.ivoclarvivadent.com.br

Ivoclar Vivadent Inc.

1-6600 Dixie Road
Mississauga, Ontario
L5T 2Y2
Canada
Tel. +1 905 670 8499
Fax +1 905 670 3102
www.ivoclarvivadent.us

Ivoclar Vivadent (Shanghai)

Trading Co., Ltd.
2/F Building 1, 881 Wuding Road,
Jing An District
200040 Shanghai
China
Tel. +86 21 6032 1657
Fax +86 21 6176 0968
www.ivoclarvivadent.com

Ivoclar Vivadent Marketing Ltd.

Calle 134 No. 7-B-83, Of. 520
Bogotá
Colombia
Tel. +57 1 627 33 99
Fax +57 1 633 16 63
www.ivoclarvivadent.co

Ivoclar Vivadent SAS

B.P. 118
F-74410 Saint-Jorioz
France
Tel. +33 450 88 64 00
Fax +33 450 68 91 52
www.ivoclarvivadent.fr

Ivoclar Vivadent GmbH

Dr. Adolf-Schneider-Str. 2
D-73479 Ellwangen, Jagst
Germany
Tel. +49 (0) 79 61 / 8 89-0
Fax +49 (0) 79 61 / 63 26
www.ivoclarvivadent.de

Wieland Dental + Technik GmbH & Co. KG

Schwenninger Strasse 13
D-75179 Pforzheim
Germany
Tel. +49 (0) 72 31 / 37 05-0
Fax +49 (0) 72 31 / 35 79 59
www.wieland-dental.com

Ivoclar Vivadent Marketing (India) Pvt. Ltd.

503/504 Raheja Plaza
15 B Shah Industrial Estate
Veera Desai Road, Andheri (West)
Mumbai, 400 053
India
Tel. +91 (22) 2673 0302
Fax +91 (22) 2673 0301
www.ivoclarvivadent.in

Ivoclar Vivadent s.r.l.

Via Isonzo 67/69
40033 Casalecchio di Reno (BO)
Italy
Tel. +39 051 611 35 55
Fax +39 051 611 35 65
www.ivoclarvivadent.it

Ivoclar Vivadent K.K.

1-28-24-4F Hongo
Bunkyo-ku
Tokyo 113-0033
Japan
Tel. +81 3 6903 3535
Fax +81 3 5844 3657
www.ivoclarvivadent.jp

Ivoclar Vivadent Ltd.

12F W-Tower, 1303-37
Seocho-dong, Seocho-gu,
Seoul 137-855
Republic of Korea
Tel. +82 (2) 536 0714
Fax +82 (2) 596 0155
www.ivoclarvivadent.co.kr

Ivoclar Vivadent S.A. de C.V.

Av. Insurgentes Sur No. 863,
Piso 14, Col. Napoles
03810 México, D.F.
México
Tel. +52 (55) 50 62 10 00
Fax +52 (55) 50 62 10 29
www.ivoclarvivadent.com.mx

Ivoclar Vivadent Ltd.

12 Omega St, Rosedale
PO Box 303011 North Harbour
Auckland 0751
New Zealand
Tel. +64 9 914 99 99
Fax +64 9 914 99 90
www.ivoclarvivadent.co.nz

Ivoclar Vivadent Polska Sp. z o.o.

Al. Jana Pawla II 78
00-175 Warszawa
Poland
Tel. +48 22 635 54 96
Fax +48 22 635 54 69
www.ivoclarvivadent.pl

Ivoclar Vivadent Marketing Ltd.

Prospekt Andropova 18 korp. 6/
office 10-06
115432 Moscow
Russia
Tel. +7 499 418-03-00
Fax +7 499 418-03-10
www.ivoclarvivadent.ru

Ivoclar Vivadent Marketing Ltd.

Qlaya Main St.
Siricon Building No.14, 2nd Floor
Office No. 204
P.O. Box 300146
Riyadh 11372
Saudi Arabia
Tel. +966 1 293 83 45
Fax +966 1 293 83 44
www.ivoclarvivadent.com

Ivoclar Vivadent Pte. Ltd.

171 Chin Swee Road
#02-01 San Centre
Singapore 169877
Tel. +65 6535 6775
Fax +65 6535 4991
www.ivoclarvivadent.com

Ivoclar Vivadent S.L.U.

c/ Emilio Muñoz N° 15
Entrada c/ Albarracín
E-28037 Madrid
Spain
Tel. +34 91 375 78 20
Fax +34 91 375 78 38
www.ivoclarvivadent.es

Ivoclar Vivadent AB

Dalvägen 14
S-169 56 Solna
Sweden
Tel. +46 (0) 8 514 93 930
Fax +46 (0) 8 514 93 940
www.ivoclarvivadent.se

Ivoclar Vivadent Liaison Office

: Tesvikiye Mahallesi
Sakayik Sokak
Nisantasi' Plaza No:38/2
Kat:5 Daire:24
34021 Sisli – Istanbul
Turkey
Tel. +90 212 343 08 02
Fax +90 212 343 08 42
www.ivoclarvivadent.com

Ivoclar Vivadent Limited

Ground Floor Compass Building
Feldspar Close
Warrens Business Park
Enderby
Leicester LE19 4SE
United Kingdom
Tel. +44 116 284 78 80
Fax +44 116 284 78 81
www.ivoclarvivadent.co.uk

Ivoclar Vivadent, Inc.

175 Pineview Drive
Amherst, N.Y. 14228
USA
Tel. +1 800 533 6825
Fax +1 716 691 2285
www.ivoclarvivadent.us