Heliobond

English

Instructions for Use
Gebrauchsinformation
Mode d'emploi
Istruzioni d'uso
Instrucciones de uso
Instruções de Uso
Bruksanvisning
Brugsanvisning
Käyttöohjeet
Bruksanvisning
Productinformatie
Οδηγίες Χρήσεως

- Light-curing primer agent
- Lichthärtender Haftvermittler
- Agent de liaison photopolymérisable
- Promotore di adesione fotoindurente
- Adhesivo fotopolimerizable
- Adesivo fotopolimerizável
- Ljushärdande bonding resin
- Lyshærdende resin
- Valokovetteinen sidosresiini
- Lichtuithardend hechtmiddel
- Φωτοπολυμεριζόμενη υγρή ρητίνη

Instructions for Use

Description

Light-curing, single-component bonding material (bonding liquid) for optimizing the enamel etching technique in combination with all light-curing restorative materials (Heliomolar®, Helio Progress®, Tetric®, Tetric EvoCeram®)

Composition

Heliobond contains:

Bis-GMA	60 %wt.
Triethylene glycol dimethacrylate	40 %wt.

Indication

- Bonding agent (bonding liquid) for adhesive restorations
- Transparent sealant for fissures and restorations
- Bonding layer for the repair of resin-, crown-, or bridge works

Contraindication

Heliobond should not be used if patients are known to be allergic to any of the ingredients.

Other interactions

To achieve an optimum bond to the enamel, the enamel etching technique (e.g. with Total Etch) has to be applied. Contamination of the etched and dried enamel surface with saliva or blood hampers this effect. Phenolic substances (e.g. eugenol) inhibit polymerization. The use of base materials containing these substances is thus contra-indicated.

Application

Preparation of enamel

- Depending on the case, bevel enamel surfaces to be treated with a diamond or clean them carefully with an abrasive (e.g. pumice powder and water). Do not use any oily or greasy paste. Possible residue of base

CE 0123

594771/1105/WE3/H

ivoclar
vivadent
clinical

materials on the enamel are removed in that way. Then, rinse and dry.

- Apply Total Etch on the enamel surfaces to be treated and allow it to react. After that, carefully rinse with water spray and dry the etched surfaces with an air blower (use only water- and oil-free air).
- The area of application must be kept dry, e.g. with cotton rolls. For difficult working situations, we recommend isolation with a rubber dam.
- The etched and dried enamel surface must not be contaminated prior to the application of Heliobond (if the surface comes in contact with saliva or blood, the etching and drying procedures must be repeated).

Dosage

- Dosage depends on the indication. When used as a bonding agent (bonding liquid), a very thin layer of Heliobond has to be applied. If necessary, use an air blower to achieve an optimally thin layer.
- When used as a transparent fissure sealant, however, Heliobond must be applied in thicker layers.

Method of application

1. Bonding agent for adhesive restorations

- Etch desired enamel surfaces (see preparation of the enamel). If necessary, cover exposed dentin surfaces with Syntac® dentin bonding agent (see Instructions for Use of Syntac).
- Apply a thin layer of Heliobond with a brush, a spherical instrument, or by means of the attachable cannula on the etched enamel surface.
- You can achieve an optimally thin layer with an air blower.
- If used in combination with light-curing composites, Heliobond does not require separate polymerization. When working on freely modelled restorations, however, we recommend separate polymerization of Heliobond for 10 s (bluephase®).
- Application of the composite
- Polymerization of the composite
- Finishing of the restoration

2. Transparent sealant for fissures and pits

- Etch desired enamel surfaces (see preparation of the enamel).
- Apply Heliobond with a suitable instrument, a brush, or by means of the attachable cannula into the fissures. Avoid entrapment of air.
- Polymerize with light for 20 s (bluephase®).
- After curing, remove inhibited layer. Check occlusion and grind off any excess material.

3. Preparation of resin repair works

- Roughen the smooth resin surfaces.
- If necessary, apply Monobond S (see Instructions for Use of Monobond S).
- Apply a thin layer of Heliobond (use an air blower, if necessary).
- Polymerize with light for 10 s (bluephase®).

Polymerization

- 10 s polymerization time is sufficient for thin layers (bluephase®).
- Thicker layers (e.g. sealings) require 20 s curing time (bluephase®)
- The light probe should be held less than 5 mm above the Heliobond surface, but it should not touch unpolymerized material

Note

- Do not expose Heliobond to intensive light during application, as this significantly shortens the working time.
- If Heliobond is used as covering layer (e.g. as a sealant), a thin layer of unpolymerized, sticky material will remain on the surface after polymerization, due to oxygen inhibition. We recommend removing this layer with cotton rolls or pellets.

Storage

- Close bottle immediately after use. Exposure to light may lead to premature polymerization.
- Storage: 2–28 °C / 36–82 °F

- Do not use Heliobond after the date of expiration
- Date of expiration: See note on the bottle/packaging

Warning

Avoid contact of Heliobond with the skin/mucous membrane, or eyes. Unpolymerized Heliobond may cause slight, reversible irritation, or lead to a general sensitisation against methacrylates.

Keep material out of the reach of children!
For use in dentistry only!

Date information prepared: 10/2005

Manufacturer

Ivoclar Vivadent AG, FL-9494 Schaan / Liechtenstein

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent AB
Ivoclar Vivadent UK Limited
Ivoclar Vivadent, Inc.