# IPS InLine® – VENEER-VERARBEITUNG

## Konventionelle Metallkeramik

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
# IPS InLine® – VENEER FABRICATION

## Conventional metal-ceramic

The following chart shows the step-by-step layering of veneers on refractory dies.

**Important:** After each working step, the master model has to be immersed in water for about 5–10 minutes, depending on the size.


**TIP:** For the Veneer fabrication smaller working steps and several intermediate firing cycles are recommended.

### Fabricating the model

Fabricate a duplicating model and subsequently cast it with commercially available refractory die material, e.g. BegoForm® from Bego, Cosmotech VEST from GC or G-CERA™ VEST from GC (observe the instructions of the manufacturer).


**Important:** Correct processing and properly degassed dies are an important prerequisite for accurately fitting veneers.

### Wash firing

After degassing the refractory dies, apply IPS InLine Add-On mixed with the IPS InLine/IPS InLine PoM Glaze and Stains Liquid in a thin layer and fire.



T	B	S	t ↗	H	V <sub>1</sub>	V <sub>2</sub>
830°C 1526°F	403°C 757°F	4 min. 4 min.	60°C 108°F	1 min. 1 min.	450°C 842°F	829°C 1524°F

### Cervical firing

Build up the marginal areas using a mixture of Dentin and, for example, Occlusal Dentin brown.

T	B	S	t ↗	H	V <sub>1</sub>	V <sub>2</sub>
940°C 1724°F	403°C 757°F	8 min. 8 min.	60°C 108°F	1 min. 1 min.	450°C 842°F	939°C 1722°F

### Dentin/Impulse firing

Internal layering is modelled on the natural characteristics and consists of a dentin build-up and various effects.


Individual layering with the Impulse materials enables mamelons, opalescence and translucent effects to be achieved.

T	B	S	t ↗	H	V <sub>1</sub>	V <sub>2</sub>
940°C 1724°F	403°C 757°F	8 min. 8 min.	60°C 108°F	1 min. 1 min.	450°C 842°F	939°C 1722°F

### Incisal firing

Subsequently, build up the outer enamel layer and fire.

T	B	S	t ↗	H	V <sub>1</sub>	V <sub>2</sub>
930°C 1706°F	403°C 757°F	8 min. 8 min.	60°C 108°F	1 min. 1 min.	450°C 842°F	929°C 1704°F

### Glaze firing

Apply IPS InLine/IPS InLine PoM glaze paste to the surface and fire.

T	B	S	t ↗	H	V <sub>1</sub>	V <sub>2</sub>
860°C 1580°F	403°C 757°F	8 min. 8 min.	60°C 108°F	1 min. 1 min.	450°C 842°F	859°C 1578°F

### Divesting of laminate veneers

Remove large amounts of die material using a grinding disk. Subsequently, clean the veneers with Ivoclar Vivadent jet polishing material at approx. 1 bar pressure.

### Pretreatment of the veneer for adhesive luting

Etch the inner aspects with hydrofluoric acid (IPS Ceramic Etching Gel) for 120 seconds, thoroughly rinse with water and dry.

These firing temperatures are guidance values. They are valid for the Programat furnaces from Ivoclar Vivadent.

**ACCURATE, FAST and BEAUTIFUL – IPS InLine!**