# Digital Denture

Materials | Materialien

ProArt CAD Try-In | IvoBase® CAD | SR Vivodent® CAD, SR Vivodent® CAD Multi |  
ProArt CAD Transfer | IvoBase® CAD Bond | ProArt CAD Wax pink

**DE Gebrauchsinformation**

Seite 2

**EN Instructions for Use**

Page 20

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Page 38

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Pagina 56

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Página 74

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Página 92

**3** ProArt CAD Try-In

**5** IvoBase® CAD

**7** SR Vivodent® CAD, SR Vivodent® CAD Multi

**10** ProArt CAD Transfer

**11** IvoBase® CAD Bond

**18** ProArt CAD Wax pink

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**29** IvoBase® CAD Bond

**36** ProArt CAD Wax pink

Please read these Instructions for Use carefully and familiarize yourself with the use of these materials.

The following materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. These regulations also apply if the materials are used in conjunction with products of other manufacturers.

## Product Description

ProArt CAD Try-In are PMMA discs for the fabrication of try-in mock-ups, impression and functional trays using CAD/CAM technology.

Composition: Polymethyl methacrylate, pigments

## Indications

- Functional try-in mock-ups for removable denture prosthetics
- Functional try-in mock-ups for fixed denture prosthetics
- Impression and functional trays for removable denture prosthetics

## Contraindications

- Permanent dentures
- Temporary restorations
- Abutment mock-ups
- If the patient is known to be allergic to any of the components

## Side effects

Systemic side effects are not known to date. In individual cases local allergic reactions to methyl methacrylate materials have been reported.

## Processing with CAD/CAM

When designing the restoration, the specified minimum wall thicknesses and parameters of the final restoration must be complied with. If the minimum values are not observed, this can lead to failure (e.g. fracture).

ProArt CAD Try-In discs can be processed in CAD/CAM milling systems with a standard holder of 98.5 mm. Only use milling tools recommended for this milling system when processing the material. Please note that the material can be subject to overheating if milling tools, polishing brushes or steam cleaners are handled incorrectly. This may result in the material being damaged.

## Processing after the milling procedure

### Finishing

After milling, separate the completed object from the disc with a tungsten carbide bur. Use a cross-cut tungsten carbide bur for finishing or shape adjustments, and use PMMA polishing equipment if necessary.

The impression tray made of ProArt CAD Try-In can be used immediately for impression taking. High gloss polishing is not required due to the short time during which the material remains in the mouth.

### **Adjustment options**

If necessary, impression trays made of ProArt CAD Try-In can be adjusted using auto-polymerizing tray materials (e.g. SR Ivolen®). To do so, roughen or sandblast the surface. Then, clean it and adjust the material.

### **Impression materials**

If the fit of the impression tray is checked on the patient before impression taking, the tray must be cleaned under running water after the try-in. After cleaning and drying of the impression tray, it is ready for impression taking.

All impression materials which do not require moist storage can be used for impression taking with the respective adhesive.

### **Cleaning and disinfection**

Clean the impression tray together with the impression material under running water. Use a suitable commercial disinfecting agent for disinfection.

## **Important processing restrictions**

The following points have to be observed for the successful working with ProArt CAD Try-In:

- Observe the required minimum thickness
- Milling the discs using only a compatible CAD/CAM system

## **General instructions**

- Failure to observe the stipulated limitations of use and processing instructions may lead to failure.
- Do not inhale grinding dust during finishing. Use extraction equipment and wear a mask.
- No special storage conditions apply.
- Keep out of the reach of children!

## Product Description

IvoBase® CAD are impact-resistant PMMA discs for the CAD/CAM technology to fabricate denture bases for removable denture prosthetics using an subtractive manufacturing process.

Composition: Polymethyl methacrylate, co-polymer for impact toughness modification, pigments

### Typical material properties

Chem. description	PMMA-based denture base resin
Flexural strength	$\geq 65$ MPa
Flexural modulus	$\geq 2000$ MPa
Fracture toughness Kmax	$\geq 1.9$ MPam <sup>1/2</sup>
Overall work of fracture Wf	$\geq 900$ J/m <sup>2</sup>
Remaining quantity MMA	$\leq 4.5$ %
Water absorption	$\leq 32$ µg/mm <sup>3</sup>
Water solubility	$\leq 1.6$ µg/mm <sup>3</sup>

## Indication

Denture bases in removable denture prosthetics.

## Contraindication

The use of IvoBase CAD is contraindicated if the patient is known to be allergic to any of its ingredients.

## Side effects

Systemic side effects are not known to date. In individual cases local allergic reactions to methyl methacrylate materials have been reported.

## Processing with CAD/CAM

During construction, a minimum wall thickness (denture base) of 2 mm must be observed. Failure to observe the minimum thickness may result in failures (e.g. fracture of the restoration).

IvoBase CAD discs can be processed in CAD/CAM milling systems with a standard holder of 98.5 mm. Only use milling tools recommended for this milling system when processing the material. Please note that the material can be subject to overheating if milling tools, polishing brushes or steam cleaners are handled incorrectly. This may result in the material being damaged.

## Processing after the milling procedure

### Finishing

After milling, separate the restoration from the disc with a suitable tungsten carbide bur. Use a cross-cut tungsten carbide bur for finishing or shape adjustments. Pay special attention to the margins of the tooth cavities and the interdental spaces. The required minimum thickness must be observed.

### Polishing

Polishing is performed as usual using polishing instruments suitable for denture base materials.

### Repair / extension/ relining

Dentures fabricated with IvoBase CAD can be repaired, adjusted or relined using conventional auto-polymerizing denture base resin (e.g. ProBase® Cold).

The practical application corresponds with that of conventionally fabricated dentures. More details can be found in the respective Instructions for Use.

### Cleaning and disinfection

Clean IvoBase CAD under running water using a brush. Use a suitable commercial disinfecting agent for disinfection.

## Important processing restrictions

The following points have to be observed for the successful working with IvoBase CAD:

- Observe the required minimum thickness (2 mm)
- Milling the discs using only a compatible CAD/CAM system
- Staining/layering using materials that are approved and/or recommended

## General instructions

- Failure to observe the stipulated limitations of use and processing instructions may lead to failure.
- Contact of solvents or monomer with polymerized denture base material may lead to white discolouration.
- Do not inhale grinding dust during finishing. Use extraction equipment and wear a mask.
- Protect the materials from direct sunlight.
- Note the storage instructions and expiry date on the packaging.
- Keep out of the reach of children!

# SR Vivodent® CAD

# SR Vivodent® CAD Multi

CE 0123

## Product Description

SR Vivodent® CAD and SR Vivodent® CAD Multi are tooth-coloured discs for the CAD/CAM technology to fabricate permanent denture teeth for removable denture prosthetics using a subtractive manufacturing process.

Composition: Highly cross-linked DCL material

### Typical material properties

Chem. description	PMMA-based DCL material
Flexural strength	≥ 80 MPa
Water absorption	≤ 40 µg/mm <sup>3</sup>
Water solubility	≤ 7.5 µg/mm <sup>3</sup>

## Indication

Artificial individual teeth and tooth segments for the removable denture prosthetics.

## Contraindication

The use of SR Vivodent CAD (Multi) is contraindicated if the patient is known to be allergic to any of its ingredients.

## Side effects

Systemic side effects are not known to date. In individual cases local allergic reactions to methyl methacrylate materials have been reported.

## Processing with CAD/CAM

During construction, a minimum wall thickness of 2.5 mm must be observed. Failure to observe the minimum thickness may result in failures (e.g. fracture of the restoration) and compromise the shade effect.

In addition to fabricate individual teeth, it is also possible to fabricate tooth segments with multi units. When finishing segments manually, make sure that the connector dimension remains as large as possible.

SR Vivodent CAD (Multi) discs can be processed in CAD/CAM milling systems with a standard holder of 98.5 mm. Only use milling tools recommended for this milling system when processing the material. Please note that the material can be subject to overheating if milling tools, polishing brushes or steam cleaners are handled incorrectly. This may result in the material being damaged.

SR Vivodent CAD Multi must be placed in the milling machine so that the printed side is pointing toward the incisal/occlusal region of the restoration. The printed side corresponds with the incisal-occlusal side.

## Processing after the milling procedure

After milling in the CAD/CAM device, the restoration is separated from its holder with a fine tungsten carbide bur or a diamond separating disk and the attachment points are smoothed out using cross-cut tungsten carbide burs to prevent possible injury by sharp edges. Make sure that the anatomical surface structure is not too seriously damaged.

### Finishing

The following procedure is recommended for finishing SR Vivodent CAD (Multi) restorations:

- Carry out shape adjustments with cross-cut tungsten carbide burs (overheating of the material must be avoided)
- Check the fit of the teeth or tooth segments before cementation in the milled denture base
- In case of interfering contacts in the tooth cavity of the denture base, the tooth should be adjusted.
- Make sure that the minimum thickness is maintained even after the minor adjustments.

### Polishing

Careful polishing is the prerequisite for an optimum esthetic result. Polishing reduces plaque accumulation and the resulting shade disturbances. Interdental areas and occlusal surfaces should be given special attention.

Manual polishing of the restoration is carried out with rotary instruments and polishing paste.

In order to achieve an optimum surface gloss, please observe the following procedure:

- Polish the restoration at a suitable rpm using the handpiece and gentle pressure to avoid heat development
- Smooth out the surface (convex areas) of the natural structure, as well as the marginal ridges and the interdental areas with rubber polishers and silicone polishing wheels
- Pre-polishing is done with the handpiece / goat hair brushes and fine pumice / Universal Polishing Paste.
- The restorations are polished to a high gloss using a goat hair brush, cotton or leather buff, as well as Universal Polishing Paste.
- Modify the goat hair brush to become star-shaped in order to optimally polish the interdental area and occlusal surfaces. Given the smaller size brush, only the desired areas are polished.
- Depending on the type of high gloss desired, leather buffs can be used to achieve a high shine, while cotton buffs are used to achieve a lesser degree of lustre.

## Important processing restrictions

The following points have to be observed for the successful working with SR Vivodent CAD (Multi):

- Observe the required minimum thickness (2.5 mm)
- Milling the discs using only a compatible CAD/CAM system
- Staining/layering using materials that are approved and/or recommended

## General instructions

- When using SR Vivodent CAD Multi, always pay attention to the placement direction (incisal and dentin).
- SR Vivodent CAD (Multi) is based on the A–D shade system.
- With SR Vivodent CAD Multi the incisal area and therefore the general esthetics vary slightly according to the position of the restoration in the CAD/CAM software.
- For reasons of esthetics, the material is not suitable for the fabrication of crowns, bridges, hybrid abutments due to its high translucency.
- Avoid overheating of the material when finishing or when using an alcohol torch.
- Failure to observe the stipulated limitations of use and processing instructions may lead to failure.
- Contact of solvents or monomer may lead to white discolouration.
- Do not inhale grinding dust during finishing. Use extraction equipment and wear a mask.
- After removing the material from the packaging protect it from direct sunlight.
- Keep out of the reach of children!

# ProArt CAD Transfer

## Product Description

ProArt CAD Transfer are discs for the CAD/CAM technology to fabricate transfer templates for the cementation of artificial teeth such as pre-fabricated or milled teeth in a CAD/CAM fabricated denture base using a subtractive manufacturing process.

Composition: Polyoxymethylene (POM)

## Indication

- Transfer templates

## Contraindication

- The use of ProArt CAD Transfer is contraindicated if the patient is known to be allergic to any of its ingredients.
- Intra-oral application
- Do not use ProArt CAD Transfer for temporary or permanent restorations.

## Processing with CAD/CAM

During construction, the suggested values of the software must not be changed as this may result in failure (e.g. fracture or deformation of the transfer template).

ProArt CAD Transfer discs can be processed in CAD/CAM milling systems with a standard holder of 98.5 mm. Only use milling tools recommended for this milling system when processing the material. Please note that the material can be subject to overheating if milling tools, polishing brushes or steam cleaners are handled incorrectly. This may result in the material being damaged.

## Processing after the milling procedure

After milling, separate the completed object from the disc with a tungsten carbide bur and smooth out the attachment points to prevent possible injury by sharp edges.

## Important processing restrictions

The following points have to be observed for the successful working with ProArt CAD Transfer:

- Milling the discs using only a compatible CAD/CAM system

## General instructions

- Do not inhale grinding dust during finishing. Use extraction equipment and wear a mask.
- No special storage conditions apply.
- Keep out of the reach of children!

## Product Description


IvoBase® CAD Bond is a specially developed PMMA-based bonding material for the cementation of pre-fabricated or CAD/CAM-fabricated resin teeth and segments made of PMMA or composite in a denture base made of IvoBase CAD.

Composition

- **Powder**  
Polymethyl methacrylate, softener, benzoyl peroxide, catalyst, pigments
- **Liquid**  
Methyl methacrylate, dimethacrylate, catalyst
- **Modelling liquid**  
Methyl methacrylate, dimethacrylate, catalyst

List of parts:

- A) IvoBase CAD Bond Polymer (10 x)
- B) IvoBase CAD Bond Monomer (10 x)
- C) IvoBase CAD Bond Modelling Liquid (1 x)
- D) Application syringes (20 x)
- E) Cannulas (10 x)

## Indication

Cementation of artificial teeth such as pre-fabricated or CAD/CAM-fabricated teeth and segments into a denture base made of IvoBase CAD.

## Contraindication

The use of IvoBase CAD Bond is contraindicated if the patient is known to be allergic to any of its ingredients. Intraoral use of unpolymerized material. Cementation of ceramic teeth.

## Side effects

In individual cases local allergic reactions to methyl methacrylate materials have been reported.

## Safety notes

- Monomer and Modelling Liquid contain methyl methacrylate (MMA).
- MMA is highly flammable and irritating (flash point: +10 °C/50 °F).
- MMA is irritating to eyes, respiratory organs and skin.
- Skin contact may cause sensitization.
- Avoid skin contact with monomer and unpolymerized material. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.
- Do not inhale vapours.
- Keep away from sources of ignition - do not smoke.
- Prevent monomer from reaching the sewage system.
- Take measures against electrostatic charge.
- Do not inhale grinding dust. Use extraction equipment and wear a mask.

## General instructions

- Not suitable for metal or ceramic materials.
- With warmer ambient temperatures, it is recommended to store the material at a temperature of 2–8 °C / 36–46 °F (refrigerator).
- Closely observe the polymerization and processing times.
- The product may only be used and processed by trained personnel.
- Store in a humid environment between the completion of the dentures and their incorporation.
- The processing of IvoBase CAD Bond with pre-fabricated or CAD/CAM-fabricated teeth or segments is identical. The use of a transfer template is recommended.
- With oversized tooth segments no transfer template is required.

## Processing

### Step 1

Select the teeth used in the design software and completely remove wax residue from the tooth card.

#### Note:

Use a tooth sieve to store the teeth in order to avoid mix-ups.

#### Important!

Do not separate the milled denture base, since another milling procedure is to be performed in the basal area after the cementation of the teeth.

### Step 2

Separate the transfer template from the ProArt CAD Transfer. For this purpose, separate the holding pins with a suitable rotary instruments, e.g. resin bur, and smooth out the attachment points.

#### Note:

The use of the transfer template enables correct positioning of the teeth in the denture cavities.

### Step 3

Check the individual teeth in sequence for the correct fit in the transfer template (with viewing window). If necessary, remove interfering contacts with the help of occlusion spray from the transfer template rather than from the tooth.

The transfer template has been designed in such a way that the occlusal impressions of the maxillary and mandibular teeth are keyed in centric position. The locked bite is approximately 6 mm.

**Step 4**

Check the fit of the teeth in the denture base. In case of interfering contact of a tooth in the tooth cavity of the denture base, the tooth should be adjusted by grinding.

**Step 5**

Perform a final examination of the tooth position in the transfer template and cavity in the denture base.

**Step 6**

Pretreat the denture teeth by blasting them with  $100\ \mu\text{m}\ \text{Al}_2\text{O}_3$  at 1–2 bar. Remove blasting residue with oil-free compressed air. Ensure that the tooth surfaces are sufficiently roughened where they are to be adhesively fixed, circumferentially around the neck of the tooth as well as in the proximal areas.

**Step 7**

Pretreat the tooth cavities of the denture base by blasting them with  $100\ \mu\text{m}\ \text{Al}_2\text{O}_3$  at 1–2 bar.

### Step 8

Locally blast the tooth cavities in the denture base in the region of the teeth to approximately 1–2 mm beyond the contact surface of the marginal seam with 100 µm Al<sub>2</sub>O<sub>3</sub> at 1–2 bar. This creates a good bonding surface for Ivobase CAD Bond. Clean the pretreated surface with oil-free compressed air.

**Important! Do not clean with steam!**

### Step 9

Wet the teeth with monomer, e.g. ProBase® Cold Monomer. The monomer evaporates within 30 seconds. After that, contamination of the teeth must be prevented.

### Step 10

Wet the denture base in the region of the blasted bonding area with monomer and allow it to react.

### Step 11

The predosed Ivobase CAD Bond is sufficient for the cementation of a complete (max. 14 teeth) maxillary or mandibular denture. Add the entire monomer (B) to the polymer jar (A) and mix for 10–20 seconds with a clean spatula. Tapping results in air bubbles coming to the surface. Close the jar with the lid.

**Step 12**

Open the jar after another 10–20 seconds. IvoBase CAD Bond has a low-viscosity consistency now, which can be easily drawn up with two of the application syringes (D) provided. It is recommended to attach a cannula (E) to the syringe to be able to apply the material precisely in the smallest areas.

**Important!**

Close the jar immediately after filling the application syringe to prevent quick evaporation of the monomer.

**Step 13**

Now, apply the material to the tooth cavities starting with one side. Dosing depends on the fit of the denture teeth in the denture base.

**Step 14**

Once IvoBase CAD Bond has been applied, place the pretreated denture teeth in the cavities one after the other. The sequence can be freely chosen.

**Note:**

The processing time is approximately 10 minutes after the beginning of mixing (at 23 °C / 73 °F).

**Step 15**

Place the transfer template on the tooth row. The transfer template is used to verify and adjust the position of all denture teeth in the denture base and, ultimately, the intercuspatation.

**Step 16**

Perform a visual examination of the position of the teeth by means of the transfer template.

**Step 17**

Remove excess material with a brush or disposable sponge. A brush is suitable to contour the material while it is ductile. The additional use of Ivobase CAD Bond Modelling Liquid (C) facilitates contouring.

**Note:**

Avoid extensive wetting of the denture base with modelling liquid as this increases the risk of surface damage caused by MMA.

**Step 18**

If there are local deficits, more material can be applied and post-contoured and smoothed with a brush and modelling liquid. Thus, post-processing is reduced to a minimum.

**Step 19**

Secure the transfer template on the teeth to prevent the teeth from changing position. Polymerization should be conducted 8–12 minutes after the start of mixing in the pressure pot, e.g. Ivomat. Polymerization time: 15 minutes in a water bath at 50 °C (122 °F) and at 2–5 bar.

**Step 20**

Remove the transfer template after the polymerization time has elapsed. The disc can now be replaced in the holder of the CAM milling unit in the correct position with the help of notches milled during the previous milling procedure. After that, milling of the basal area is performed.

**Step 21, 22**

Separate the denture with the teeth from the disc at the holding pins with a suitable rotary instrument.

Finishing is carried out with cross-cut tungsten carbide burs and polishing media as it is usual for denture base materials.

**Storage Instructions**

- Store material in a cool, dark, and well ventilated place.
- Protect from direct sunlight.
- Storage temperature: 2–28 °C / 36–82 °F.
- Do not use the products after the indicated expiration date.
- Keep out of the reach of children!

## Product Description

ProArt CAD Wax pink are wax discs for the fabrication of prosthetic try-ins and individualized wax bite rims for the functional impression taking including wax bite registration using CAD/CAM technology.

### Properties

Chem. composition	Synthetic wax
Shade	Pink
Odour	Peppermint
Drop melting point	70–90 °C (158–194 °F)
Density	0.97 g/cm <sup>3</sup>
Water solubility	Insoluble

## Indication

- Try-ins
- Wax bite rims

## Contraindication

- The use of ProArt CAD Wax pink is contraindicated if the patient is known to be allergic to any of its ingredients.
- Not suitable for use as lost form in the casting technique.

## Side effects

Systemic side effects are not known to date.

## Processing with CAD/CAM

During construction, a minimum wall thickness (denture base) of 2.5 mm must be observed. Failure to observe the minimum thickness may result in failures (e.g. fracture).

ProArt CAD Wax pink discs can be processed in CAD/CAM milling systems with a standard holder of 98.5 mm. Only use milling tools recommended for this milling system when processing the material. The ambient temperature during processing should not exceed 35 °C (95 °F).

## Processing after the milling procedure

### Finishing

After milling, separate the completed object from the disc with a coarse-cut saw or a cross-cut tungsten carbide bur suitable for wax. Use a suitable cross-cut tungsten carbide bur for finishing or shape adjustments.

The material can be adapted with conventional set-up wax. Adhesive wax must not be used (risk of deformation).

**Use as wax try-in**

For prosthetic try-ins, the denture teeth are fixed in the cavities of the manufactured denture base using set-up wax. For this purpose, use a wax knife to make the hot set-up wax flow into the gap between the denture tooth and the cavity.

After the try-in on the patient, the wax-up can be transferred into a denture base resin using the conventional, shrinkage-compensated injection technique (e.g. IvoBase® System).

**Use in combination with impression materials**

If the fit of the impression tray is checked on the patient before impression taking, the tray must be cleaned under running water after the try-in. After cleaning and drying of the impression tray, it is ready for impression taking. Addition-reaction silicones are suitable for impression taking if used in combination with a respective tray adhesive.

**Cleaning and disinfection**

Clean ProArt CAD Wax pink after use on the patient under running water. Use a suitable commercial disinfecting agent for disinfection.

**Important processing restrictions**

The following points have to be observed for the successful working with ProArt CAD Wax pink:

- If the used pre-fabricated teeth penetrate the denture base on the basal side, a try-in with ProArt Wax pink is not recommended. Instead do a functional try-in with ProArt CAD Try-In.
- Observe the required minimum thickness (2.5 mm)
- Milling the discs using only a compatible CAD/CAM system

**General instructions**

- Failure to observe the stipulated limitations of use and processing instructions may lead to failure.
- Protect from direct sunlight.
- Store at 2–28 °C (36–82 °F).
- Keep out of the reach of children!

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© Ivoclar Vivadent AG, Schaan / Liechtenstein
**Ivoclar Vivadent Limited**
**Ivoclar Vivadent, Inc.**
Ivoclar Vivadent AG, 9494 Schaan/Liechtenstein