**Käyttöohjeet**  
**Bruksanvisning**  
**Productinformatie**  
**Οδηγίες Χρήσεως**  
**Kullanım Kılavuzu**  
**Инструкция**  
**Instrukcja stosowania**

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For dental use only.  
Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed dentist.

Made in Liechtenstein  
Ivoclar Vivadent AG,  
FL-9494 Schaan/Liechtenstein

**Instructions for Use**  
**Gebrauchsinformation**  
**Mode d'emploi**  
**Istruzioni d'uso**  
**Instrucciones de uso**  
**Instruções de Uso**  
**Bruksanvisning**  
**Brugsanvisning**

- Self-curing temporary c&b material
- Selbsthärtendes temporäres K&B Material
- Matériau autopolymérisant pour la réalisation au cabinet de couronnes, bridges et facettes provisoires
- Composito autoindurente per ponti e corone provisorio
- Material provisional C&B autopolimerizable
- Material autopolimerizável para coroas e pontes provisórias
- Kemisk hærdende provisorisk krone- & bromateriale
- Kemiallisokovetteinen kruunu-ja siltamateriaali

Complies with:  
ISO 4049, EN ISO 4049

**CE 0123**

**ivoclar**  
**vivadent**  
*clinical*

## Instructions for Use

### Description

Telio CS C&B is a self-curing composite material in paste form used for the fabrication of temporary crowns, bridges, inlays, onlays, core build-ups in combination with temporary posts and veneers. It is supplied in specially designed cartridges that provide a 10:1 base/catalyst ratio.

### Shades

Telio CS C&B is available in five shades: Bleach BL3, A1, A2, A3, A3.5. Telio CAD and Telio Lab are available in the same shades. They are coordinated and compatible with each other.

### Composition

Telio CS C&B is composed of polyfunctional methacrylates (48 wt %) and inorganic fillers (47 wt %). Additional ingredients include additives, initiators, stabilizers and pigments (5 wt %).

### Indications

Telio CS C&B is a temporary crown and bridge material suitable for the fabrication of crowns, bridges, inlays, onlays, core build-ups in combination with temporary posts and veneers. Because it is compatible with all the other products of the Telio range (both in terms of shade and chemistry), it is an ideal lining material for temporary restorations made of Telio CAD and/or Telio Lab.  
In addition, Telio CS C&B may be used for lining prefabricated polycarbonate resin crowns.

### Contraindication

Do not use if a patient is known to be allergic or intolerant to any of the ingredients of Telio CS C&B.

### Side effects

In rare cases, components of Telio CS C&B may lead to a sensitizing reaction. In this case, the use of the product should be discontinued.

### Application

#### 1. Preliminary impression

##### Fabrication of an alginate or silicone matrix.

If the unprepared teeth exhibit a suitable occlusal morphology that may be transferred to the provisional restoration or even the planned permanent restoration, and if abutment teeth in suitable position are available, a preliminary impression can be made to serve as a matrix for the fabrication of the temporaries. Prior to the preparation or extraction of the teeth, an impression of the pre-operative situation is made using addition-cured silicone impression material (e.g. Virtual®; dimensionally stable impressions) or alginate (e.g. Vival NF; impressions that are dimensionally stable for a limited period of time). To facilitate repositioning of the impression in the oral cavity, the interdental gingival areas should be trimmed and undercut

removed in the impression, if necessary. If the occlusion is defective, the fabrication of a wax-up in the laboratory may additionally be required prior to taking the impression.

**Plastic template.** When using a vacuum-formed plastic matrix care should be taken that the plastic material does not bond chemically to Telio CS C&B. Suitable plastic sheets (e.g. polyethylene sheets) should be selected, so that the vacuum-formed matrix can be reused.

If a bond between the plastic and Telio CS C&B is desired, plastic sheets should be chosen that are especially designed for this purpose. The interdental aspect of the matrix can be conditioned with AdheSE® Bond prior to applying the temporary restoration composite for adhesion.

## 2. Cartridge system

Fig. 1

1. Press the black release lever located below the plunger on the back of the dispenser and pull the plunger as far back as possible (Fig. 1).

Fig. 2

2. Lift the cartridge lock and insert the cartridge. Lower the cartridge lock (Fig. 2). Push slide forward until it touches the cartridge.

Fig. 3

3. Remove the sealing cap or used mixing tip by turning it a 1/4 turn counter clockwise and discard it, do not reuse it (Fig. 3).

Fig. 4

4. Insert a new mixing tip in a clean motion, so that the "V"-shaped notch on the colored base is aligned with the "V"-shaped notch on the housing of the cartridge. While gripping the colored base, and not the mixing tip, turn the mixing tip base 1/4 turn clockwise (Fig. 4).

5. The material is mixed and dispensed through the mixing tip by pulling the trigger and can be applied directly into the matrix.

## Removal of the cartridge

Lift the release lever below the plunger on the back of the dispenser and pull the plunger. Lift the cartridge lock on the front side on top of the dispenser and remove the cartridge.

**Leave the used mixing tip on the syringe as a cap until the next application.**

## 3. Application

Telio CS C&B is automatically mixed when it is dispensed with slight pressure directly into the impression or vacuum-formed matrix.

Fig. 5

Prior to each application, extrude and discard a pea-sized amount of material (Fig 5). To prevent voids, apply Telio CS C&B to the occlusal surfaces first, then proceed by filling the gingival areas. Keep the mixing tip immersed in the material to prevent bubbles.

If required, Telio CS C&B may also be syringed around the prepared teeth to avoid bubble formation along the preparation margin. The working time of the material is approx. 1 minute at room temperature (23 °C / 73 °F).

## 4. Shaping of temporaries

### a) Fabrication in the mouth

The prepared teeth should be slightly moist. This can be achieved by painting them with water-soluble glycerine gel (e.g. Liquid Strip). Undercut areas of the tooth preparation or adjacent teeth may have to be blocked out with appropriate materials (e.g. with wax). Carefully reposition the impression filled with Telio CS C&B on the prepared teeth. The setting time in the mouth (37 °C / 98 °F) is approximately 1 to 2 minutes. By then, Telio CS C&B has attained a hard-elastic consistency and can be removed from the mouth together with the impression.

### b) Fabrication on the model

Lubricate the preparation areas on the model with separator (e.g. Vaseline). Carefully reposition the impression filled with Telio CS C&B on the model. After approx. 3 minutes (at 23 °C / 73 °F), Telio CS C&B attains a hard-elastic, partially-set consistency and can be removed from the model together with the impression.

## 5. Post-curing and finishing

After removing the temporary restoration from the impression (or from the tooth preparation), remove excess material using rotary instruments once the material is cured completely (after 4 to 5 minutes). Suitable finishing instruments are fissured carbide burs. Remove the oxygen-inhibited layer with e.g. alcohol or by polishing the surface of the restoration. Polish using e.g. silicon carbide rubber polishes (e.g. Astropol®).

## 6. Cementation of the temporary

In order to reduce the sensitivity of the prepared dentin surfaces during the temporization phase, Telio CS Desensitizer may be applied prior to cementing the temporary restoration. Brush Telio CS Desensitizer into the dentin surfaces for 10 seconds using a suitable instrument (brush, application aid). Cautiously disperse excess to a thin layer with a stream of air. Subsequently seat the temporary restoration, preferably using a non-eugenol temporary cement (e.g. Telio CS Link)

## 7. Relining, repairs, addition of material

If Telio CS C&B temporaries are relined, repaired or adjusted, the following procedure is recommended:

Telio CS C&B temporaries can be repaired with freshly mixed Telio CS C&B material after having conditioned them with AdheSE Bond / Heliobond.

1. Grind the areas of the Telio CS C&B temporary which need to be relined, repaired or adjusted with a coarse diamond.

2. Apply AdheSE Bond / Heliobond.
3. Disperse AdheSE Bond / Heliobond with a very weak stream of air. Avoid pooling! Care should be taken not to blow away the bonding agent. All the prepared surfaces must be adequately coated with AdheSE Bond.  
**Note:** AdheSE Bond / Heliobond does not contain a solvent that needs to be evaporated!
4. Polymerize AdheSE Bond / Heliobond for 10 s using a LED or halogen curing light with a light intensity of more than 500 mW/cm<sup>2</sup> (eg bluephase in the LOP mode). If other curing lights are used (eg plasma arc or laser curing units), the recommendations of the respective manufacturer must be followed to obtain fully polymerized results. The curing time depends on the intensity and the wavelength range of the light emitted by the different curing units. Curing times of less than 5 s are not recommended.
5. Apply Telio CS C&B
6. If necessary, finish the relined, repaired and adjusted areas of the Telio CS C&B temporary with eg cross-cut tungsten carbide burs and subsequently polish them with silicon carbide rubber polishers (eg OptraPol®).

Alternatively, Telio CS C&B temporaries can be repaired or adjusted directly with a light-curing Ivoclar Vivadent composite (eg Telio Add-On Flow, Tetric EvoFlow®). In this case, the respective Instructions for Use or Telio Chair-side Instructions of Use must be observed.

If Telio CS C&B, Telio CAD and/or Telio Lab restorations are repaired, relined or adjusted at a later date, the respective Telio Chairside / Telio Labside Instructions for Use must be observed.

#### **Working time**

Approx. 1 minute at room temperature (23°C / 73°F)

#### **Setting time**

Approx. 1 to 2 minutes at 37 °C / 98 °F, approx. 3 minutes at room temperature (23 °C / 73 °F).

#### **Special notes**

- Uncured material can be removed using alcohol or similar solvents.

#### **Precautions**

- Prevent any contact of uncured Telio CS C&B with the skin/mucous membrane and eyes.
- Uncured Telio CS C&B may cause slight irritation and may lead to a sensitization against methacrylates.
- Commercially available medical gloves do not provide protection against the sensitizing effect of methacrylates.

#### **Storage and shelf life**

- Shelf life: see expiration date on label
- Storage temperature 2–28 °C / 36–82 °F
- Only use Telio CS C&B at room temperature (cooled material is more viscous and exhibits a slower setting behaviour).

**Keep out of the reach of children!  
For use in dentistry only!**

**Date information prepared:** 06/2010, Rev. 1

#### **Manufacturer**

Ivoclar Vivadent AG, FL-9494 Schaan/Liechtenstein

The product has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes.

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## Ivoclar Vivadent UK Limited
## Ivoclar Vivadent, Inc.