# Plaque Test

## Instructions for Use Gebrauchsinformation Mode d'emploi Istruzioni d'uso Instrucciones de uso Instruções de Uso Bruksanvisning Brugsanvisning

- Disclosing liquid to make bacterial biofilm visible
- Indikatorflüssigkeit zum Sichtbarmachen des bakteriellen Biofilms
- Liquide révélateur du biofilm bactérien
- Liquido indicatore per rendere visibile il biofilm batterico
- Líquido indicador de biofilm bacteriano
- Líquido evidenciador para mostrar un bio-filme bacterial, líquido e visível.
- Indikatorvätska för att göra bakteriell biofilm synlig

## Käyttöohjeet Bruksanvisning Productinformatie Οδηγίες Χρήσεως Kullanım Kılavuzu Инструкция Instrukcja stosowania

- Væske til synlig påvisning af bakteriel biofilm
- Indikaattoriiliuos bakteeribiofilmin paljastamiseen
- Indikatorvæske for synliggjøring av bakteriell biofilm
- Vloeistof om bacteriële biofilm zichtbaar te maken
- Υγρό αποκάλυψης που καθιστά ορατή τη βακτηριακή βιομεμβράνη
- Bakteriyel biyofilmî görünür hale getirmek için ortaya çıkarcı sıvı
- Жидкость для обнаружения бактериальной биопленки
- Płyn do wybarwienia płytki nazębnej

For dental use only.

Caution: U.S. Federal Law restricts this device to sale by or on the order of a licensed dentist.

Made in Liechtenstein  
Ivoclar Vivadent AG,  
FL-9494 Schaan/  
Liechtenstein

555092/WE3/0611/REV1

**ivoclar**  
**vivadent**  
clinical

## Step by Step

1. Apply Plaque Test generously on all tooth surfaces using a cotton swab.
2. Instruct the patient to rinse. Plaque Test adheres to areas covered with plaque.
3. Shine a polymerization light, e.g. bluephase on the teeth. Plaque covered with the disclosing liquid has a fluorescent yellow appearance.
4. Remove plaque by brushing the teeth or by professional tooth cleaning using Proxyl<sup>®</sup> pastes.

## Note

The glycerine gives Plaque Test a sweet taste. The disclosing liquid is sugarless.

## Storage

- Close bottle immediately after use.
- Shelf life: see label / packaging.
- Do not use Plaque Test after the indicated date of expiration.
- Storage temperature: 2–28 °C/36–82 °F.

## Store out of the reach of children!

## For use in dentistry only!

Date information prepared: 06/2011, Rev. 1

## Manufacturer

Ivoclar Vivadent AG, FL-9494 Schaan/Liechtenstein

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

## Deutsch

## English

## Instructions for Use

### Description

Plaque Test is a disclosing liquid used to reveal plaque on the surfaces of natural and restored teeth as well as implants. Under a polymerization light, e.g. bluephase<sup>®</sup>, plaque appears fluorescent yellow after Plaque Test has been applied.

### Composition

Fluorescein disodium salt < 1 wt.%, glycerine and distilled water.

### Indication

Exposure of plaque on tooth surfaces.

### Contraindication

If patients are known to be allergic to any of the ingredients in Plaque Test, the material should not be applied.

### Side effects

None known to date.

### Application

Use Plaque Test to selectively stain biofilm on teeth. This measure is particularly recommended to reveal retention niches for plaque during orthodontic treatment or if restorations or implants are present. The fluorescent pigment does not stain teeth, soft tissue and restorations. Only the plaque-coated areas glow under the light of a polymerization unit. This method can also be used in the determination of plaque indices.

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**Ivoclar Vivadent Marketing Ltd.**
**Ivoclar Vivadent AB**
**Ivoclar Vivadent Limited**
**Ivoclar Vivadent, Inc.**