# Variolink® II

## Try-In

### Instructions for Use

### Gebrauchsinformation

### Mode d'emploi

### Istruzioni d'uso

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### Bruksanvisning

### Brugsanvisning

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### Productinformatie

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For dental use only.  
Caution: U.S. Federal Law restricts this device to sale by or on the order of a licensed dentist.

Made in Liechtenstein  
Ivoclar Vivadent AG, FL-9494 Schaan/Liechtenstein

561838/0803/WE3/G

**ivoclar**  
**vivadent**  
clinical

English

## Instructions for Use

### Description

Variolink II Try-In pastes are water-soluble glycerine pastes. They permit shade simulation of indirect composite, and ceramic restorations (inlays, onlays, crowns, metal-free adhesive bridges, veneers) prior to final cementation using Variolink II luting composite. The ideal luting composite shade is determined with the help of the Try-In pastes. The translucency and shades of these pastes correspond to those of the thoroughly polymerized Variolink II luting composite (see Note).

Like Variolink II, the Try-In pastes are available in 6 shades and 3 degrees of translucency:

### Shades

	Chromascope	A–D
<b>Variolink II</b>		
<b>Base paste</b>		
bleach XL	010	–
transparent	T	–
white	110	A1
opaque white	–	–
yellow (universal)	210	A3
brown	340	A4
<b>Catalyst paste-high/low viscosity</b>		
transparent	T	–
yellow (universal)	210	A3
<b>Variolink II Try-In</b>		
bleach XL	010	–
transparent	T	–
white	110	A1
opaque white	–	–
yellow (universal)	210	A3
brown	340	A4

### Composition

Try-In pastes are composed of glycerine, mineral fillers, and colorants.

## Indication

Shade simulation of the final luting composite. Checking accuracy of fit of the final restoration. Protection of the etching pattern during try in when the ceramic has been etched by the laboratory.

## Contraindication

Do not use the Try-In pastes if the patient is known to be allergic to any of the material's ingredients.

## Side Effects

None known to date

## Interactions

Carefully remove Try-In paste from the restoration and the tooth before final seating. Otherwise, the bond between the restoration and Variolink II may be impaired.

## Application

Select the Try-In paste whose shade corresponds most accurately to that of the restoration. Then, apply a small amount of paste on the restoration and carefully seat it without pressure. If the restoration cannot be seated easily, check the contact surface as well as the proximal regions for premature contacts. Do not force the restoration into place (risk of fracture). Do not check the occlusion before the final cementation (risk of fracture).

## Important

Variolink II luting composite permits light-curing or dual-curing seating (light-curing and self-curing). When using the dual-curing method, simulation of the exact shade can be achieved by mixing Try-In paste shade 210/A3 (simulates Catalyst paste) with the second Try-In paste (simulates the selected Base paste). Pastes are mixed in equal parts.

Try-In pastes are water soluble. They are easily removed by holding the restoration under running water or placing it in an ultrasonic bath. Clean and dry the restoration thoroughly. Finally, the restoration is cemented using Variolink II (see Variolink II Instructions for Use).

## Note

Try-In paste shades and the corresponding Variolink II luting composite shades vary slightly in the uncured state. The shade of all luting composites changes slightly during and after polymerization. **Try-In paste shades are coordinated with those of the thoroughly polymerized luting composite.**

## Storage

- Do not use Variolink II Try-In after expiry date.
- Store the material at room temperature (2–28 °C/36–82 °F).
- Shelf life: see expiry date

**Keep out of the reach of children!  
For use in dentistry only!**

## Date information prepared

08/2003

## Manufacturer

Ivoclar Vivadent AG  
FL-9494 Schaan/Liechtenstein

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the materials for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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### **Ivoclar Vivadent S.A.**
### **Ivoclar Vivadent AB**
### **Ivoclar Vivadent, Inc.**