# AdheSE® / AdheSE® DC Activator

## Brugsanvisning

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## Brugsanvisning

## Productinformatie

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CE 0123

For dental use only.

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed dentist.

Made in Liechtenstein  
Ivoclar Vivadent AG,  
FL-9494 Schaan/Liechtenstein

ivoclar®  
vivadent®  
clinical

English

## Instructions for Use

### Description

**AdheSE** is a self-etching light-curing filled two-component adhesive system for enamel and dentin. It comprises

1. AdheSE Primer
2. AdheSE Bond

**AdheSE DC Activator** is an optional component, allowing the system to be used for dual/self-curing adhesive applications. The activator is mixed with AdheSE Bond in a 1:1 ratio (activation of the bonding component).

### Composition

**AdheSE Primer:** dimethacrylate, phosphonic acid acrylate, initiators and stabilizers in an aqueous solution

**AdheSE Bond:** HEMA, dimethacrylate, silicon dioxide, initiators and stabilizers

**AdheSE DC Activator:** initiators, solvents

### Indications of AdheSE (without DC Activator)

- Adhesive for direct light-curing composite and compomer restorations
- After light-curing of the bonding component for build-ups with dual or self-curing composites

### Indications of AdheSE with DC Activator

- Adhesive cementation of indirect metal-free restorations (inlays, onlays, crowns, bridges made of eg IPS Empress, IPS e.max CAD/Press) or of metal-free posts (eg FRC Postec Plus) with dual or self-curing composites (eg Variolink II)
- Direct restorations and core build-ups with dual and self-curing composites (eg MultiCore)

### Contraindication

If the patient is allergic to any of the components of AdheSE or AdheSE DC or if the stipulated working technique cannot be used.

Direct pulp capping.

## Side effects

Components of AdheSE may cause a sensitizing reaction in patients with a hypersensitivity to any of the material's ingredients. In these cases, the material should not be used.

## Interactions

Substances containing eugenol/clove oil inhibit the polymerization of AdheSE and AdheSE DC.

## Application

### A: Direct restorations – light-curing materials (composites, compomers, eg Tetric EvoCeram, Heliomolar, Compoglass)

If purely light-curing composites are used, AdheSE Bond must be light cured prior to the placement of a restorative. In this case, AdheSE is used without the DC Activator.

1. A dry operating field must be ensured. Ideally, a rubber dam should be used.
2. Apply an adequate amount of AdheSE Primer with a brush. All the cavity surfaces, beginning at the enamel, should be adequately wetted. If necessary, fresh primer should be applied. Once the cavity surfaces have been completely coated, the primer must be brushed into the entire surface for another 15 s. The total reaction time should not be shorter than 30 s.
3. Disperse excess AdheSE Primer with a strong stream of air until the mobile liquid film is no longer visible.
4. Apply AdheSE Bond, beginning at the dentin.
5. Disperse AdheSE Bond with a very weak stream of air. Avoid pooling! Care should be taken not to blow away the bonding agent. All the prepared surfaces must be adequately coated with AdheSE Bond.  
**Note:** AdheSE Bond does not contain a solvent that has to be evaporated!
6. Polymerize AdheSE Bond for 10 s using a LED or halogen curing light with a light intensity of more than 500 mW/cm<sup>2</sup> (eg bluephase in the LOP mode).

If other than halogen lights (eg plasma arc or laser polymerization units) are used, the recommendations of the manufacturer of the device used must be followed to obtain fully polymerized results. The light-curing time depends on the intensity and the wavelength range of the light emitted by the different curing units.

**Note:** Light-curing times of less than 5 s are not recommended.

7. The restorative material must be applied according to the instructions of the manufacturer.

### B: Indirect restorations / Posts:

1. Prepare the surface of the restorations according to the instructions of the manufacturer. The following conditioning steps apply to IPS Empress restorations after they have been tried in: Etch the inner surface of the restoration with IPS Ceramic Etching Gel (5% hydro-fluoric acid) for 60 s (IPS e.max CAD/Press for 20 s). Rinse thoroughly with water and dry with an air syringe. Silanize the inner surfaces with Monobond Plus for 60 s and dry with an air syringe.  
Do not touch the surfaces once they have been silanized!  
**Note:** Do not apply activated AdheSE Bond to the inner restoration surfaces!
2. A dry operating field must be ensured. Ideally, a rubber dam should be used.
3. Apply an adequate amount of AdheSE Primer with a brush. All the cavity surfaces, beginning at the enamel, should be adequately wetted. If necessary, fresh primer should be applied. Apply to all prepared surfaces for 30 s.
4. Disperse excess AdheSE Primer with a strong stream of air until the mobile liquid film is no longer visible. Remove excess from root canals with paper points.
5. Mix equal amounts of AdheSE Bond and AdheSE DC Activator. The two components must be mixed thoroughly.
6. Apply the activated AdheSE Bond to all prepared surfaces.
7. Remove excess material in such a way that the activated AdheSE Bond completely covers the enamel and dentin without pooling. This can be achieved by
  - a) applying a weak stream of air
  - b) holding a strong suction device over the cavity
  - c) drying root canals with paper points.

When indirect restorations or posts are placed in conjunction with dual-curing cement, activated AdheSE Bond must not be light cured, as a cured bond layer may compromise the fit of the restoration. If purely light-curing composites are used for cementation, the activated AdheSE Bond must be light cured prior to the placement of the restoration.

8. Next, coat the restoration with a light- / dual-curing cement (eg Variolink II) and place the restoration. Remove excess cement and light cure the restoration from all sides according to the instructions of the manufacturer, using a LED or halogen polymerization unit. Make sure that all the surfaces are fully cured. Check the light performance of the curing light on a regular basis.

**Note:** If cement is applied to cavity surfaces that have been coated with activated AdheSE Bond, premature polymerization may occur which could adversely affect the fit of the restoration.

**Note:** For the cementation of veneers, the total etch technique should be applied.

### **C: Direct restorations – self or dual-curing composites (eg core build-up composites)**

1. AdheSE Primer is applied, allowed to react and dispersed with air and the bonding agent is activated and applied as described in items B 3 to B7a.
2. Activated AdheSE Bond does not have to be light cured if self or dual-curing filling materials are used. If light-curing of the adhesive is preferred, AdheSE Bond may be polymerized for 10 s using a LED or halogen curing unit with a light intensity of more than 400 mW/cm<sup>2</sup> (eg bluephase in the LOP mode).
3. Mix the composite (eg MultiCore) according to its Instructions for Use and place it in the cavity.

### **Warning**

AdheSE is irritating. Avoid contact of the material with the skin, mucous membrane and eyes. If the product comes in contact with the skin, wash it off with plenty of water. In case of accidental contact with the eyes, wash immediately with plenty of water and seek medical attention. Commercial medical gloves do not protect against sensitivity towards methacrylate.

### **Storage:**

- See label / packaging for the expiry date
- Storage temperature 2–28 °C / 36–82 °F

**Keep out of the reach of children.  
For use in dentistry only.**

**Date of preparation:** 09/2010, Rev. 2

### **Manufacturer**

Ivoclar Vivadent AG  
9494 Schaan Liechtenstein

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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## Ivoclar Vivadent Limited
## Ivoclar Vivadent, Inc.