Brugsanvisning

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Bruksanvisning

Productinformatie

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CE 0123

For dental use only.
Caution: U.S. Federal Law restricts this device to sale by or on the order of a licensed dentist.

Made in Liechtenstein
Ivoclar Vivadent AG, FL-9494 Schaan /
Liechtenstein

ivoclar
vivadent®
clinical

633298/0410/W/E3/G

Instructions for Use

Product description

Telio Activator is a methyl methacrylate (MMA)-based liquid for extra-oral use, which is applied to activate the surfaces of restorations made of Telio CS C&B, Telio CAD and Telio Lab that are subsequently adjusted, relined or customized.

Composition

Methyl methacrylate (86-87 %), dimethacrylate and triethylene glycol dimethacrylate (13%) and catalyst (<1 %).

Indications

Extraoral activation of the surfaces of restorations made of

- Telio CS C&B
- Telio CAD
- Telio Lab

Contraindication

- Direct application in the oral cavity

Important processing notes

The following processing errors should be avoided in order to ensure successful results with Telio Activator:

- Insufficient surface conditioning and/or surface cleaning
- Reaction time too short
- Material applied too sparingly
- Combination with materials that are not approved and/or recommended

Side effects

In individual cases allergic reactions to MMA materials have been reported. In rare cases, components of Telio Activator may lead to sensitization. Telio Activator should not be used in such cases.

Application in the dental practice

- **Relining of Telio CAD / Telio Lab restorations with Telio CS C&B**
Roughen the area to be relined using coarse diamonds or sandblast (Al₂O₃, 100 µm, 1–2 bar/15–29 psi). Thoroughly rinse with water and dry with oil-free compressed air. Then

wet the areas to be bonded to extraorally with Telio Activator. Brush Telio Activator into the entire surface for at least 30 seconds using an application brush in order to ensure an even distribution and quick penetration of Telio Activator. Subsequently, allow the Activator to react for another 30–60 seconds (total reaction time 1-2 minutes). Now apply the Heliobond bonding agent, thinly disperse it with blown air and polymerize for 10 seconds (650 mW/cm², e.g. bluephase® using LOP).

Apply Telio CS C&B directly to the conditioned surface for the relining.

Please refer to the corresponding Instructions for Use (Telio Chairside or Telio Labside) for the detailed working steps regarding the relining with Telio CS C&B.

• Subsequent addition of material / repairs of Telio CS C&B restorations with Telio Add-On Flow

Roughen the area to be supplemented or repaired using coarse diamonds or sandblast (Al₂O₃, 100 µm, 1–2 bar/ 15–29 psi). Thoroughly rinse with water and dry with oil-free compressed air. Then wet the areas to be bonded to extraorally with Telio Activator. Brush Telio Activator into the entire surface for at least 30 seconds using an application brush in order to ensure an even distribution and quick penetration of Telio Activator. Subsequently, allow the Activator to react for another 30-60 seconds (total reaction time 1-2 minutes). Now apply the Heliobond bonding agent, thinly disperse it with blown air and polymerize for 10 seconds (650 mW/cm², e.g. bluephase using LOP). Apply Telio Add-On Flow or a light-curing Ivoclar Vivadent composite (e.g. Tetric EvoFlow, Tetric EvoCeram) and cure according to the Instructions for Use.

Application in the dental laboratory

• Customization of / addition of material to Telio CAD / Telio Lab restorations using

- Telio Lab
- Telio Stains
- Telio Lab LC Massen
- Telio Add-On Flow

Roughen the area to be supplemented by grinding or sandblast (Al₂O₃, 100 µm, 1–2 bar/15–29 psi). Then thoroughly clean with steam and dry with oil-free compressed air. After that, condition with Telio Activator. For that purpose, distribute Telio Activator evenly and generously on the conditioned surface and allow it to react for at least 2 minutes to a maximum of 4 minutes.

After the reaction time has passed, Telio Lab can be directly applied. If the other materials mentioned are used (Telio Stains, Telio Lab LC materials and Telio Add-On Flow), apply SR Composit to the activated surfaces according to the Instructions for Use in a first step and then polymerize. Subsequently, the Telio Stains, Telio Lab LC materials are layered to customize the restoration. Please refer to the corresponding Instructions for Use (Telio Lab/CAD) for the detailed working steps regarding the individualization or add-ons.

Overview of material combinations and processing

Restorative material	Surface Preparation	1) Surface Activation 2) Bonding Agent	Add-on material
Telio CS C&B	Blasting (Al ₂ O ₃ , 100 µm, 1–2 bar/ 15-29 psi pressure), or roughening with a coarse diamond	1) Telio Activator 2) Heliobond	Telio CS C&B Telio Add-On Flow Telio Stains
Telio Lab / Telio CAD	Blasting (Al ₂ O ₃ , 100 µm, 1–2 bar/ 15-29 psi pressure), or roughening with rotary instruments	1) Telio Activator 2) Heliobond	Telio CS C&B Telio Add-On Flow Telio Stains
		1) Telio Activator or 1) Telio Lab Cold Liquid	Telio Lab
Ivoclar Vivadent denture teeth made of PMMA and composite (excluding ceramic teeth)	Blasting (Al ₂ O ₃ , 100 µm, 1–2 bar/ 15-29 psi pressure), or roughening with rotary instruments	1) Telio Activator or 1) Telio Lab Cold Liquid 2) SR Composit	Telio Stains Telio Lab LC Telio Add-On Flow
		1) Telio Activator or 1) Telio Lab Cold Liquid	Telio Lab

Warning

- Telio Activator contains methyl methacrylate (MMA).
- MMA is highly flammable and irritating (flash point: +10 °C/50 °F).
- Do not inhale vapours.
- Irritating to eyes, respiratory organs and skin
- Skin contact may cause sensitization. Customary medical

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gloves do not provide protection against the sensitizing effect of methacrylates.

- Keep away from sources of ignition. Do not smoke.
- Prevent contamination of the sewage system.
- Take measures against electrostatic charge.
- Telio Activator is only suitable for extraoral application.

Storage

- Storage temperature 2-28 °C (36-82 °F).
- Shelf life: see expiry date on label/packaging.
- Do not use Telio Activator after the indicated expiry date.
- Close bottle immediately after use.

Keep out of the reach of children!

For use in dentistry only!

Date information prepared: 02/2010, Rev. 0

Manufacturer

Ivoclar Vivadent AG, FL-9494 Schaan/Liechtenstein

The product has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes.

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Ivoclar Vivadent UK Limited
Ivoclar Vivadent, Inc.