Multilink® Automix

Multilink® Primer A and B

Brugsanvisning
Käyttöohjeet
Brugsanvisning
Produktinformatie
Οδηγίες Χρήσεως
Kullanım bilgileri
Инструкция по применению
Instrukcja stosowania

Instructions for Use
Gebrauchsinformation
Mode d'emploi
Istruzioni d'uso
Instrucciones de Uso
Instruções de Uso
Bruksanvisning

- Selvhærdende dentalt kompositement med mulighed for lyspolymerisering
- Selvhærdende og selvetsende primer
- Itsekøvetten ja valinnaisesti valokosteitten yhdistelmämuokkuinnityksestä hammaslääketieteelliseen käyttöön
- Itsekøvetten ja itse-etsäava primer
- Kjemisk herdende dentalt kompositement med mulighet for lysherding
- Kjemisk herdende og selvetsende primer
- Zelf uit hardend tandheelkundig bevestigingscomposit met de optie van licht uit harding
- Zelf uit hardende en zelfetsende primer
- Αυτοπυρολιζόμενο οδοντιατρική συγκολλητική σύνθετη ρητίνη με προεραιτική δυνατότητα φωτοπολυμερισμού
- Αυτοπυρολιζόμενο και αυτοεπιπορευόμενο Primer
- Dig hekimiglindne kullannila amaçi, kendiliginde sertlegen, şikla sertlegme opsiyonu yaparirma kompoziti
- Kendiliginde sertlegen ve kendiliginde astileyen (self-etching) primer
- Самовердеющий композитный цемент с возможностью фотополимеризации
- Самовердеющий самоотравляющийся праймер
- Cement adhezivny na bazie materialow zlozonych o podwojnym mechanizmie polimeryzacji
- Samowytwarzający, chemoutwardzalny Primer

- Self-curing resin-based dental luting material with light-curing option
- Self-curing and self-etching primer
- Selbsthärtendes zahnrärztliches Befestigungscomposit mit optionaler Lichthärtung
- Selbsthärtender und selbstätzender Primer
- Composite de collage dentaire autopolymérisant, avec option photopolymérisation
- Adhészil automordáncant et autopolymérisant
- Cemento composito autoindurente con fotopolimerizzazione facultativa
- Primer automordenzante ed autoindurente
- Material de cimentação dental autopolimerizável em base a resina com opção de fotopolimerização
- Primer autopolimerizável y autograbante
- Composito de cimentação autopolimerizável, com opção fotopolimerizável.
- Primer auto-condicionante e autopolimerizável.
- Självhärdande resinbaserat dentalt cementeringsmaterial med ljushärdning som tillval
- Självhärdande och självettsande primer

CE 0123
ivoclar
vivadent
clinical

62.8305/081.010/WE

English

Instructions for Use

Description

Multilink® Automix is a self-curing luting composite with light-curing option for the adhesive luting of indirect restorations made of metal, metal-ceramic, all-ceramic and composite. Multilink Automix is applied together with the self-etching and self-curing Multilink Primer.

Monobond Plus is recommended as a coupling agent to achieve a stronger bond to precious and non-precious alloys, as well as to all-ceramics made of zirconium and aluminium oxide and silicate ceramics.

Shades

Multilink Automix is available in three shades with different degrees of translucency:

- Transparent (high translucency)
- Yellow (high translucency)
- Opaque (low translucency)

Working time

The working and setting times depend on the ambient temperature. Once Multilink Automix has been dispensed from the automix syringe, the following times apply:

	At room temperature 23 °C ± 1 °C (73 °F ± 2 °F)	Intraorally (in combination with Multilink Primer) approx. 37 °C ± 1 °C (99 °F ± 2 °F)
Working time	180 ± 30 sec	120 ± 30 sec
Curing time (without working time)	300 ± 30 sec	180 ± 30 sec

Mixing ratio

Multilink Automix is always dispensed from the automix syringe in the optimum ratio. Multilink Primer A and Multilink Primer B are mixed in a 1:1 ratio (e.g. 1 drop of Primer A + 1 drop Primer B) or a multiple thereof.

Composition

Multilink Automix

The monomer matrix is composed of dimethacrylate and HEMA. The inorganic fillers are barium glass, ytterbium trifluoride, spheroid mixed oxide. The particle size is 0.25–3.0µm. The mean particle size is 0.9µm. The total volume of inorganic fillers is approximately 40%.

Multilink Primer A and B

Multilink Primer A contains an aqueous solution of initiators. Multilink Primer B contains HEMA, phosphonic acid and acrylic acid monomers.

Indication

Multilink Automix and Multilink Primer are used for the permanent cementation of indirect restorations where a strong bond is desired:

- Inlays, onlays, crowns, adhesive bridges, bridges and endodontic posts made of
 - metal and metal-ceramics
 - all-ceramics, in particular opaque zirconium oxide ceramics
 - composites and fiber-reinforced composites

Contraindication

The use of Multilink Automix and Multilink Primer is contraindicated

- if a dry working field cannot be established or the stipulated working technique cannot be applied
- if a patient is known to be allergic to any of the ingredients of Multilink Automix or Multilink Primer A/B.

Side effects

Systemic side effects are not known to date. In individual cases, allergic reactions to single components have been reported.

Interactions

Phenolic substances (e.g. eugenol) inhibit polymerization. Consequently, application of materials containing these substances is to be avoided. Disinfectants with an oxidative effect (e.g. hydrogen peroxide) may interact with the initiator system, which in turn can impair the curing process. Thus, the automix syringe must not be disinfected using oxidative agents. The syringe can be disinfected by wiping it with e.g. medical alcohol. Alkaline jet media (e.g. Cojet) impair the effect of Multilink Primer.

Application

1. Removal of the temporary restoration and cleaning of the cavity

Remove possible residue of the temporary luting cement from the cavity or preparation with a polishing brush and an oil- and fluoride-free cleaning paste. Rinse with water spray. Subsequently, dry with water- and oil-free air. Avoid overdrying.

2. Try-in of the restoration and drying

Next, the shade, fit and occlusion of the restoration can be checked. Care should be taken when checking the occlusion of fragile and brittle ceramic objects before they are permanently cemented as there is a risk of fracture. If necessary, make corrections with fine diamonds at medium speed and with slight pressure using copious amounts of water. Polish ground surfaces. When applying an adhesive luting technique with composites, absolute isolation of the operative area – preferably with a rubber dam such as OptraDam, alternatively with cotton rolls and saliva ejector – is required. If the dental hard tissue is contaminated with saliva, it needs to be cleaned again as indicated in section 1.

3. Pre-treatment of the restoration

3.1.1 Glass-ceramic restorations (e.g. IPS Empress®)

- Etch with 5% hydrofluoric acid (e.g. IPS® Ceramic etching gel) for 60 seconds or according to the instructions of the manufacturer of the restorative material.
- Thoroughly rinse the restoration with water spray and dry with oil-free air.

3.1.2 Lithium disilicate glass-ceramic restorations

(e.g. IPS e.max® Press, IPS e.max CAD)

- Etch with 5% hydrofluoric acid (e.g. IPS® Ceramic etching gel) for 20 seconds or according to the instructions of the manufacturer of the restorative material.
- Thoroughly rinse the restoration with water spray and dry with oil-free air.

3.1.3 Zirconium oxide (e.g. IPS e.max ZirCAD) and aluminium oxide ceramic restorations

- Sandblast the inner restoration surfaces (sandblasting parameters according to the instructions of the manufacturer of the restorative material).
- If necessary, clean the restoration in an ultrasonic unit for about 1 minute.
- Thoroughly rinse with water spray and dry with oil-free air.
- **IMPORTANT!** In order to create a strong bond, do not clean the zirconium oxide surfaces with phosphoric acid.

3.1.4 Metal or metal-supported restorations

- Sandblast the inner restoration surfaces (sandblasting parameters according to the instructions of the manufacturer of the restorative material) until a uniformly mat surface is achieved.
- If necessary, clean the restoration in an ultrasonic unit for about 1 minute.
- Thoroughly rinse with water spray and dry with oil-free air
- **IMPORTANT!** In order to create a strong bond, do not clean the metal surfaces with phosphoric acid.

3.1.5 Composite or fibre-reinforced composite restorations

- Sandblast the inner restoration surfaces (sandblasting parameters according to the instructions of the manufacturer of the restorative material).
- If necessary, clean the restoration in an ultrasonic unit for about 1 minute.
- Thoroughly rinse with water spray and dry with oil-free air.

3.2 Apply Monobond Plus with a brush or microbrush to the pre-treated surfaces.

Allow the material to react for 60 seconds. Disperse with a strong stream of air.

Special notes

- Contamination with saliva, blood or water must absolutely be avoided during the entire procedure.
- In case of contamination, the entire treatment of the restorations is to be repeated (see section 3).

4. Mixing of Multilink Primer A and Multilink Primer B

Mix the two Multilink Primer liquids A and B in a 1:1 mixing ratio (e.g. 1 drop of Primer A and 1 drop of Primer B). The mixed Primer A/B is solely self-curing and does not need to be protected against light, but it must be applied within 10 minutes.

5. Application of the mixed Multilink Primers A/B to the enamel and the dentin

Apply the mixed Multilink Primer A/B to all preparation surfaces (prepared enamel and dentin) using a microbrush – starting from the enamel if there is any enamel left and scrubbing with slight pressure for 15 seconds. If necessary, fresh primer should be applied. A reaction time of 30 seconds is recommended on the enamel and 15 seconds on the dentin. Disperse excess Multilink Primer with a strong stream of air until the mobile liquid film is no longer visible. As the Primer is solely self-curing, no light-curing is necessary.

6. Application of Multilink Automix to the restoration

For each application place a new automix tip on the syringe. Dispense Multilink from the automix syringe and apply the desired quantity directly onto the restoration. As the luting material will cure in the partially used mixing tip, it may serve as a seal for the contents of the syringe until needed once again (replace with a new tip just before the next use).

Note

Multilink Automix should be used quickly after it has been dispensed. As soon as Multilink Automix comes in contact with the mixed Multilink Primer A/B, the self-curing reaction is accelerated so that the curing time is shortened.

Apply Multilink Automix directly to the inner surface of the restoration. The application of Multilink directly to the preparation or the cavity, which has been pre-treated with Multilink Primer, is not indicated as this would lead to a considerable acceleration of the curing process and thus may cause improper fit.

7. Placement of the restoration and removal of excess cement

a) solely self-curing

Seat the restoration in place and fix/hold. Remove excess material immediately with a microbrush/brush/foam pellet/dental floss or scaler. Be sure to remove excess material in time, especially in areas that are difficult to reach (proximal or gingival margins). Due to the reaction between Multilink Automix and Multilink Primer A/B, a high bond strength and a high degree of cure is achieved only a few minutes after the restoration has been seated in place.

b) self-curing with additional light-curing

(quarter technique, indicated for cases with up to 2 bridge abutments = 3- to 4-unit bridges)

Seat the restoration in place and fix/hold. Excess material is light-cured with the polymerization light (approx. 650 mW/cm², e.g. bluephase, LOP mode) for **2-4 seconds per quarter surface** (mesio-oral, disto-oral, mesio-buccal, disto-buccal) at a distance of approx. 0-10 mm. Excess material can be easily removed with a scaler. Be sure to remove excess in time, especially in areas that are difficult to reach (proximal or gingival margins). Subsequently, light-cure all margins again for 20 seconds (approx. 1200 mW/cm², e.g. bluephase, HIP mode).

c) self-curing with additional light-curing

(quarter technique, indicated for cases with up to 6 bridge abutments = circular bridges)

Light-curing and subsequent removal of excess:

Seat the restoration in place and fix/hold. Excess material on all abutments is light-cured with the polymerization light (approx. 650 mW/cm², e.g. bluephase, LOP mode) for **1 second per quarter surface** (mesio-oral, disto-oral, mesio-buccal, disto-buccal) at a distance of approx. 0-10 mm. Excess material can be easily removed with a scaler. Be sure to remove excess in time, especially in areas that are difficult to reach (proximal or gingival margins). Subsequently, light-cure all margins again for 20 seconds (approx. 1200 mW/cm², e.g. bluephase, HIP mode).

Note

As all composites, Multilink Automix is subject to oxygen inhibition. This means that the surface layer (approximately 100 µm) does not polymerize during curing, as it comes in contact with atmospheric oxygen. In order to avoid this, we recommend covering the restoration margins with glycerine gel/air block (e.g. Liquid Strip) immediately after the removal of excess. After complete polymerization, the glycerine gel/air block is rinsed off with water.

8. Finishing of the restoration

- Apply finishing and polishing strips in the proximal regions
- Check occlusion and function and make corrections if necessary.
- Polish restoration with polishers (Astropol®) or discs.

Special instructions for the cementation of endodontic posts

1. For the cementation of endodontic posts, thoroughly clean the root canal (Residue of eugenol-based sealers in particular may inhibit the polymerization of the luting composite). Ideally, the mixed Multilink Primer A/B should be applied to the root canal and occlusal surface of the preparation for about 15 seconds using the thin (violet) microbrushes. Remove excess material from the root canal using paper points.
2. Coat the endodontic post, which has been prepared according to the instructions of the post manufacturer, with the mixed Multilink cement. **IMPORTANT!** Do not lute the Multilink cement into the root canal that has been wetted with Multilink Primer A/B. Otherwise, premature curing might occur which possibly prevents the post from being placed into its required position.

3. The endodontic post is placed, thus producing excess cement.
4. Excess Multilink cement can be distributed across the occlusal preparation surface for total coverage. Then light-cure Multilink Automix for 20 seconds, maintaining the post in place with the polymerization light.
5. Core build-up material (e.g. MultiCore®) can then be directly applied onto Multilink Automix, which serves as bonding agent, and cured according to the instructions of the manufacturer (light-curing).

Warning

Unpolymerized Multilink Automix pastes and the Multilink Primer liquids may cause slight irritations. Avoid contact with the skin, mucous membrane or eyes. If Multilink comes into contact with the eyes, immediately rinse with water and, if required, seek medical advice. If the material comes into contact with the skin, rinse with copious amounts of water. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Storage

- Do not use Multilink Automix and Multilink Primer after the expiry date.
- Cool storage (2–8 °C/ 36–46 °F) for Multilink Automix and Multilink Primer.
- Multilink Primer bottles should be tightly closed after use.
- In order to ensure a tight seal of the Multilink Automix syringe, leave the used mixing tip on the syringe.
- Shelf life: see expiration date

Keep material out of the reach of children!

For use in dentistry only!

Date information prepared: 09/2010; Rev. 1

Manufacturer:

Ivoclar Vivadent AG
FL-9494 Schaan/Liechtenstein

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.