**Bruksanvisning**  
**Brugsanvisning**  
**Käyttöohjeet**  
**Bruksanvisning**  
**Productinformatie**  
**Οδηγίες Χρήσεως**  
**Kullanım Kılavuzu**

**Instructions for Use**  
**Gebrauchsinformation**  
**Mode d'emploi**  
**Istruzioni d'uso**  
**Instrucciones de uso**  
**Instruções de Uso**

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## Instructions for Use

### Description

Variolink Veneer is a micro-filled, light-curing luting composite system for the adhesive luting of ceramic and composite restorations with low thickness (up to 2.0 mm), e.g. veneers, inlays, onlays.

Variolink Veneer exhibits the following characteristics:

- Optimized consistency and easy removal of excess
- Very high colour stability
- High translucency
- Excellent adhesive strength
- High resistance to wear

The following Variolink Veneer shades are radiopaque (>100% relative AI standard):

- Low Value -1
- Low Value -2
- Low Value -3
- High Value +3

As a result of their high translucency, the following Variolink Veneer shades are not radiopaque (<100% relative AI standard):

- Medium Value 0
- High Value +1
- High Value +2

### Shades

The gradation of Variolink Veneer shades is based on the effect a luting composite exerts on the brightness value of the final restoration. Variolink Veneer Medium Value 0 affects the brightness value only slightly or not at all. At the same time, it demonstrates the highest translucency and is thus the most neutral in shade. While the High Value pastes +1, +2, +3 let the restoration appear brighter, Low Value -1, -2, -3 impart a "warmer" or darker character to the restoration. In order to evaluate the overall effect of the various shades prior to permanent cementation, the use of the Variolink Veneer Try-In pastes is recommended.

The light-curing Variolink Veneer shades are available in the following shade gradations:

- |                  |                           |                           |
|------------------|---------------------------|---------------------------|
| - High Value +3  | Translucency approx. 5 %  | ↑ Brighter / more whitish |
| - High Value +2  | Translucency approx. 12 % |                           |
| - High Value +1  | Translucency approx. 24 % |                           |
| - Medium Value 0 | Translucency approx. 50 % |                           |
| - Low Value -1   | Translucency approx. 14 % | ↓ Darker / more yellowish |
| - Low Value -2   | Translucency approx. 9 %  |                           |
| - Low Value -3   | Translucency approx. 9 %  |                           |

CE 0123

5972260210WEE3H

ivoclar  
vivadent®  
clinical

## Working time

The Variolink Veneer pastes are to be used exclusively in light-curing procedures.

## Composition

Variolink Veneer is composed of dimethacrylates, silicon dioxide and ytterbium trifluoride. Additional contents are catalysts, stabilizers and pigments.

- The total content of inorganic filler is approx. 40 vol%.
- Filler particle sizes range from 40 nm to 300 nm.

## Indication

Adhesive luting of ceramic and composite restorations (fabricated at the chair-side or in the laboratory) with a low layer thickness (< 2.0 mm) that enable the use of a purely light-curing technique due to their high translucency:

- Preferably veneers, inlays and onlays made of e.g. IPS Empress®, IPS Empress® 2, IPS e.max® Ceram

## Contraindication

The application of Variolink Veneer is contraindicated:

- If a dry working field cannot be established, or the stipulated working procedures cannot be applied.
- If the patient is known to be allergic to any of the ingredients of Variolink Veneer.
- For opaque restorations made of e.g. metal

## Side effects

Systemic side effects are not known to date. In rare cases, allergic reactions to individual components have been reported. The use of the product should not be continued in these cases.

## Interactions

Phenolic substances (e.g. eugenol) inhibit polymerization. Consequently, materials containing this type of ingredient (e.g. eugenol-containing provisional luting materials) should not be used in conjunction with Variolink Veneer.

## Application (step by step)

Please also refer to the Instruction for Use of the products mentioned for more detailed information.

### Try-in and adhesive luting of tooth-coloured restorations

#### 1. Removal of the temporary restoration

After removing the provisional restoration, clean cavity or the prepared tooth using a polishing brush and non-fluoride, oil-free paste. Make sure to remove all temporary cement residues. Following this, rinse with water spray.

## Note:

Use only eugenol-free temporary luting materials, as eugenol may inhibit the polymerization of Variolink Veneer and other luting composites.

#### 2. Try-in and shade adjustment

For optimum aesthetic results, evaluate the effect of Variolink Veneer on the shade of the restoration with the help of the Variolink Veneer Try-In pastes at this stage. It is not advisable to verify occlusion at this point, as this might lead to fracture. Subsequently, adjust proximal contacts and re-polish using ceramic or composite polishing instruments.

#### 3. Isolation

In conjunction with the adhesive luting technique, the use of rubber dam isolation is recommended.

#### 4. Cleansing of the restoration

Use water spray to wash the Variolink Veneer Try-In paste off the restoration and dry it with moisture- and oil-free air. Clean the surface of the restoration by etching it with phosphoric acid for 15 seconds. Rinse off and dry.

#### 5. Pre-treatment of the restoration

- Make sure that the restoration was pre-treated in the laboratory according to the instructions of the manufacturer (e.g. check whether the ceramic material was etched with hydrofluoric acid).
- Apply silane to the inner surfaces of the restoration (e.g. Monobond Plus) and leave to react for 60 seconds. Subsequently, dry the surface with oil-/moisture-free air. Do not rinse.

#### 6. Pre-treatment of the cavity

- Apply Total Etch (37% phosphoric acid) to the prepared enamel first, and then to the dentin surfaces. Allow phosphoric acid to react for 15–30 seconds on the enamel, and a maximum 15 seconds on the dentin.  
Using a matrix band in the proximal regions prevents accidental etching of the adjacent teeth and facilitates the subsequent removal of excess Variolink Veneer. Secure the matrix bands in place using suitable interdental wedges.
- Then thoroughly rinse with copious amounts of water for at least 5 seconds. Remove excess water, leaving the dentin surface slightly, but visibly, moist for wet bonding. This can be done with either:
  - a) compressed air, or
  - b) a high-volume evacuation tip held directly over the prepared surface for 1–2 seconds, or
  - c) a dry brush, foam pellet or other lint-free absorbent material

**Do not overdry dentin!**

## 7. Application of the bonding agent (see Instructions for Use of the bonding agent)

## 8. Variolink Veneer luting techniques

Variolink Veneer is especially recommended for highly translucent, aesthetic restorations in the anterior region (e.g. veneers). After selecting a suitable shade, apply Variolink Veneer paste directly to the inner surfaces of the restoration. For controlled, precise placement use a Variolink Veneer application tip.

### Note:

Variolink Veneer is a light-curing material and thus sensitive to blue light (operatory light, ambient light). Therefore, Variolink Veneer should not be dispensed from the syringe until right before use. Exposure to intensive light should be avoided during application.

## 9. Placement

### a) Inlays/onlays

Use a brush or spatula to apply Variolink Veneer to the cavity and/or the inner surface of the restoration (to avoid the inclusion of air bubbles if concavities are present).

### b) Veneers

Apply Variolink Veneer directly to the inner surface of the restoration using a brush or spatula.

## 10. Removal of excess

### a) Conventional technique

- Seat the restoration exerting light pressure and remove gross excess with a suitable instrument (e.g. spatula, brush).
  - Make sure to remove excess material in time, particularly in areas that are difficult to reach (proximal or gingival margins).
  - Increase pressure and maintain for several seconds.
  - Maintain the pressure and light-cure a small area with a small light probe for a few seconds to tack the restoration into place.
- Important: Do not light-cure any proximal or marginal regions!**
- Remove further excess with a suitable instrument. When doing this, take care not to wipe any material out of the marginal areas.
  - Thoroughly light-cure each section for 10-30 seconds, depending on the light intensity of the curing unit used.

### b) "Complete curing technique"

- Seat the restoration exerting light pressure.
- Maintain the pressure and light-cure a small area with a small light probe for 10 seconds to tack the restoration into place
- Thoroughly light-cure each section for 10–30 seconds.
- Variolink Veneer excess may be removed using suitable instruments (e.g. chipping it away with a probe or scaler).

### c) "Wave technique"

- Seat the restoration exerting light pressure.
- Maintain pressure on the restoration and wave the light guide of a curing light from a distance of 2–3 cm across the gingival areas of the restorations for 2–3 seconds, initiating the cement to a "semi-set" stage. Use an explorer to touch the cement and verify the "semi-set" stage which demonstrates a solid, rubbery consistency. Additional cure time may be required to achieve the desired consistency.
- Larger pieces of excess Variolink Veneer can be easily removed from the margins and interproximally. Remove all excess cement.
- Thoroughly light-cure from all aspects for 10–30 seconds.
- **Note:** The Wave Technique is time, intensity, and distance dependant. It is critical to find the combination of variables that work using your own experience.

## Notes on oxygen inhibition

As all composite resins, Variolink Veneer is subject to oxygen inhibition. This means that the surface layer (approx. 50 µm in height), which is exposed to atmospheric oxygen during the curing process, does not polymerize. There are two different ways to overcome this problem:

- 1) Leave excess Variolink Veneer in place. Thus, only the excess is affected by oxygen inhibition and can easily be removed after polymerization when finishing of the margin is performed.
- 2) Cover restoration margins with Liquid-Strip (glycerine gel) after the removal of excess, but prior to polymerization. In this way, oxygen inhibition is prevented.

## Notes on the curing unit

The curing time per section is strongly dependent on the light intensity achieved by the curing unit. As many different types of curing lights are currently available (quartz halogen lights, laser lights, LED lights) we ask you to refer to the manufacturer's operating manual for the appropriate curing times.

## 11. Finishing and polishing

- Remove any excess material using finishing diamonds (grit-size below 25 µm) and flexible polishing discs.
- Apply finishing and polishing strips in the proximal regions.
- Verify occlusion and adjust if required. Any adjustments can be made with finishing diamonds operated at a speed of 3000–7000 rpm, exerting light pressure. Ground areas need to be repolished.
- Polish restoration margins with silicone polishers (Politip®-F-P, Astropol®) or polishing discs.

**Warning**

Avoid contact of uncured Variolink Veneer with the skin, mucous membranes and eyes. Uncured Variolink Veneer may have a slight irritating effect and cause sensitization to methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates!

**Storage and shelf life**

- Do not use Variolink Veneer after the indicated date of expiration.
- Variolink Veneer can be stored at room temperature.
- Important: Close syringes immediately after use. Exposure to light leads to premature polymerization.
- Shelf life: See expiration date

**Keep material out of the reach of children!****For use in dentistry only!**

**Date information prepared:** 09/2009, Rev. 1

**Manufacturer:**

Ivoclar Vivadent AG, FL-9494 Schaan / Liechtenstein

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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## **Ivoclar Vivadent Limited**
## **Ivoclar Vivadent, Inc.**