


FRONTZAHNKRONE | ANTERIOR CROWN

Variante: Standard Charakterisierung | Variant: Standard characterization




Ausgangssituation
Starting situation



1. Schritt / 1st Step:
Dentin: Shade A2

2. Schritt / 2nd Step:
Incisal: Shade I1

3. Schritt und Brennen / 3rd Step and firing:
Halo: Shade A2




Endsituation glasiert
Completed restoration glazed

Variante: Individuelle Charakterisierung | Variant: Individual characterization




Ausgangssituation
Starting situation



1. Schritt / 1st Step:
Cervical/Incisal: Stains orange
Dentin: Shade A2

2. Schritt / 2nd Step:
Incisal: Shade I1

3. Schritt und Brennen / 3rd Step and firing:
Halo: Shade A2
Approximal / proximal: Stains white



Endsituation glasiert
Completed restoration glazed

Brennparameter Malfarben- und Charakterisierungsbrand
Firing parameters Stains and characterization firing

B	S	t _↑	T	H	V ₁	V ₂
403°C	6 min	100°C	790°C	1 min	450°C	789°C
757°F	6 min	180°F	1454°F	1 min	842°F	1453°F


Brennparameter Glanzbrand
Firing parameters Glaze firing

B	S	t _↑	T	H	V ₁	V ₂
403°C	6 min	100°C	790°C	1–2 min	450°C	789°C
757°F	6 min	180°F	1454°F	1–2 min	842°F	1453°F


Bitte entsprechende Verarbeitungsanleitung beachten. / Please read the corresponding Instructions for Use.

SEITENZAHNKRONE | POSTERIOR CROWN

Variante: Standard Charakterisierung | Variant: Standard characterization




Ausgangssituation
Starting situation



1. Schritt / 1st Step:
Occlusal: Shade 130

2. Schritt / 2nd Step:
Incisal: Shade I1

3. Schritt und Brennen / 3rd Step and firing:
Fissur / fissure: Stains orange



Endsituation glasiert
Completed restoration glazed

Variante: Individuelle Charakterisierung | Variant: Individual characterization




Ausgangssituation
Starting situation



1. Schritt / 1st Step:
Occlusal: Shade 130
Occlusal Fossa: Stains orange

2. Schritt / 2nd Step:
Incisal: Shade I1

3. Schritt und Brennen / 3rd Step and firing:
Fissur / fissure: Stains mahogany
Höcker / cusp: Stains white



Endsituation glasiert
Completed restoration glazed

Brennparameter Malfarben- und Charakterisierungsbrand
Firing parameters Stains and characterization firing

B	S	t st	T	H	V ₁	V ₂
403°C	6 min	100°C	790°C	1 min	450°C	789°C
757°F	6 min	180°F	1454°F	1 min	842°F	1453°F

Brennparameter Glanzbrand
Firing parameters Glaze firing

B	S	t st	T	H	V ₁	V ₂
403°C	6 min	100°C	790°C	1–2 min	450°C	789°C
757°F	6 min	180°F	1454°F	1–2 min	842°F	1453°F

Bitte entsprechende Verarbeitungsanleitung beachten. / Please read the corresponding Instructions for Use.