

# IVOCLAR VIVADENT REMOVABLE

*More than dentures.*



## The Stratos® System

Simple articulation with a full range of options.

## STRATOS 100

The **AVERAGE-VALUE** articulator



The average-value articulator for easy and exact tooth set-up in removable denture prosthetics.

## STRATOS 200

The **SEMI-ADJUSTABLE** articulator



The semi-adjustable articulator for universal application in the field of fixed and removable denture prosthetics.







## STRATOS 300

The **INDIVIDUAL** articulator



The individual articulator for functional simulations in conjunction with the high-quality crown and bridge technique as well as implant prosthetics.

## COMPARISON OF THE STRATOS 100, STRATOS 200 AND STRATOS 300 JOINTS

		Stratos 100	Stratos 200	Stratos 300
Condylar path curvature		Radius 12.5	Radius 12.5	Radius 12.5
Protrusive movement		30°, fixed	30°, [15°, 20°, 25°, 35°, 40°, 45°, 60°]*	0 – 60°, adjustable
Bennett movement		15° / 30°, fixed	15° / 30°, fixed	0 – 30°, adjustable
Retrusive movement		35°	35°	35° / 0 – 2 mm, adjustable
Side-shift setting		–	0 – 1.5 mm, adjustable	0 – 1.5 mm, adjustable
Protrusive centric shift		0 – 4 mm, adjustable	0 – 4 mm, adjustable	0 – 4 mm, adjustable

\*Protrusion inserts are available individually or as part of the "personalized" Accessories Assortment.

## UTS 3D TRANSFERBOW

The Universal Transferbow System (UTS) is very efficient, economic and precise. In less than three minutes, an experienced user is able to transfer the arbitrary position of a patient's maxilla into the Stratos.



Indirect method with CE registration joint holder (without UTS 3D transferbow)



With height-adjustable support pins (Type 2) and UTS 3D transferbow





THE AVERAGE-VALUE ARTICULATOR    SIMPLE    PRACTICAL    UNCOMPLICATED

\* included as standard

# STRATOS 200

## THE SEMI-ADJUSTABLE ARTICULATOR    VERSATILE    PRACTICAL    UNIVERSAL

The Stratos 200 is a semi-adjustable articulator with various setting options due to exchangeable joint inserts. It has been especially designed for use in fixed and removable denture prosthetics.

### Indication:

Removable denture prosthetics, fixed crowns and bridges

#### Opening stop

This device prevents the upper frame from tilting forward and creates an optimum working position with the articulator open.



#### Centric locking catch

The centric locking catch can be applied precisely and quickly during the working procedure.

#### Practical separation

By loosening the centric fixation screw, the upper and lower frames of the Stratos 200 may easily be separated and re-fixed.



#### Protrusion screw

This screw enables a centric shift in a protrusive direction.



#### Side-shift movement

The ISS screw of the Stratos 200 permits an immediate side-shift movement of 0 – 1.5mm.

#### 45° incline support holder

This component enables a stable inclined position of the articulator.

#### Incline supports

These incline supports provide ample space when working on the models and permit ergonomic handling of the articulator.

#### Exchangeable joint inserts

With the exchangeable joint inserts, the protrusion and Bennett angles can be individually adjusted in the Stratos 200.

#### Accessories

The system can be individually extended with practice-oriented, well thought-out accessories.

#### Incisal guide rod

The height of the rod can be adjusted and it features millimetre markings and an initial set position.

#### Retrusion

With the Stratos 200, a retrusive movement of 35° can be executed.

#### Splitcast systems\*

The Stratos 200 can be coordinated with other articulators due to its compatibility with Splitcast systems.

#### Base plates

The Stratos 200 is suitable for use with retentive base plates or magnetic base plates.

#### Incisal plate 0°\*

It can be replaced with the 15° plate, the 30° plate or the individual anterior guidance plate.

#### Connection

For CP and FH registration joint holders.

### DELIVERY FORM

1 Basic articulator. 1 Incline support holder. 1 Plaster protection. 1 Base block (retentive). 1 Incisal indicator. 2 Protrusion inserts 30°. 1 Incisal plate 0°

### THE HIGHLIGHTS

#### Exchangeable joint inserts




#### Centric locking catch



#### Protrusion screw



#### Side-shift setting



\* included as standard


# STRATOS 300

THE INDIVIDUAL ARTICULATOR HIGH-QUALITY VERSATILE FLEXIBLE

The Stratos 300 is a state-of-the-art precision articulator with individual setting options. With its wealth of possible variations, it represents the optimum basis for exact functional simulations for demanding, high-quality dentures.

## Indication:

Removable denture prosthetics, fixed crowns and bridges, implant-supported dentures, functional diagnostics



## DELIVERY FORM

1 Basic articulator. 1 Incline support holder. 1 Plaster protection. 1 Base block (retentive). 1 Incisal indicator. 1 Incisal plate 0°

## THE HIGHLIGHTS

### Bennett angle setting



### Centric locking catch



### Protrusion screw



### 120° opening stop



### Side-shift setting



\* included as standard



# Stratos® Accessories



## Mounting Plate Options:

- Screw retained plastic mounting plates (white)
- Black plastic screw retained split-cast base for magnetic mounting with metal discs (rubber ring included)
- Magnetic mounting plates for use with calibrated or non-calibrated bases



**Calibration Kit:** provides precise calibration of one articulator to another for accurate transfer of mounted models.



**Incisal tables:** fixed 15° and 30°, and custom incisal guide tables for personalized anterior guidance.



**Instrument carrier:** to attach the horizontal guide, setup table, 2-D setup template, and bite fork support.

## Mounting and Setup Accessories




**Horizontal guide:** for the orientation of dentulous and edentulous models according to Bonwill's triangle. (Attached with instrument carrier)



**Flat setup table:** for easy, average model orientation of models according to Bonwill's triangle. (Attached with instrument carrier)



**2-D setup template:** for posterior tooth positioning for average valve setup, taking the curves of Spee and Wilson into account. (Attached with instrument carrier)



**3-D setup template:** for correct setup when using personalized model mounting with the transferbow, taking the curves of Spee and Wilson into account.

## UTS Transferbow and Accessories



**UTS Transferbow:** the ideal instrument for easy transfer of personalized patient data according to the reference planes (FH/CP).



**Registration joint holder (Campers Plane):** the entire Transferbow no longer needs to be sent to the laboratory. The registration joint holder ensures accurate model transfer.



**Registration joint holder (Frankfort Horizontal):** for orientation of the UTS Transferbow registrations according to the Frankfort horizontal (FH).



**Plane indicator:** indicates the reference planes during personalized model orientation.



**Bite fork support:** supports the bite fork during model mounting. (Attached with instrument carrier)



**UTS Transferbow support pins:** for direct connection of the UTS.



**Centric Tray:** records centric and vertical jaw relationship for complete, partial, and implant removable prosthetics.