Indications
Recommended for inlays, onlays, 3/4 crowns, posts, short and long span bridges, PFM crowns.

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Yellow, gold-based dental metal ceramic casting alloy, type 4

Produktformate
Instruktioner för bruk
Instrucciones de uso
Modo d'emploi
Gebruiksaanwijzing
Brugsanvisning
Instrucciones de uso
Brugsanvisning
Instrucciones de uso
Modo d'emploi
Instrucciones de uso
Brugsanvisning
Käyttöohjeet

Manufacturer USA
Ivoclar Vivadent Inc.
175 Pineview Drive
Amherst, NY 14228
Tel. +1 800 533 6825
Fax +1 716 691 2285
www.ivoclarvivadent.us

Manufacturer Canada
Ivoclar Vivadent Inc.
175 Pineview Drive
Amherst, NY 14228
Tel. +1 800 533 6825
Fax +1 716 691 2285
www.ivoclarvivadent.us

Manufacturer Europe
Ivoclar Vivadent GmbH
Brenschtr. 16
Postfach 223
A-6706 Burs
Tel. +43 5552 624 49
www.ivoclarvivadent.com

Manufacturer Asia
Ivoclar Vivadent Pty. Ltd.
1 - 5 Overseas Drive
P.O. Box 367
Noble Park, Vic. 3174
Tel. +61 3 979 595 99
Fax +61 3 979 596 45
www.ivoclarvivadent.com.au

Manufacturer Australia
Ivoclar Vivadent Pty. Ltd.
1 - 5 Overseas Drive
P.O. Box 367
Noble Park, Vic. 3174
Tel. +61 3 979 595 99
Fax +61 3 979 596 45
www.ivoclarvivadent.com.au

Manufacturer Mexico
Ivoclar Vivadent S.A. de C.V.
139025 Nahuatl (82)
Tel. +52 (55) 50 62 10 29
Fax +52 (55) 50 62 10 29
www.ivoclarvivadent.com.mx

Manufacturer Colombia
Ivoclar Vivadent Marketing Ltd.
Calle 134 No. 7-8-83, Of. 520
Bogotá
Tel. +57 (1) 627 33 99
Fax +57 (1) 627 33 99
www.ivoclarvivadent.com

Manufacturer Brazil
Ivoclar Vivadent Ltda.
Rua Geraldo Flausino Gomes,
78 - 6ª andar, Cj. 61/62
Bairro: Brooklin Novo
CEP: 04575-060 - São Paulo - SP
Tel. +55 (11) 3466 0800
Fax +55 (11) 3466 0840
www.ivoclarvivadent.com.br

Manufacturer India
Ivoclar Vivadent Marketing Ltd.
503504 Raheja Plaza
15 B Shah Industrial Estate
Veera Desai Road, Andheri (West)
Mumbai, 400 053
Tel. +91 (22) 2673 0302
Fax +91 (22) 2673 0301
www.ivoclarvivadent.ru

Manufacturer Russia
Ivoclar Vivadent Marketing Ltd.
171 Chin Sweet Road
Singapore 169877
Tel. +65 6535 6175
Fax +65 6535 4991

Manufacturer Spain
Ivoclar Vivadent S.A.
c/ Emilio Muñoz nº 15
E-28037 Madrid
Tel. +34 913 75 78 20
Fax +34 913 75 78 38
www.ivoclarvivadent.com

Manufacturer Sweden
Ivoclar Vivadent AB
Dahlgren 14
S-169 56 Solna
Tel. +46 (0) 8 514 93 930
Fax +46 (0) 8 514 93 940
www.ivoclarvivadent.se

Manufacturer Turkey
Ivoclar Vivadent Liason Office
Ahi Evran Caddesi no 1
Polaris Is Merkezi Kat: 7
80670 Maslak
Istanbul
Tel. +90 536 577 12 62

Manufacturer UK
Ivoclar Vivadent UK Limited
Ground Floor Compass Building
Feldspar Close
Warrens Business Park
Enderby
Leicester LE19 4SE England
Tel. +44 116 284 78 80
Fax +44 116 284 78 81
www.ivoclarvivadent.co.uk

Manufacturer Mexico
Ivoclar Vivadent S.A. de C.V.
139025 Nahuatl (82)
Tel. +52 (55) 50 62 10 29
Fax +52 (55) 50 62 10 29
www.ivoclarvivadent.com.mx

Manufacturer Colombia
Ivoclar Vivadent Marketing Ltd.
Calle 134 No. 7-8-83, Of. 520
Bogotá
Tel. +57 (1) 627 33 99
Fax +57 (1) 627 33 99
www.ivoclarvivadent.com

Manufacturer Brazil
Ivoclar Vivadent Ltda.
Rua Geraldo Flausino Gomes,
78 - 6ª andar, Cj. 61/62
Bairro: Brooklin Novo
CEP: 04575-060 - São Paulo - SP
Tel. +55 (11) 3466 0800
Fax +55 (11) 3466 0840
www.ivoclarvivadent.com.br

Manufacturer India
Ivoclar Vivadent Marketing Ltd.
503504 Raheja Plaza
15 B Shah Industrial Estate
Veera Desai Road, Andheri (West)
Mumbai, 400 053
Tel. +91 (22) 2673 0302
Fax +91 (22) 2673 0301
www.ivoclarvivadent.ru

Manufacturer Russia
Ivoclar Vivadent Marketing Ltd.
171 Chin Sweet Road
Singapore 169877
Tel. +65 6535 6175
Fax +65 6535 4991

Manufacturer Spain
Ivoclar Vivadent S.A.
c/ Emilio Muñoz nº 15
E-28037 Madrid
Tel. +34 913 75 78 20
Fax +34 913 75 78 38
www.ivoclarvivadent.com

Manufacturer Sweden
Ivoclar Vivadent AB
Dahlgren 14
S-169 56 Solna
Tel. +46 (0) 8 514 93 930
Fax +46 (0) 8 514 93 940
www.ivoclarvivadent.se

Manufacturer Turkey
Ivoclar Vivadent Liason Office
Ahi Evran Caddesi no 1
Polaris Is Merkezi Kat: 7
80670 Maslak
Istanbul
Tel. +90 536 577 12 62

Manufacturer UK
Ivoclar Vivadent UK Limited
Ground Floor Compass Building
Feldspar Close
Warrens Business Park
Enderby
Leicester LE19 4SE England
Tel. +44 116 284 78 80
Fax +44 116 284 78 81
www.ivoclarvivadent.co.uk

Manufacturer Mexico
Ivoclar Vivadent S.A. de C.V.
139025 Nahuatl (82)
Tel. +52 (55) 50 62 10 29
Fax +52 (55) 50 62 10 29
www.ivoclarvivadent.com.mx

Manufacturer Colombia
Ivoclar Vivadent Marketing Ltd.
Calle 134 No. 7-8-83, Of. 520
Bogotá
Tel. +57 (1) 627 33 99
Fax +57 (1) 627 33 99
www.ivoclarvivadent.com

Manufacturer Brazil
Ivoclar Vivadent Ltda.
Rua Geraldo Flausino Gomes,
78 - 6ª andar, Cj. 61/62
Bairro: Brooklin Novo
CEP: 04575-060 - São Paulo - SP
Tel. +55 (11) 3466 0800
Fax +55 (11) 3466 0840
www.ivoclarvivadent.com.br

Manufacturer India
Ivoclar Vivadent Marketing Ltd.
503504 Raheja Plaza
15 B Shah Industrial Estate
Veera Desai Road, Andheri (West)
Mumbai, 400 053
Tel. +91 (22) 2673 0302
Fax +91 (22) 2673 0301
www.ivoclarvivadent.ru

Manufacturer Russia
Ivoclar Vivadent Marketing Ltd.
171 Chin Sweet Road
Singapore 169877
Tel. +65 6535 6175
Fax +65 6535 4991

Manufacturer Spain
Ivoclar Vivadent S.A.
c/ Emilio Muñoz nº 15
E-28037 Madrid
Tel. +34 913 75 78 20
Fax +34 913 75 78 38
www.ivoclarvivadent.com

Manufacturer Sweden
Ivoclar Vivadent AB
Dahlgren 14
S-169 56 Solna
Tel. +46 (0) 8 514 93 930
Fax +46 (0) 8 514 93 940
www.ivoclarvivadent.se

Manufacturer Turkey
Ivoclar Vivadent Liason Office
Ahi Evran Caddesi no 1
Polaris Is Merkezi Kat: 7
80670 Maslak
Istanbul
Tel. +90 536 577 12 62

Manufacturer UK
Ivoclar Vivadent UK Limited
Ground Floor Compass Building
Feldspar Close
Warrens Business Park
Enderby
Leicester LE19 4SE England
Tel. +44 116 284 78 80
Fax +44 116 284 78 81
www.ivoclarvivadent.co.uk

ISO 13485
Quality Management System Certified

MADE IN U.S.A.
unles otherwise indicated
Made in Canada
(fabrique au Canada)

LOT

Rev. Date 3/10
Rev. 3

Ivoclar Vivadent
passion vision innovation

EN INSTRUCTIONS FOR USE

MODELLATION
Design the framework in a reduced anatomical shape taking the planned veneer into consideration. Single crowns require a thickness of minimum 0.3 mm; abutment crowns thickness a minimum of 0.5 mm. Make sure the framework demonstrates adequate stability of shape. Avoid sharp angles. Design the connector areas to be adequate for the position and alloy being used.

SPRUNG
Provide the modeled bridge framework or coping with sprues of a suitable size. Use the direct or indirect technique being sure that the reservoir is positioned in the heat center. The connection sprues between the reservoir and the coping should be 2.5-3.0 mm in length and width.

INVESTMENT
Weigh the wax pattern including the sprue to determine the quantity of the alloy to be used. (See wax conversion sheet/formula: weight x density = gr. of alloy). Use investment following the manufacturer's instructions.

BURN-OUT
The suggested burnout temperature: 750-820C/1380-1510F

MELTING AND CASTING
Use a separate carbon/ceramic crucible for each alloy and preheat the crucible (ceramic) in the burnout furnace. Used and new alloys must be in a ratio of 1:1. Depending on the type of casting machine, follow the manufacturers instructions for use. If you use Ivoclar Vivadent Magic Wand, set propane to 0.35 bar/5 psi and oxygen to 0.7 bar/10 psi pressure. Keep the alloy in the reducing atmosphere of the flame between the inner and outer cones. Do not use flux. After casting bench cool to room temperature.

Casting Temperature: 1260-1320C/2300-2410F

METAL PREPARATION
Carefully divest and clean the object with Al₂O₃. Do not use a hammer for divesting the object to prevent deformation. Finish the framework with carbide burs or with ceramic-bonded grinding instruments. Avoid inhalation of dust during grinding!

OXIDATION
Blast surface with 50-100 micron Al₂O₃ at max. 4.5 bar/65 psi pressure before oxidation. Subsequently, clean with ultrasonic or steam. Place the object on the firing tray and provide adequate support. Place the tray in a porcelain furnace set at a low temperature of 650C/1200F and increase the furnace temperature to 1010C/1850F without vacuum with 5 min. hold time at top temperature.

HEAT TREATMENT
Hardening: 480C/900F for 15 minutes; air cool.

SOLDERS AND FLUXES
Design the soldering patty as small as possible and preheat it in the furnace at approximately 600C/1112F. The soldering gap should be the same thickness as the soldering strip. Allow the object to cool slowly after soldering.

Pre-Solder: Y-2 Ceramic Solder, High Fusing Bondal Flux
Post-Solder: .615 Fine Solder, Low Fusing White Gold Solder, .585 Fine Solder, Bondal Flux
Laser weld material: Laser Ceramic Yellow

POLISHING
After glazing remove oxide and flux residue and finish/polish the framework with rubber finishers/polishers.

INDICATIONS
Recommended for inlays, onlays, 3/4 crowns, crowns, telescope crowns, conus crowns, posts, short and long span bridges, PFM crowns.

CONTRAINDICATIONS
For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended.

SIDE EFFECTS
In individual cases, sensitivity or allergies to components of this alloy may occur.

INTERACTIONS
Galvanic effects may occur between different alloys in the same oral environment.
For additional information look into the alloy property chart.

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