Inlays, onlays, 3/4 crowns, telescopes crowns, conus crowns, short and long span bridges, implants, superstructure, partial dentures.

Table with 10 columns: Au, Ag, Cu, Sn, Zn, In, Ga, Pb, Ni, Other. Values range from 88.3 to 1.0.


Composition

Instructions for Use, Gebrauchsinformation, Instrucciones de uso, Bruksanvisning, Modé d'emploi, Instrucciones de uso, Produktionsformate, Instrucciones d'uso.

Yellow, gold-based dental casting alloy, Type 4 Harmony® Extra Hard

- USA: Ivoclar Vivadent Inc., 175 Pineview Drive, Amherst, NY 14228. Tel: +1 800 533 6825.
Brazil: Ivoclar Vivadent Ltda., Rua Gerardo Flausino Gomes, 78 - 6.º andar, Cx. 61162, Bairro: Brooklin Novo, CEP: 04575-060 - São Paulo - SP. Tel: +55 (11) 3466 0800.
Canada: Ivoclar Vivadent Inc., 2785 Skyway Avenue, Unit 1, Mississauga, Ontario L4W4Y3. Tel: +1 905 238 5700.
China: Ivoclar Vivadent Marketing Ltd., Rm 603 Xuen Tang International Business Plaza No. 798 Zhao Jia Bang Road, Shanghai 200030. Tel: +86 21 5456 0776.
India: Ivoclar Vivadent Marketing Ltd., 503/504 Raheja Plaza, Veera Desai Road, Andheri (West), Mumbai, 400 053. Tel: +91 (22) 2673 0302.
Italy: Ivoclar Vivadent s.r.l. & C. s.a.s., Via Gustav Flora 32, I-39025 Naturno (BZ), Italy. Tel: +39 0473 67 01 11.
Japan: Ivoclar Vivadent K.K., 1-28-24-4F Horogo Bunkyo-ku, Tokyo 113-0033. Tel: +81 3 6903 3535.
Mexico: Ivoclar Vivadent S.A. de C.V., Av. Mazatlan No. 61, Piso 2, Col. Condesa, 06170 Mexico City, D.F., Mexico. Tel: +52 (55) 50 62 10 29.
New Zealand: Ivoclar Vivadent Ltd., 12 Omega St, Albany PO Box 5243 Wellesley St, Auckland. Tel: +64 9 630 02 06.
Poland: Ivoclar Vivadent Sp. z o.o., ul. Jana Pawla 1 78, Leica Center, 01-74319 Warszawa, Poland. Tel: +48 22 635 54 96.
Russia: Ivoclar Vivadent Marketing Ltd., Derbenevskaja Nabershnaja, 11, Geb. W, 115114 Moscow. Tel: +7 495 913 66 19.
Singapore: Ivoclar Vivadent Marketing Ltd., 171 Chin Swee Road, Singapore 169877. Tel: +65 6535 6775.
Spain: Ivoclar Vivadent S.A., c/ Emilio Muñoz Nº 15, E-28037 Madrid, Spain. Tel: +34 913 75 78 20.
Sweden: Ivoclar Vivadent AB, Dahlvägen 14, S-169 56 Solna, Sweden. Tel: +46 (0) 8 514 93 940.
Switzerland: Ivoclar Vivadent Ltd., Polaris is Merket Kat: 7, 80670 Maslik, Istanbul, Turkey. Tel: +90 536 577 12 62.
USA: Ivoclar Vivadent UK Limited, Ground Floor Compass Building, Feldspar Close, Warrens Business Park, Enderby, Leicestershire, LE19 4SE England, UK. Tel: +44 116 284 78 80.

ISO 13485 Quality Management System Certified. Made in U.S.A. unless otherwise indicated. Made in Canada (fabrique au Canada). Rev. Date 2/09 Rev. 4.

EN INSTRUCTIONS FOR USE

MODELLATION Wax to full contour for crown and bridge frame design. For composite, build up the framework in a reduced anatomic shape taking the planned veneer into consideration.

SPRUNG Provide the modeled bridge framework or coping with sprues of a suitable size. Use the direct or indirect technique being sure that the reservoir is positioned in the heat center.

INVESTMENT Weigh the wax pattern including the sprue to determine the quantity of the alloy to be used. (See wax conversion sheet/formula: weight x density = gr. of alloy).

BURN-OUT The suggested burnout temperature: High heat temperature investment: 650-760C/1200-1400F Low heat temperature investment: 480-540C/900-1000F

MELTING AND CASTING Use a separate carbon/ceramic crucible for each alloy. Used and new alloy must be in a ratio of 1:1. Depending on the type of casting machine, follow the manufacturers instructions for use.

METAL PREPARATION Carefully divest and clean the object with Al2O3 glass beads, or a pickling agent (such as Prevoxx®). Do not use a hammer for divesting the object to prevent deformation.

HEAT TREATMENT Annealing: 675C/1250F for 10 minutes; quench immediately (water) Hardening: 345C/650F for 30 minutes; air cool.

SOLDERS AND FLUXES Design the soldering patty as small as possible and preheat it in the furnace at approximately 600C/1112F. The soldering gap should be the same thickness as the soldering strip.

LASER WELD MATERIAL Laser C&B Yellow

POLISHING After soldering or heat treatment, remove oxide and flux residue and finish and polish the framework with rubber finishers and polishers.

INDICATIONS Recommended for inlays, onlays, 3/4 crowns, telescopes crowns, conus crowns, short and long span bridges, posts, implant superstructure, partial dentures.

CONTRAINDICATIONS For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended.

SIDE EFFECTS In individual cases, sensitivity or allergies to components of this alloy may occur.

INTERACTIONS Galvanic effects may occur between different alloys in the same oral environment.

For additional information look into the alloy property chart.

SPRUNG METHOD

RECOMMENDATIONS DIRECT: single crowns, inlays and onlays INDIRECT: multiple units and multiple single crowns

INSTRUCTIONS: 1. Select a sprue with a reservoir equal to or larger than the thickest cross-section of the restoration. 2. Maintain the reservoir(s) within the heat center of the investment; keep the restoration(s) approximately 5 mm from the end of the investment and no closer than 5 mm from the sides.

SUGGESTIONS: 1. Thickness and design of wax pattern: follow your manufacturers instructions. 2. Use chill-vents when casting heavy and/or large restorations.

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