

Ivoclar Vivadent All-Ceramics

Preparing for cementation

IPS e.max®			
Material	IPS e.max® Press IPS e.max® CAD		IPS e.max® ZirCAD
	Lithium disilicate glass-ceramic (LS ₂)		Zirconium oxide (ZrO ₂)
Indication	Thin Veneers, veneers, occlusal veneers, inlays, onlays, partial crowns, minimally invasive crowns	Crowns, 3-unit bridges up to the 2 nd premolar	Crowns and bridges
Cementation method	adhesive	self-adhesive/conventional ¹	adhesive self-adhesive/conventional
Blasting	–		Cleaning with Al ₂ O ₃ at max. 1 bar
Etching	Option 1: Agitate Monobond Etch & Prime® for 20 s and allow it to react for another 40 s	Option 2: 20 s with IPS® Ceramic Etching Gel	–
Conditioning		60 s with Monobond® Plus ²	60 s with Monobond® Plus –
Cementation system	Variolink® Esthetic, Multilink® Automix ³	SpeedCEM® Plus, Vivaglass® CEM	Multilink® Automix SpeedCEM® Plus, Vivaglass® CEM

The range of products on offer may vary from country to country.

¹ Crown layer thickness of at least 1.5 mm

² Conventional cementation is done without conditioning.

³ Not recommended for veneers.



Please read the corresponding Instructions for Use.



more information
www.cementation-navigation.com

Date Information prepared: 2017-05-24
627461/en

Ivoclar Vivadent All-Ceramics

Preparing for cementation

	IPS e.max®	IPS Empress®								
Material	IPS e.max® Ceram IPS e.max® ZirPress	IPS Empress® Esthetic IPS Empress® CAD								
	Fluorapatite glass-ceramic	Leucite glass-ceramic								
Indication	Veneers	Veneers, inlays, onlays, partial crowns, crowns								
Cementation method	adhesive	adhesive								
Blasting	–	–								
Etching	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%; text-align: center;">Option 1: Agitate Monobond Etch & Prime® for 20 s and allow it to react for another 40 s</td> <td style="width: 50%; text-align: center;">Option 2: 20 s with IPS® Ceramic Etching Gel</td> </tr> <tr> <td></td> <td style="text-align: center;">60 s with Monobond® Plus</td> </tr> </table>	Option 1: Agitate Monobond Etch & Prime® for 20 s and allow it to react for another 40 s	Option 2: 20 s with IPS® Ceramic Etching Gel		60 s with Monobond® Plus	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%; text-align: center;">Option 1: Agitate Monobond Etch & Prime® for 20 s and allow it to react for another 40 s</td> <td style="width: 50%; text-align: center;">Option 2: 60 s with IPS® Ceramic Etching Gel</td> </tr> <tr> <td></td> <td style="text-align: center;">60 s with Monobond® Plus</td> </tr> </table>	Option 1: Agitate Monobond Etch & Prime® for 20 s and allow it to react for another 40 s	Option 2: 60 s with IPS® Ceramic Etching Gel		60 s with Monobond® Plus
Option 1: Agitate Monobond Etch & Prime® for 20 s and allow it to react for another 40 s	Option 2: 20 s with IPS® Ceramic Etching Gel									
	60 s with Monobond® Plus									
Option 1: Agitate Monobond Etch & Prime® for 20 s and allow it to react for another 40 s	Option 2: 60 s with IPS® Ceramic Etching Gel									
	60 s with Monobond® Plus									
Conditioning										
Cementation system	Variolink® Esthetic	Variolink® Esthetic, Multilink® Automix *								

The range of products on offer may vary from country to country.

* Not recommended for veneers.



Please read the corresponding Instructions for Use.



more information
www.cementation-navigation.com