ivoclar vivadent
 Patented

 Recommended for inlays, onlays, 3/4 crowns, crowns, telescopes, crowns, posts, short and long span bridges, PFM crowns, implant superstructure, partial dentures.

Indications
 Recommended for inlays, onlays, 3/4 crowns, crowns, telescopes, crowns, posts, short and long span bridges, PFM crowns, implant superstructure, partial dentures.

Composition

Au	48.7	Pt	39.6	Ag	0.4	Sn	<1.0	Zn	10.6	Bi	<1.0	Ir	0.1	Other	<1.0
----	------	----	------	----	-----	----	------	----	------	----	------	----	-----	-------	------

Instructions for Use
 Instrucciones de uso
 Gebrauchsinformation
 Instruções de uso
 Bruksanvisning
 Modé d'emploi
 Instrucciones de uso
 Gebrauchsinformation
 Instruções de uso
 Bruksanvisning
 Modé d'emploi
 Instrucciones de uso
 Gebrauchsinformation
 Instruções de uso
 Bruksanvisning
 Modé d'emploi

White, gold-based dental metal ceramic casting alloy, Type 4

 Patented

Brazil Ivoclar Vivadent Ltda. Rua Geraldo Flausino Gomes, 78 - 6º andar Cjs. 61/62 Bairro: Brooklin Novo CEP: 04575-060 - São Paulo - SP Tel. +55 (11) 3466 0800 Fax +55 (11) 3466 0840 www.ivoclarvivadent.com.br	Canada Ivoclar Vivadent Inc. 2765 Skymark Avenue, Unit 17 Mississauga, Ontario L4W4Y3 Tel. +1 905 238 5700 Fax +1 905 238 5711 www.ivoclarvivadent.ca	China Ivoclar Vivadent Marketing Ltd. Rm 602 Kuen Tang No. 798 Zhao Jia Bang Road Shanghai 200030 Tel. +86 21 5456 0776 Fax +86 21 6445 1561 www.ivoclarvivadent.com.cn	Colombia Ivoclar Vivadent Marketing Ltd. Calle 134 No. 7-8-83, Of. 520 Bogotá Tel. +57 11 627 33 99 Fax +57 11 633 16 63 www.ivoclarvivadent.com	Austria Ivoclar Vivadent GmbH Brenschtr. 16 Postfach 223 A-6706 Bius Tel. +43 5552 624 49 www.ivoclarvivadent.com	Germany Ivoclar Vivadent GmbH Dr. Adolf-Schneider-Str. 2 D-73479 Ellwangen, Jagst Tel. +49 (0) 79 61 7 8 89-0 Fax +49 (0) 79 61 63 26 www.ivoclarvivadent.de	USA Ivoclar Vivadent Inc. 175 Pineview Drive Amherst, NY 14228 Tel. +1 800 533 6825 Fax +1 716 691 2285 www.ivoclarvivadent.us	France Ivoclar Vivadent SAS B.P. 118 F-73410 Saint-Jorioz Tel. +33 450 88 64 00 Fax +33 450 88 61 52 www.ivoclarvivadent.fr	Italy Ivoclar Vivadent s.r.l. & C. s.a.s. Via Gustav Flora 32 I-39025 Naturno (BZ) Tel. +39 0473 67 01 11 Fax +39 0473 66 77 80 www.ivoclarvivadent.it	Spain Ivoclar Vivadent S.A. c/ Emilio Muñoz Nº 15 E-28037 Madrid Tel. +34 913 75 78 20 Fax +34 913 75 78 38 www.ivoclarvivadent.com	Sweden Ivoclar Vivadent AB Dahlvägen 14 S-169 56 Solna Tel. +46 (0) 8 514 93 930 Fax +46 (0) 8 514 93 940 www.ivoclarvivadent.se	Turkey Ivoclar Vivadent Liaison Office Ahî Evran Caddesi no 1 Polaris Is Merkezi Kat: 7 80670 Maslak Istanbul Tel. +90 536 577 12 62	UK Ivoclar Vivadent UK Limited Ground Floor Compass Building Feldspar Close Warrens Business Park Enderby Leicester LE19 4SE England Tel. +44 116 284 78 80 Fax +44 116 284 78 81 www.ivoclarvivadent.co.uk	Russia Ivoclar Vivadent Marketing Ltd. Derbenevskaya Naberezhnaya 11, Geb. V. 115114 Moscow Tel. +7 495 913 66 19 Fax +7 495 913 66 15 www.ivoclarvivadent.ru	Singapore Ivoclar Vivadent Marketing Ltd. 171 Chin Swee Road #02-01 San Centre Singapore 169877 Tel. +65 6535 6775 Fax +65 6535 4991	Japan Ivoclar Vivadent K.K. 1-28-24 4F Hongo Bunkyo-ku, Tokyo 113-0033 Tel. +81 3 6903 3535 Fax +81 3 5844 3657	Mexico Ivoclar Vivadent S.A. de C.V. Av. Mazatlan No. 61, Piso 2 Col. Condesa 06170 México City, D.F. Tel. +52 (55) 50 62 10 00 Fax +52 (55) 50 62 10 29 www.ivoclarvivadent.com.mx	New Zealand Ivoclar Vivadent Ltd 12 Omega St, Albany PO Box 5243 Wellesley St Auckland Tel. +64 9 630 52 06 Fax +64 9 914 99 90 www.ivoclarvivadent.co.nz	Poland Ivoclar Vivadent Sp. z o.o. ul. Jana Pawła 1 78 D-73479 Ellwangen, Jagst Tel. +48 22 635 54 96 Fax +48 22 635 54 69 www.ivoclarvivadent.pl
---	--	---	---	--	---	---	--	---	--	---	---	---	---	---	---	---	---	--

ISO 13485
 Quality Management System Certified
Made in U.S.A. unless otherwise indicated
 Made in Canada (fabrique au Canada)
 Caution: US Federal Law restricts this device to sale by or on the order of a licensed dentist.
 Rev. Date 9/10
 Rev. 6

 passion vision innovation

INSTRUCTIONS FOR USE
EN
MODELLATION
 Design the framework in a reduced anatomical shape taking the planned veneer into consideration. Single crowns require a thickness of minimum 0.3 mm; abutment crowns thickness a minimum of 0.5 mm. Make sure the framework demonstrates adequate stability of shape. Avoid sharp angles. Design the connector areas to be adequate for the position and alloy being used.
SPRUNG
 Provide the modeled bridge framework or coping with sprues of a suitable size. Use the direct or indirect technique being sure that the reservoir is positioned in the heat center. The connection sprues between the reservoir and the coping should be 2.5-3.0 mm in length and width.
INVESTMENT
 Weigh the wax pattern including the sprue to determine the quantity of the alloy to be used. (See wax conversion sheet/formula: weight x density = gr. of alloy). Use investment following the manufacturer's instructions.
BURN-OUT
 The suggested burnout temperature: 750-820C/1380-1510F
MELTING AND CASTING
 Use a separate ceramic crucible for each alloy and preheat the crucible in the burnout furnace. Used and new alloys must be in a ratio of 1:1. Depending on the type of casting machine, follow the manufacturer's instructions for use. If you use Ivoclar Vivadent Magic Wand, set propane to 0.35 bar/5 psi and oxygen to 0.7 bar/10 psi pressure. Keep the alloy in the reducing atmosphere of the flame between the inner and outer cones. Do not use flux. After casting bend cool to room temperature.
CASTING TEMPERATURE: 1350-1410C/2460-2570F
METAL PREPARATION
 Carefully divest and clean the object with Al₂O₃. Do not use a hammer for divesting the object to prevent deformation. Finish the framework with carbide burs or with ceramic-bonded grinding instruments. Avoid inhalation of dust during grinding!
OXIDATION
 Blast surface with 50-100 micron Al₂O₃ at max. 4.5 bar/65 psi pressure before oxidation. Subsequently, clean with ultrasonic or steam. Place the object on the firing tray and provide adequate support. Place the tray in a porcelain furnace set at a low temperature of 400C/750F and increase the furnace temperature to 950C/1740F without vacuum with 1 min. hold time at top temperature.
HEAT TREATMENT
Hardening: 595C/1100F for 15 minutes; air cool.
SOLDERS AND FLUXES
 Design the soldering pattern as small as possible and preheat it in the furnace at approximately 600C/1112F. The soldering gap should be the same thickness as the soldering strip. Allow the object to cool slowly after soldering.
Pre-Solder: High Fusing White Ceramic Solder, High Fusing Bondal Flux
Post-Solder: .615 Fine Solder, Low Fusing White Gold Solder, .585 Fine Solder, Bondal Flux
Laser weld material: Laser Ceramic White
POLISHING
 After glazing remove oxide and flux residue and finish/polish the framework with rubber finishers/polishers.
INDICATIONS
 Recommended for inlays, onlays, 3/4 crowns, crowns, telescopes, crowns, posts, short and long span bridges, PFM crowns, implant superstructure, partial dentures, conus crowns, posts.
CONTRAINDICATIONS
 For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended.
SIDE EFFECTS
 In individual cases, sensitivity or allergies to components of this alloy may occur.
INTERACTIONS
 Galvanic effects may occur between different alloys in the same oral environment.
For additional information look into the alloy property chart.

SPRUNG METHOD
RECOMMENDATIONS
DIRECT: single crowns, inlays and onlays
INDIRECT: multiple units and multiple single crowns
INSTRUCTIONS:
 1. Select a sprue with a reservoir equal to or larger than the thickest cross-section of the restoration.
 2. Maintain the reservoir(s) within the heat center of the investment; keep the restoration(s) approximately 5 mm from the end of the investment and no closer than 5 mm from the sides.
 3. Make sure that the sprue is connected to the thickest part of the restoration.
 4. The connection between the sprue and the restoration should be flared ("trumpet" configuration), eliminating turbulence of the alloy (causing investment erosion) and facilitating an undisturbed flow of alloy during casting and solidification.
 5. Maintain proper alloy weight in order to eliminate the negative effect of a too large button during solidification. The rule of thumb for proper alloy weight determination is: total wax weight x specific gravity of the alloy = proper alloy weight.
SUGGESTIONS:
 1. Thickness and design of wax pattern: follow your manufacturer's instructions.
 2. Use chill-vents when casting heavy and/or large restorations.

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