


d.SIGN[®] 59

Pd-based dental metal-ceramic alloy, Type 4

<p>EN Instructions for use</p> <p>Pd-based dental metal-ceramic alloy, Type 4</p> <p>DE Gebrauchsinformation</p> <p>Pd-baltige Dentalkeramik-Legierung, Typ 4</p> <p>FR Mode d'emploi</p> <p>Aliaje ceramo-metalic dentarie, à base de Pd, Type 4</p> <p>IT Istruzioni d'uso</p> <p>Legia dentaria per metalceramica à base di Pd, Tipo 4</p> <p>ES Instruciones de uso</p> <p>Aliajeón dentario para metal cerámica basada en Pd, Tipo 4</p> <p>PT Instruções de uso</p> <p>Liga dentária para metalcerâmica à base de Pd, Tipo 4</p> <p>SV Bruksanvisning</p> <p>Pd-baserad dental metal-keramisk legering, Typ 4</p>	<p>DA Brugsanvisning</p> <p>Pd-baseret dental metal-keramik legering, Type 4</p> <p>FI Käyttöohjeet</p> <p>Pd-pohjainen metallkeramiainnennäköislegiirinä, Tyyppi 4</p> <p>NO Bruksanvisning</p> <p>Pd-basert dentalkeramisk legering, Type 4</p> <p>NL Gebruiksaanwijzing</p> <p>Pd-gebaseerd metaalkeramische legering, Type 4</p> <p>EL Οδηγίες Χρήσεως</p> <p>odontotικό μετάλλοκεραμικό κράμα για Pd-based, Τύπος 4</p> <p>TR Kullanma Talimatı</p> <p>Pd tabanlı dental seramik alaşımı, Tip 4</p> <p>RU Инструкция по применению</p> <p>Pd-основанный стоматологический сплав, тип 4</p> <p>PL Instrukcja stosowania</p> <p>Pd-Stopni metalceramiczny powłoczony, Typ 4</p>
---	--

Condition / Treatment	Type	Indications
Porcelain Fired	4	Onlays, 3/4 Crowns, Crowns, Telescope Crowns, Conus Bars, Attachments, Implant Retained Superstructures, Partial Dentures

Indications for Use

Density	<p>10.9 (g/cm³)</p>
EC REP	<p> ivoclar Vivadent AG Biedersteinstrasse 2 9494 Schaan Liechtenstein Tel. +423 235 35 35 Fax +423 235 33 60 www.ivoclarvivadent.com</p>
REPR	<p> Ivoclar Vivadent Inc. 175 Pineview Drive Amherst, NY 14228 USA Tel. +1 800 533 6825 Fax +1 716 691 2285 www.ivoclarvivadent.com</p>
<p>Date information prepared: 585409/2019-05-23 Rev. 13</p>	 ivoclar vivadent passion vision innovation

Composition (mass %)	<table> <tbody><tr> <td>Au</td> <td>Pt</td> <td>Pd</td> <td>Ag</td> <td>Cu</td> <td>Ga</td> <td>In</td> <td>Ir</td> <td>Re</td> <td>Ru</td> <td>Sn</td> <td>Zn</td> <td>Other</td> </tr> <tr> <td>-</td> <td><1.0</td> <td>59.2</td> <td>27.9</td> <td>-</td> <td>-</td> <td>2.7</td> <td>-</td> <td><1.0</td> <td><1.0</td> <td>8.2</td> <td>1.3</td> <td>Li<0.1</td> </tr> </tbody></table>	Au	Pt	Pd	Ag	Cu	Ga	In	Ir	Re	Ru	Sn	Zn	Other	-	<1.0	59.2	27.9	-	-	2.7	-	<1.0	<1.0	8.2	1.3	Li<0.1
Au	Pt	Pd	Ag	Cu	Ga	In	Ir	Re	Ru	Sn	Zn	Other															
-	<1.0	59.2	27.9	-	-	2.7	-	<1.0	<1.0	8.2	1.3	Li<0.1															

  see Instructions

ivoclar vivadent
passion vision innovation

Instructions for Use

<p>II PRODUCT DESCRIPTION</p> <p>Pd-based dental metal-ceramic alloy, Type 4</p> <p>II INDICATIONS*</p> <p>Onlays, 3/4 Crowns, Crowns, Telescope Crowns, Conus Bars, Attachments, Implant Retained Superstructures, Partial Dentures</p> <p>II WAXING/MODELLATION</p> <p>Design the framework in a reduced anatomical shape taking the planned veneer into consideration. Single crowns require a minimum thickness of 0.3 mm. Abutment crowns require a minimum thickness of 0.5 mm. Ensure the framework design provides adequate support for the veneering material. Avoid sharp angles. Connectors must have the required dimensions to provide resistance to deformation. Create large surface areas for plated soldering, with a gap of 0.05–0.2 mm.</p> <p>III SPRING</p> <p>Provide the modelled single-tooth restoration or bridge framework with sprues of a suitable size. In general the reservoir, sprue leads, and connector sprues, whether peak shaped or traditional, must be sized according to the specific technique used. When using the direct or indirect technique be sure that the reservoir is positioned in the indirect denture. The connector sprues between the resin and wax casting should be a maximum of 2.5–3.0 mm in length and width. The wax pattern including the sprues must be weighed in grams or in order to determine the needed amount of alloy. Wax conversion formulas: wax weight (gram) x alloy density = grams of alloy required.</p> <p>III INVESTING</p> <p>Use a phosphate-bonded investment material. Follow the manufacturer's instructions.</p> <p>III PREHEATING / BURN-OUT</p> <p>Recommended burn-out temperature: 750–820 °C / 1380–1510 °F</p> <p>III MELTING AND CASTING</p> <p>Torch: Propane 0.35 bar / 5 psi, Oxygen 0.7 bar / 10 psi Other specifics may be required by the type of casting machine. It is recommended to use a separate and clean ceramic crucible for each alloy. Preheat the ceramic brick in the burnout furnace. The recommended ratio of used material to new material is 1:1. Do not use flux.</p> <p>Casting Temperature: 1350–1400 °C / 2460–2550 °F</p> <p>III FRAMEWORK FINISHING</p> <p>After bench cooling, carefully dust and clean the casting with aluminum oxide (Al₂O₃). Do not use a hammer for dusting. Finish the casting with carbide burs and/or with ceramic-bonded grinding instruments. Blast the surface with 50–100 micron aluminum oxide (AlO₃) 5.2 bar / 75 psi pressure. Subsequently, steam clean or ultrasonic clean with distilled water or ethanol and dry the framework.</p> <p>III OXIDATION</p> <p>Place the framework on the firing tray providing adequate support. To achieve a uniform result follow the oxidation cycle.</p> <p>Temperature: 1010 °C / 1850 °F; Holding time: 10 min; Vacuum: No If the oxide layer is stained, grind and blast the surface again. Repeat the oxide firing. Use the appropriate ceramic veneering material, following the manufacturer's instructions. Highly recommended firing temperature: 1090 °C / 1995 °F</p> <p>III HEAT TREATMENT</p> <p>Hardening: 540 °C / 1005 °F for 15 min; bench cool</p> <p>III SOLDERING AND LASER WELDING</p> <p>The golding gap should not be wider than the thickness of the soldering material. Allow the soldered casting to cool slowly. Use flux sparingly.</p> <p>Pre Solder: SHWFC Post Solder: 585 Fine Gold Solder, LWFG Laser Welding Wire: Laser Ceramic Yellow</p> <p>III POLISHING</p> <p>Carefully remove any oxide and flux residue. Smooth the metal surfaces with rubber polishes. Polish to a high gloss finish using polishing paste. Subsequently, clean using ultrasonic cleaning equipment or careful steam cleaning.</p>

Condition / Treatment	Type	Indications
Porcelain Fired	4	Onlays, 3/4 Crowns, Crowns, Telescope Crowns, Conus Bars, Attachments, Implant Retained Superstructures, Partial Dentures

Density	<p>10.9 (g/cm³)</p>
EC REP	<p> ivoclar Vivadent AG Biedersteinstrasse 2 9494 Schaan Liechtenstein Tel. +423 235 35 35 Fax +423 235 33 60 www.ivoclarvivadent.com</p>
REPR	<p> Ivoclar Vivadent Inc. 175 Pineview Drive Amherst, NY 14228 USA Tel. +1 800 533 6825 Fax +1 716 691 2285 www.ivoclarvivadent.com</p>
<p>Date information prepared: 585409/2019-05-23 Rev. 13</p>	 ivoclar vivadent passion vision innovation

Composition (mass %)	<table> <tbody><tr> <td>Au</td> <td>Pt</td> <td>Pd</td> <td>Ag</td> <td>Cu</td> <td>Ga</td> <td>In</td> <td>Ir</td> <td>Re</td> <td>Ru</td> <td>Sn</td> <td>Zn</td> <td>Other</td> </tr> <tr> <td>-</td> <td><1.0</td> <td>59.2</td> <td>27.9</td> <td>-</td> <td>-</td> <td>2.7</td> <td>-</td> <td><1.0</td> <td><1.0</td> <td>8.2</td> <td>1.3</td> <td>Li<0.1</td> </tr> </tbody></table>	Au	Pt	Pd	Ag	Cu	Ga	In	Ir	Re	Ru	Sn	Zn	Other	-	<1.0	59.2	27.9	-	-	2.7	-	<1.0	<1.0	8.2	1.3	Li<0.1
Au	Pt	Pd	Ag	Cu	Ga	In	Ir	Re	Ru	Sn	Zn	Other															
-	<1.0	59.2	27.9	-	-	2.7	-	<1.0	<1.0	8.2	1.3	Li<0.1															

  see Instructions

ivoclar vivadent
passion vision innovation

Density	<p>10.9 (g/cm³)</p>
EC REP	<p> ivoclar Vivadent AG Biedersteinstrasse 2 9494 Schaan Liechtenstein Tel. +423 235 35 35 Fax +423 235 33 60 www.ivoclarvivadent.com</p>
REPR	<p> Ivoclar Vivadent Inc. 175 Pineview Drive Amherst, NY 14228 USA Tel. +1 800 533 6825 Fax +1 716 691 2285 www.ivoclarvivadent.com</p>
<p>Date information prepared: 585409/2019-05-23 Rev. 13</p>	 ivoclar vivadent passion vision innovation

Condition / Treatment	Type	Indications
Porcelain Fired	4	Onlays, 3/4 Crowns, Crowns, Telescope Crowns, Conus Bars, Attachments, Implant Retained Superstructures, Partial Dentures

Condition / Treatment	Type	Indications
Porcelain Fired	4	Onlays, 3/4 Crowns, Crowns, Telescope Crowns, Conus Bars, Attachments, Implant Retained Superstructures, Partial Dentures

ADDITIONAL SAFETY CONCERNS AND INSTRUCTIONS

III CONTRAINDICATIONS
For patients with known allergy/sensitivity to any major or minor elements of this alloy, consultation with a physician is recommended. Alloy is not to be used for any application not included within the indications.

III SIDE EFFECTS
In individual cases, sensitivity or allergies to elements of this alloy may occur. Ivoclar Vivadent makes no claims regarding the MRI-compatibility of its dental alloys. It is recommended that the patient be made aware of the possibility for dental alloys to affect MRI results and to disclose the presence of dental alloys to the MRI technician prior to conducting a test.

III INTERACTIONS
Galvanic effects may occur between different or dissimilar alloys in the same oral environment.

III CAUTION
Metal vapors and metal dust are harmful if inhaled. Therefore, the use of extraction equipment and/or suitable protective masks is advised!

III STORAGE CONDITION
Store in a dry environment at room temperature.

III DISCLAIMER
This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and is responsible for any purpose not explicitly stated in the Instructions. These regulations also apply if the materials are used in conjunction with products of other manufacturers.

PROCESSING DATA	
Investment Material:	Phosphate-bonded
Preheating/Burn-out Temperature:	750–820 °C / 1380–1510 °F
Crucible:	Ceramic Crucible
Casting Temperature:	1350–1400 °C / 2460–2550 °F
Oxidation:	Temperature: 1010 °C / 1850 °F; Holding time: 10 min; Vacuum: No
CTE:	(25–500 °C): 14.5 x 10 ⁻⁶ (K) (20–600 °C): 14.8 x 10 ⁻⁶ (K)
Recommended Ceramic Material:	IPS Style®/ IPS InLine®/ One, IPS InLine®, IPS Classic®, IPS d.SIGN®
Hardening:	540 °C / 1005 °F for 15 min; bench cool
Pre Solder / Flux:	SHWFC, High Fusing Bondal Flux
Post Solder / Flux:	585, LWFG, Bondal Flux
Laser Welding Wire:	Laser Ceramic Yellow
TECHNICAL DATA (ISO 22674:2016 & ISO 9693-1:2012)	
Type/Color:	4 White
Density (g/cm ³):	10.9
Melting Range (Solidus/Liquidus):	1190–1290 °C / 2175–2355 °F
Elastic Modulus (GPa/psi):	134 19,435,000
Wetters Hardness:	Procetain Fired Hardened
Vickers Strength (MPa/psi):	740 / 107,330 780 / 113,130
0.2% Proof Stress (MPa/psi):	490 / 71,070 640 / 92,820
Elongation (%):	10 5
*See TYPE CLASSIFICATION DUE TO PHYSICAL PROPERTIES	

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