


Description

FRC Postec® Plus is a radiopaque root canal post made of glass-fibre-reinforced composite with a conicity of 5°18'. The adhesive technique used for cementation provides a durable bond between FRC Postec Plus, the core build-up and tooth. The elastic behaviour is comparable to that of dentin. The high translucency allows the fabrication of highly esthetic, metal-free restorations.

Size 0**Size 1****Size 3****Reamer**

High-quality rotary instruments substantially contribute to the success of the treatment. The dimensions of the reamers correspond to the sizes of the FRC Postec Plus posts and thus ensure accuracy of fit in the root canal.

For easier identification, the instruments are colour-coded:

blue → Size 0

white → Size 1

red → Size 3

Composition

FRC Postec Plus is a light-conducting, radiopaque root canal post made of glass fibres. The polymer matrix is composed of aromatic and aliphatic dimethacrylates. It also contains ytterbium trifluoride.

Indications

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Size 0: 32–42

Size 1: 12, 22, 32–42

Size 3: 11, 13, 21, 23, 33, 43

- **Posterior region**

Depending on the clinical situation, Size 0, Size 1 and Size 3 can be used in the posterior region.

Contraindications

- Patients suffering from bruxism or suspected bruxism
- Deep overbite
- Less than 2–3 mm supragingival tooth structure
- Circular, isogingival damage
- If patients are known to be allergic to any of the ingredients, FRC Postec Plus should not be used.

Side effects and other interactions

None known to date.

FRC Postec® Plus

EN	Instructions for Use Glass-fibre reinforced composite root canal post	FI	Käyttöohjeet Lasikuituvahvisteinen komposiittijuurikanavanasta
DE	Gebrauchsinformation Wurzelstift aus Glasfaserverstärktem Composite	NO	Bruksanvisning Røtsift av glasfiberforsterket komposit
FR	Mode d'emploi Tenon radiculaire en composite renforcé aux fibres de verre	NL	Gebruiksaanwijzing Wortelstift van glasvezelversterkt composiet
IT	Istruzioni d'uso Perno radicolare in composito rinforzato con fibre di vetro	EL	Οδηγίες Χρήσεως Ενδορριζικός αξονας από σύνθετη ρητίνη ενισχυμένη με ίνες υάλου
ES	Instrucciones de uso Espiga s radicales de composite reforzado con fibra de vidrio	TR	Kullanma Talimatları Cam lifyle güçlendirilmiş kompozitten yapılmış bir radyoopak kök kanal postu
PT	Instruções de Uso Poste intra-radicular de composito reforçado com fibra de vidro	RU	Инструкция Композитный корневой штифт, армированный стекловолокном

Clinical instructions

For more detailed information, please refer to the separate Instructions for Use of the mentioned products.

1. Preparing the tooth

- Successfully completed endodontic treatment.
- Select the FRC Postec Plus suitable for the tooth to be restored.
- Prepare the coronal residual tooth structure to accommodate the crown (minimum dentin wall thickness ≥ 1.5 mm) and determine the reconstruction method (direct or indirect technique by means of a post and core build-up) depending on the residual dentin structure.
- Determine the required preparation depth and mark the distance on the corresponding instruments by using silicone stoppers.
 - **The remaining root canal filling should demonstrate an apical length no shorter than 4 mm.**
 - **The length of the post within the root canal must at least correspond to the coronal length of the prosthetic restoration.**
- Remove the root canal filling (e.g. Gates-Glidden or Peeso reamer)
- 1,000–5,000 rpm (instrument speed, green-angled handpiece)
- For posts of Size 0, use instruments with a maximum diameter of 0.5 mm.
- For posts of Size 1, use instruments with a maximum diameter of 0.7 mm.
- For posts of Size 3, use instruments with a maximum diameter of 0.9 mm.
- Prepare the canal down to the required depth using the reamer.
- 1,000–5,000 rpm (instrument speed, green-angled handpiece)
- Rinse and clean the reamer during preparation.
- Final preparation of the tooth to accommodate the post and core build-up. If necessary, apply antirotation protection (e.g. retention box, canal inlay).
- Rinse the root canal.

2. Try-in and conditioning of FRC Postec Plus

- Check proper fit of the post.
- Determine the coronal length and shorten the post using rotary diamond grinders (extraorally), if necessary. Do not use scissors or cutting pliers!
- After try-in, clean the post with phosphoric acid etching gel (e.g. Total Etch). Allow the etching gel to react for 60 seconds, then thoroughly rinse with water and dry.
- Silanate the post (e.g. Monobond® Plus). Allow a reaction time of 60 seconds, then carefully dry with an air syringe. Do not touch the surface with your fingers after that.

3. Adhesive cementation with a dual- or self-curing luting composite

- Use for example Multilink® Automix, MultiCore® Flow or SpeedCEM® Plus.
- For conditioning the dentin and cementing the post, please refer to the respective Instructions for Use of the cementation system in use. (If a dual-curing luting composite is used, polymerize from the occlusal surface for at least 20 seconds holding the curing light close to the post).
- **Important: When a luting composite is used in combination with an adhesive system, the luting composite should not be applied to the root canal which has previously been coated with adhesive. The luting composite and adhesive excess may intermix, which can lead to premature curing. As a result, it may not be possible to move the post into the desired position.**

4. Core build-up

- Directly sculpt the build-up using a highly filled composite such as MultiCore Flow or Tetric® EvoCeram Bulk Fill.
- For conditioning the dentin and applying the composite, please refer to the Instructions for Use of the corresponding composite.

5. Preparing the tooth core according to the planned crown

- Observe circular embrasure of the core.
- The preparation margin must be located approx. 2–3 mm in the natural tooth structure.

6. Impression taking and provisional restoration

- Only use eugenol-free cements to place the temporary, since eugenol may inhibit the polymerization of the luting composite.
- If resin-based temporary cements are used, isolate the preparation with Vaseline to prevent the cement from chemically bonding with the composite core build-up.
- During long temporization periods, the provisional should be re-cemented several times to ensure an effective seal at all times.

7. Next appointment

- Incorporate the permanent restoration (e.g. IPS e.max® System) according to the corresponding Instructions for Use.

8. Removing the post in case of a revision

- If the removal of the FRC Postec Plus root canal post is indicated, this can be done using rotary instruments.
- Due to the similar moving behaviour of the rotary cutting instruments within the glass-fibre-reinforced root canal post and the root dentin, this procedure has to be carried out with caution to control the risk of root perforation. Check the individual steps by means of X-rays.
- Working steps:
 - Make an X-ray in the right angle technique to determine the length and size of the post.
 - Remove the build-up component and determine the position of the root canal post
 - Drill a preliminary hole in the centre of the post using an extra-long round bur with a diameter of 0.6 mm.
 - After that, remove the root canal post using the FRC Postec Plus reamer.

Important notes

- If the stipulated clinical indications cannot be observed, the FRC Postec Plus should not be used. Any use for indications other than those stipulated may jeopardize the success of the subsequent prosthetic restoration.
- Do not grind FRC Postec Plus and do not apply mechanical retentions → predetermined breaking points.
- Use only adhesive luting composite systems or self-adhesive resin cements for cementation. We recommend using Multilink Automix, MultiCore Flow or SpeedCEM Plus from Ivoclar Vivadent.
- The cementation procedure requires isolation of the surfaces to be bonded. Adequate relative or absolute isolation using suitable auxiliaries, such as OptraGate® or OptraDam® Plus, is required.

Warnings

Grinding FRC Postec Plus generates glass-fibre dust, which may cause itching

of the skin. Therefore, the use of gloves is recommended. Avoid the inhalation of grinding dust.

Disinfection, cleaning and sterilization

General notes:

Please also observe the country-specific regulations and guidelines regarding the hygiene standards and reprocessing of medical devices in dental practices. Protective gloves and goggles complying with the requirements of directive 89/686/EEC must be worn when handling any used and contaminated instruments.

- Reamers

New reamers are non-sterile and must be subjected to a complete preparation cycle before they are used for the first time. Subsequently, the reamers must be cleaned and disinfected immediately after every use. Additional sterilization in the steam sterilizer is explicitly recommended. Reprocessing in a washer disinfectant (WD) unit is generally preferable to manual reprocessing.

Only a WD (compliant with EN ISO 15883) with proven effectiveness should be used to ensure reliable machine cleaning and disinfection. The operator is responsible for ensuring that the validation, performance re-qualification and periodic routine checks of the reprocessing procedures in the washer disinfectant are defined, documented and carried out appropriately.

Machine cleaning and disinfection

- Cleaning and disinfection of the reamers in a WD according to the instructions of the equipment manufacturer; e.g. Miele with Vario-TD program or Melag with universal program in a washing box 9955 (Komet); cleaning agents and disinfectants: 1.5 g/l Komet DCTherm 9869 mildly alkaline.
- After the cycle has ended, dry with oil-free compressed air.

Manual cleaning and disinfection:

- Clean reamers immediately after use with a brush (nylon bristles) under running water.
- Additionally, manual cleaning can be complemented with ultrasonic cleaning using a suitable cleaning and disinfecting agent (e.g. Komet DC1, 9826 /alkaline, aldehyde-free). The complete immersion time stipulated by the manufacturer must be observed.
Note: Do not exceed 45°C/113°F to prevent protein coagulation.
- For disinfection, immerse the reamers in a disinfecting solution for endodontic instruments, which is deemed appropriate by the manufacturer (e.g. Komet DC1, 9826 /alkaline, aldehyde-free).
- Finally, rinse the reamers with water (ideally demineralized) and dry with oil-free compressed air.

After cleaning and disinfection, visually examine the reamers for intactness (chipping, corroded surfaces) and cleanliness. If macroscopically visible residual contamination is present, repeat the cleaning and disinfection process.

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ivoclar
vivadent[®]
clinical

Sterilization:

- Steam sterilization with vacuum process at 134°C/273.2°F in a device compliant with the stipulations of DIN EN 13060.
- Use a suitable insert (e.g. endodontic sterilization container 556 (Komet)).
- Fractionated pre-vacuum (type B)

The reamers do not show any performance loss due to frequent reprocessing.

- FRC Postec Plus:

The root canal post can be disinfected with medical alcohol.

The root canal post can be sterilized in the autoclave for 5 minutes at 134°C/273.2°F and 2.1 bar.

Keep out of the reach of children!

For use in dentistry only.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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