
# IPS® Investment Ring System

300 g

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## Instructions for Use

### Product description

The IPS Investment Ring System 300g represents the latest investment ring generation for the PRESS technology. It can be used for pressing on long-span frameworks, as well as for the pressing on several smaller restorations in one working step. The Investment Ring System is ideally compatible with existing Ivoclar Vivadent press ingots and press furnaces.

### Indications

- *IPS e.max ZirPress*
  - Pressing on zirconium oxide frameworks up to a bridge span width of 10 units (depending on the dimensions of the dental arch)
  - Pressing on several smaller zirconium oxide frameworks
- *IPS Inline PoM*
  - Pressing on frameworks fabricated of a ceramic alloy up to a bridge span width of 10 units (depending on the dimensions of the dental arch)
  - Pressing on several smaller frameworks fabricated of a ceramic alloy

### Limitations of use

- Casting of alloys
- All applications not explicitly approved as an indication

### Delivery form

- *IPS Investment Ring System 300 g – Starter Kit*
  - 1x Ring Base 300 g
  - 1x Ring Gauge 300 g
  - 1x IPS Silicone Ring 300 g
  - 1x IPS Ingot Tong
  - 1x IPS Sprue Guide 300 g
  - 10x IPS One-Way Plunger 300 g

The components are also individually available.

### General note

- Only ingots of the same shade may be processed in one press cycle.
- Failure to observe the stipulated contraindications and processing instructions may lead to compromised press results.

## Instructions for Use

### Sprueing

Always attach the sprues in the direction of flow of the ceramic and at the thickest part of the wax-up so that smooth flowing of the viscous ceramic during pressing is enabled.

We recommend the following procedure to determine the accurate wax weight:

- Weigh the framework (with the fired liner/opaquer).
- Weigh the IPS ring base 300 g (seal the opening of the ring base with wax).
- Position the contoured restoration on the ring base and attach it with wax. Weigh again.
- The wax weight of the contoured restoration is calculated by deducting the weight of the ring base and the weight of the framework from the total wax weight.

The following sprueing guidelines must be observed:

	Single crowns	Bridges
<b>Wax wire ø</b>	3 mm	3–4 mm
<b>Length of the wax wire</b>	min. 3 mm, max. 12 mm	min. 3 mm, max. 12 mm
<b>Length of the wax wire incl. waxed-up objects</b>	max. 16 mm	max. 16 mm
<b>Sprue attachment point at the waxed-up object</b>	thickest part of the wax-up	thickest part of the wax-up; bridge pontics and each bridge unit
<b>Sprue angle to the waxed-up object</b>	axial	axial
<b>Sprue angle to the ring base</b>	45–60°	45–60°
<b>Design of the attachment points</b>	round and slightly flared, no angles or sharp edges	round and slightly flared, no angles or sharp edges
<b>Distance between the objects</b>	min. 3 mm	min. 3 mm
<b>Distance to the IPS Silicone Ring</b>	min. 10 mm (Check with the IPS Sprue Guide)	min. 10 mm (Check with the IPS Sprue Guide)

### Note:

By marking the position of the restoration on the rings, subsequent divestment may be facilitated.


## Investing

Investing is carried out using either IPS PressVEST (e.g. overnight) or with IPS PressVEST Speed. Please refer to the Instructions for Use of the corresponding investment material for the detailed processing parameters. Please observe the following basic guidelines:

- Do not use a debubbler on the wax objects.
- Mix the investment material. (**Caution:** The investment material contains quartz powder. Therefore, avoid the inhalation of dust)
- Use a suitable instrument for the fine investment of the cavity. Make sure that the delicate wax margins are not damaged.
- Carefully place the IPS Silicone Ring 300 g on the IPS Ring Base 300 g without damaging the wax objects. The silicone ring must sit flush on the ring base.
- Carefully fill the investment ring with investment material up to the marking on the silicone ring. Position the ring gauge with a hinged movement to allow entrapped air to escape through the eccentric hole.
- Allow the investment ring to set without manipulating it.

## Investment material mixing ratio

		Liquid : Distilled Water
IPS PressVEST (300 g)	Single tooth and bridge frameworks	39 ml : 27 ml
IPS PressVEST Speed (300 g)	Single tooth and bridge frameworks	48 ml : 33 ml

## Preheating

After the stipulated setting time of the respective investment material (IPS PressVEST or IPS PressVEST Speed), the investment ring is prepared for preheating as follows:

- Remove the ring gauge and ring base with a turning movement.
- Carefully push the investment ring out of the IPS Silicone Ring.

- Remove rough spots on the bottom surface of the investment ring with a plaster knife. Check the 90° angle. Investment material residue must not enter the sprues. Blow into the sprues if necessary.

	IPS PressVEST	IPS PressVEST Speed
<b>Setting time</b>	at least 60 min. (max. 24 hours)	at least 30 min. / max. 45 min.
<b>Preheating furnace temperature upon placing the investment ring</b>	Room temperature (max. heating rate 5 °C/min)	850 °C
<b>Position of the investment ring in the preheating furnace</b>	at the rear wall, with the opening facing down	at the rear wall, with the opening facing down
<b>IPS e.max ZirPress ingots or IPS InLine PoM ingots</b>	no preheating	no preheating
<b>IPS One-Way Plunger 300 g</b>	no preheating	no preheating
<b>Holding time at final temperature of 850 °C</b>	at least 90 min	at least 90 min

### Note:

If several Speed investments are to be conducted, they should be invested consecutively and placed into the preheating furnace at an interval of approximately 20 minutes. Make sure that the furnace temperature does not drop too much when placing the investment rings into the preheating furnace. The stipulated holding time counts from the point when the preheating temperature has been reached again.

## Pressing

Carry out the following preparatory steps for pressing before the preheating cycle for the investment ring has been completed:

- Switch on and preheat the press furnace in time. Observe the special press parameters for the IPS Investment Ring System 300 g.
- Provide a cold IPS One-Way plunger
- Provide cold ingots in the required shade and pick them up with the IPS Ingot Tongs.

## Ingot selection

- 3–6 ingots should always be used per press cycle
- with IPS InLine PoM either only small or only large ingots must be used per layer.

### IPS e.max ZirPress

Number of ingots	3 ingots	6 ingots
Wax weight	max. 2.2 g	max. 5.1 g

## IPS InLine PoM

Number of Ingots	3x XS	6x XS	3x S	6x S	3x XS + 3x S
Wax weight	max. 1.4 g	max. 3.0 g	max. 2.1 g	max. 5.1 g	max. 4.1 g

Remove the investment ring from the preheating furnace immediately after completion of the preheating cycle. This step may take max. 1 minute to prevent the investment ring from cooling down too much.

- Place the cold ingots in the hot investment ring using the IPS Ingot Tongs.
- Place the IPS One-Way plunger in the hot investment ring.
- Place the completed investment ring in the center of the hot press furnace using the investment tongs.
- The selected press program is started by pressing "START"
- After the press cycle has ended, immediately remove the investment ring and place it on a cooling grid.

IPS Investment Ring System 300 g	
maximum number of ingots	6 ingots
IPS e.max ZirPress ingots or IPS InLine PoM ingots	cold ingots
IPS One-Way-Plunger 300 g	cold plunger

## Press parameters for IPS e.max ZirPress

(Please observe the parameters for your press furnace!)

Press Furnaces	B	t <sup>∞</sup> °C	T °C/min	H °C	V <sub>1</sub> min	V <sub>2</sub> °C	A °C
EP500 (V2.9)	700	60	940	40	500	940	Progr. 31-51
EP600 / EP600 Combi	700	60	930	40	500	930	250 µm/min
EP 3000	700	60	930	40	500	930	250 µm/min
EP 5000	700	60	930	40	500	930	250 µm/min

## Press parameters for IPS InLine PoM

(Please observe the parameters for your press furnace!)

Press Furnace	B	t <sup>∞</sup> °C	T °C/min	H °C	V <sub>1</sub> min	V <sub>2</sub> °C	A °C
EP500 (V2.9)	700	60	960	40	500	960	Progr. 31-51
EP600 / EP600 Combi	700	60	950	40	500	950	50 µm/min.
EP 3000	700	60	950	40	500	950	50 µm/min.
EP 5000	700	60	950	40	500	950	50 µm/min.

## Divesting


Once the investment ring has cooled to room temperature (approx. 60 minutes), divest as follows:

- Mark the length of the plunger on the cooled investment ring.
- Separate the investment ring using a separating disk. This predetermined breaking point enables reliable separation of the investment material and the ceramic.
- Break the investment ring at the predetermined breaking point using a plaster knife.
- Always use polishing beads to divest the pressed objects. Do not use Al<sub>2</sub>O<sub>3</sub> for rough and fine divestment.
- Rough divestment is carried out at 4 bar (60 psi) pressure, while fine divestment is done at 2 bar (30 psi) pressure.
- Observe the blasting direction and distance to prevent damage to the object margins during divestment.

## Note

For long-span restorations,...

- it is recommended to cut the sprues while the restoration is still in the investment ring in order to relieve tensions before final divestment.
- make sure not to damage the restoration when cutting the various sprues.

## Date information prepared:

09/2008

## Manufacturer:

Ivoclar Vivadent AG, FL-9494 Schaan / Liechtenstein

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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### Ivoclar Vivadent UK Limited
### Ivoclar Vivadent, Inc.