Instructions for Use

EN

■ PRODUCT DESCRIPTION

Ag-based dental casting /metal-ceramic alloy, Type 4

■ INDICATIONS*

Inlays, Onlays, Crowns, Telescopic Crowns, Conus Crowns, Bridges, Wide Bridges, Cast Posts / Cores, Bars, Attachments, Implant Retained Superstructures, Partial Dentures

■ WAXING/MODELLATION

Design the framework in full contour or in a reduced anatomical shape taking the planned veneer into consideration. Use the composite veneering material as recommended. Simple crowns require a minimum thickness of 0.3 mm. Abutment crowns require a minimum thickness of 0.5 mm. Ensure the crown thickness provides adequate support for the veneering material. Avoid sharp angles. Connectors must have the required dimensions to provide resistance to deformation. Create large surface areas for planned soldering, with a gap of 0.05 – 0.2 mm.

■ SPRING

Use the modeled single-tooth restoration or bridge framework with spurs of a suitable size. In general the reservoir, spur leads, and connector spurs, whether pair shaped or traditional, must be sized according to the spring technique used. When using the direct or indirect technique be sure that the reservoir is positioned so that the coil provides adequate support for the veneering material. Avoid sharp angles. Connectors must be 2.5–3.0 mm in length and width. The wax pattern including the spurs must be weighed in grams in order to determine the needed amount of alloy. Wax conversion formula: wax weight (gram) x alloy density = grams of alloy required.

■ INVESTING

Use a phosphate-bonded investment material. Follow the manufacturer's instructions.

■ PREHEATING / BURN OUT

Recommended burn-out temperature: 850 °C

■ MELTING AND CASTING

Temp.: Propane 0.15 bar, Oxygen 0.35 bar
Use the casting machine for the type of casting machine. It is recommended to use a separate and clean carbon /ceramic crucible for each alloy. Preheat the ceramic crucible in the burnout furnace. The recommended temperature is similar to new material 1.1. Do not use flux.
Casting Temperature: 1300 °C

■ FRAMEWORK FINISHING

After bench finishing, carefully clean and cast the casting with aluminum oxide (Al₂O₃). Do not use a hammer or other tool. Clean the framework with carbon bars and/or with ceramic-bonded grinding instruments. Blast the surface with 50–110 micron aluminum oxide (Al₂O₃) at 2.0 bar. For veneering with Lab Composite castings, the framework must be finished and polished. Subsequently, steam clean with clean water and distilled water or ethanol and dry the framework.

■ OXIDATION

Place the framework on the firing tray providing adequate support. To achieve a uniform result follow the firing schedule below.
Temperature: 800 °C, Holding time: 5 min; Vacuum: No
In case of oxide layer is stained, grind and blast the surface again. Repeat the oxide firing. In the case of Zn-containing alloys, avoid splashing of the oxide layer in 10% sulfuric acid or similar pickling solutions is recommended. Note: Always observe proper safety procedures when handling acid. Before application clean framework carefully with water. Use the appropriate ceramic veneering material, following the manufacturer's instructions.
Highest recommended firing temperature: 1000 °C

■ HEAT TREATMENT

Hardening: 500 °C for 15 min; bench cool 850 °C

■ SOLDERING AND LASER WELDING

Flux: High Fusing Bondal Flux
Flux: Bondal Flux
Post Solder: Porta Optimum Lot 710
Flux: Bondal Flux
Laser Welding Wire: Porta P6

■ POLISHING

Carefully remove any oxide and flux residue. Smooth the metal surfaces with rubber polishers. Polish to a high gloss finish using polishing paste. Subsequently, clean using ultrasonic cleaning equipment or careful to clean using.

ADDITIONAL SAFETY CONCERNS AND INSTRUCTIONS

■ CONTRAINDICATIONS

For patients with known allergy/sensitivity to any major or minor elements of this alloy, consultation with a physician is recommended. Alloy is not to be used for any application not included within the indications.

■ SIDE EFFECTS

In certain individuals, sensitivity or allergies to elements of this alloy may occur. Ivoclar Vivadent makes no claims regarding the MRI-compatibility of its dental alloys. It is recommended that the patient be aware of the possibility for dental alloys to affect MRI results and to disclose the presence of dental alloys to the MRI technician prior to conducting a test.

■ INTERACTIONS

Galvanic effects may occur between different or dissimilar alloys in the same oral environment.

■ CAUTION

Metal vapors and metal dust are harmful if inhaled. Therefore, the use of extraction equipment and /or suitable protective masks is advised!

■ STORAGE CONDITION

Store in a dry environment at room temperature.

■ DISCLAIMER

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and for any purpose not explicitly stated in the Instructions. These regulations also apply if the material is used in conjunction with products of other manufacturers.

PROCESSING DATA	phosphate-bonded
Investment Material:	phosphate-bonded
Preheating / Burn-out Temperature:	850 °C
Crucible:	Carbon crucible / Ceramic crucible
Casting Temperature:	1300 °C
Oxidation:	Temperature: 800 °C; Holding time: 5 min; Vacuum: No
CTE:	(25–500 °C): 16.2 x 10 ⁻⁶ /K
Hardening/Veneering Materials:	low-fusing, high-expanding, leucite-containing metal ceramics
Hardening:	500 °C for 15 min; bench cool
Pre Solder / Flux:	Porta Lot 1090 W High Fusing Bondal Flux
Post Solder / Flux:	Porta Optimum Lot 710 Bondal Flux
Laser Welding Wire:	Porta P6
TECHNISCHE DATA (ISO 22674:2016 & ISO 9693-1:2012)	
Type / Color:	4 White
Density (g/cm ³):	10,7
Schmelztemperatur (Solidus/Liquidus):	1100–1170 °C
Elastic Modulus (GPa):	95
Vickers Hardness:	190 220 240 240
Tickers-Härte:	190 220 240 240
Tickers-Härte (MPa):	640 660 660 660
0,2% Dehnungs-Stress (MPa):	350 450 450 450
Elongation (%):	20 11 10 10
* Siehe TYP-KLASSIFIZIERUNG DUE TO PHYSICAL PROPERTIES	

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