Recommended for the fabrication of jewelry and other ornamental items.

Composition

Instructions for Use
Istruzioni d'uso
Gebrauchsinformation
Instruções de uso
Modo d'emploi
Brugsanvisning
Käyttöohjeet

Instructions for Use
Istruzioni d'uso
Gebrauchsinformation
Instruções de uso
Modo d'emploi
Brugsanvisning
Käyttöohjeet

Casting Gold

24K

Ivoclar Vivadent Worldwide

- Australia, Brazil, Canada, China, Colombia, France, Germany, India, Italy, Japan, Mexico, New Zealand, Poland, Russia, Sweden, Switzerland, Taiwan, Thailand, Turkey, UK, USA, Vietnam

ISO 13485 Quality Management System Certified Made in U.S.A. www.ivoclarvivadent.com

INSTRUCTIONS FOR USE

MODELLATION Wax to full contour for design. Avoid sharp angles. Design the connector areas to be adequate for the position and alloy being used.

INVESTMENT Weigh the wax pattern including the sprue to determine the quantity of the alloy to be used. (See wax conversion sheet/formula: weight x density = gr. of alloy).

BURN-OUT The suggested burnout temperature: High heat temperature investment: 650-760C/1200-1400F Low heat temperature investment: 480-540C/900-1000F

MELTING AND CASTING Use a separate carbon/ceramic crucible for each alloy. Used and new alloy must be in a ratio of 1:1. Depending on the type of casting machine, follow the manufacturers instructions for use.

METAL PREPARATION Carefully divest and clean the object with AL2O3, glass beads, or a pickling agent (such as Prevox®). Do not use a hammer for divesting the object to prevent deformation.

SOLDERS AND FLUXES Design the soldering patty as small as possible and preheat it in the furnace at approximately 600C/1112F. The soldering gap should be the same thickness as the soldering strip.

POLISHING After soldering or heat treatment, remove oxide and flux residue and finish and polish the framework with rubber finishers and polishers.

INDICATIONS Recommended for the fabrication of jewelry and other ornamental items.

CONTRAINDICATIONS For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended.

SIDE EFFECTS In individual cases, sensitivity or allergies to components of this alloy may occur.

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