

Media Release

Schaan, Liechtenstein – 7 October 2010

Designed for systems thinkers

Telio CS Inlay and Telio CS Onlay are the new additions to Ivoclar Vivadent's Telio system

Telio CS Inlay and Telio CS Onlay from Ivoclar Vivadent are light-cured single-component materials which allow dentists to temporize both Class I and II cavities – without having to additionally use a temporary cement. Telio CS Inlay and Onlay round off the Telio system, which is the most comprehensive system for the temporary restoration of teeth available to date.

The soft-elastic final consistency of Telio CS Inlay makes it particularly suitable for use in the temporization of deep, parallel-walled inlay preparations, even if slight undercuts are present. However, dentists may also use Telio CS Inlay to reline prefabricated polycarbonate crowns or seal implant screw channels. Telio CS Onlay features a hard-elastic consistency and is thus the material of choice if larger or less retentive preparations such as onlay preparations are temporized.

Strong features

Telio CS Inlay and Telio CS Onlay are characterized by their ability to wet dental hard tissue and their low polymerization shrinkage. Moreover, the two filling materials which also contain an antibacterial and cariostatic agent, show very low water absorption. These properties reduce the risk of marginal gap formation, caries, discolouration and bacterial growth. Consequently, bacteria-tight, colour-stable temporary restorations with very little odour build-up can be created.

Quick and easy

Telio CS Inlay and Telio CS Onlay are quick and easy to apply as they show little stickiness to the modelling instrument. The removal of the temporary filling materials is just as easy: due to their soft-elastic or hard-elastic consistency they can be easily removed in one piece.

The Telio system

Telio is the most comprehensive system for the temporary restoration of teeth available to date. All the Telio materials are compatible with each other and feature coordinated shades. Consequently, the laboratory products can be combined with the dentistry products, which renders their application much more reliable.

For more information please contact:

Ivoclar Vivadent AG
Bendererstr. 2
9494 Schaan
Principality of Schaan
Tel.: +423 235 35 35
Fax.: +423 235 33 60
E-Mail: info@ivoclarvivadent.com
www.ivoclarvivadent.com

Captions:

(Telio CS Inlay and Onlay.jpg)

Fig. 1: The new temporary filling materials Telio CS Inlay and Onlay

(Telio System_deutsch.jpg)

Fig. 2: The Telio system: The compatible systems solution for the temporary restoration of teeth from Ivoclar Vivadent

Media contact:

Lorenzo Rigliaco
Public Relations Manager
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan
Principality of Liechtenstein
Tel.: +423 235 36 98
Fax: +423 235 36 33
E-Mail: lorenzo.rigliaco@ivoclarvivadent.com