

Programat® P500 & P500/G2

Software Revision Log

Version	Changes from previous version	Status
V8.1 25.08.2017	Software optimized <ul style="list-style-type: none"> • New program structure • New programs (e.g. IPS e.max Ceram) 	Currently delivered version
V8.0 01.04.2016	Software optimized <ul style="list-style-type: none"> • New programs for IPS Style and IPS Ivocolor 	No longer available!
V7.1 01.10.2014	Software optimized <ul style="list-style-type: none"> • New programs for IPS InLine powder opaquer • Various software improvements 	No longer available!
V7.0 05.05.2014	Software optimized <ul style="list-style-type: none"> • New service interval indicator • Various software improvements • New, special speed crystallization program for IPS e.max CAD (P85): The new program replaces the existing IPS e.max CAD Speed Crystallization. The program may only be used in conjunction with the IPS e.max CAD Speed Crystallization Tray (article code 647633AN). Please observe the processing information for the IPS e.max CAD Speed Crystallization. 	No longer available!
V6.0 04.04.2011	Software optimized <ul style="list-style-type: none"> • New Ivoclar Vivadent Programs for IPS InLine One (Programs 121 – 126) • New Ivoclar Vivadent Programs for IPS e.max CAD-on (Programs 83 – 84) 	No longer available!
V5.0 15.04.2010	New Software P500/G2: Software optimized and new features: <ul style="list-style-type: none"> • New icons and optimized user interface • Power Saving Technology • TSP: Thermo Shock Protection • CSP: Cooling Shock Protection • New program-structure: Ivoclar Vivadent and individual programs • Ivoclar Vivadent firing parameters were changed 	No longer available!

Programat® P500 & P500/G2

V3.1 15.01.2008	Software optimized and new functions <ul style="list-style-type: none"> • Updated Ivoclar standard firing programs P1-P120 • New SW update via USB memory stick possible • New function: Preparing the memory stick as external program memory 	No longer available!
V2.2 05.02.2007	Software optimized <ul style="list-style-type: none"> • Adjustment of the parameters in the programs (P23 and P24) 	No longer available!
V2.0 14.06.2006	Software optimized and new functions <ul style="list-style-type: none"> • New representation of group and program names • Information line regarding the current program status • Extended information in the footer • New symbols and new writing functions • Improved sound reproduction • New error description • Improved memory stick recognition 	No longer available!
V1.6 23.02.2006	Software optimized and new functions <ul style="list-style-type: none"> • Editing of numeric values by means of the +/- keys will speed up if the key is pressed longer • Adjustment of the parameters in the programs (P15 and P16) 	No longer available!
V1.5 23.12.2005	Software optimized	No longer available!
V1.03 16.11.2005	Software optimized	No longer available!
V1.02 20.10.2005	Software optimized <ul style="list-style-type: none"> • Touch function improved 	No longer available!
V1.01 13.10.2005	Software optimized <ul style="list-style-type: none"> • Firing parameters were changed 	No longer available!
V1.00 22.09.2005	Basic version	blocked

Version 25.08.2017