

Programat P700 & P700/G2

Software Revision Log

Version	Changes from previous versions	Status
V6.0 04.04.2011	Software optimized <ul style="list-style-type: none">- New Ivoclar Vivadent Programs for IPS InLine One (Programs 121 – 126)- New Ivoclar Vivadent Programs for IPS e.max CAD-on (Programs 83 – 84)	Currently delivered version
V5.0 01.05.2010	Software optimized and new functions: <ul style="list-style-type: none">- Ivoclar Vivadent Firing parameters were changed- Coloured firing curve display	No longer available!
V4.1 01.12.2009	Software optimized and new functions: Improved voice output and optimized remaining time indicator	No longer available!
V4.0 25.08.2009	New Software P700/G2: Software optimized and new features <ul style="list-style-type: none">- Power Saving Technology- CSP: Cooling Shock Protection- Acoustic voice output in 12 languages- New icons and optimized user interface- and many more!	No longer available!
V3.0 15.11.2007	Software optimized and new functions <ul style="list-style-type: none">- Updated Ivoclar standard firing programs P1-P120- Newly added languages: Chinese simplified, Chinese traditional, Swedish, Russian, Polish, Turkish;	No longer available!
V1.04 03.09.2007	Software optimized and new functions Extended USB memory stick recognition	No longer available!

V1.03 04.06.2007	Software optimized	No longer available!
V1.02 01.06.2007	New Function Improved opening and closing movement with smooth start function	No longer available!
V1.01 22.05.2007	New Function Protocols of firing programs include SW information	No longer available!
V1.00 14.05.2007	Basic Version	No longer available!

Version: 08.05.2011