

Programat® P700 & P700/G2

Software Revision Log

Version	Changes from previous version	Status
V8.0 01.04.2016	Software optimized <ul style="list-style-type: none"> New programs for IPS Style and IPS Ivocolor 	Currently delivered version
V7.1 01.10.2014	Software optimized <ul style="list-style-type: none"> New programs for IPS InLine powder opaquer Various software improvements 	No longer available!
V7.0 05.05.2014	Software optimized <ul style="list-style-type: none"> New service interval indicator Various software improvements New, special speed crystallization program for IPS e.max CAD (P85): The new program replaces the existing IPS e.max CAD Speed Crystallization. The program may only be used in conjunction with the IPS e.max CAD Speed Crystallization Tray (article code 647633AN). Please observe the processing information for the IPS e.max CAD Speed Crystallization. 	No longer available!
V6.0 04.04.2011	Software optimized <ul style="list-style-type: none"> New Ivoclar Vivadent Programs for IPS InLine One (Programs 121 – 126) New Ivoclar Vivadent Programs for IPS e.max CAD-on (Programs 83 – 84) 	No longer available!
V5.0 01.05.2010	Software optimized and new functions <ul style="list-style-type: none"> Ivoclar Vivadent Firing parameters were changed Coloured firing curve display 	No longer available!
V4.1 01.12.2009	Software optimized and new functions <ul style="list-style-type: none"> Improved voice output and optimized remaining time indicator 	No longer available!

Programat® P700 & P700/G2

V4.0 25.08.2009	New Software P700/G2: Software optimized and new features <ul style="list-style-type: none"> • Power Saving Technology • CSP: Cooling Shock Protection • Acoustic voice output in 12 languages • New icons and optimized user interface • and many more! 	No longer available!
V3.0 15.11.2007	Software optimized and new functions <ul style="list-style-type: none"> • Updated Ivoclar standard firing programs P1-P120 • Newly added languages: Chinese simplified, Chinese traditional, Swedish, Russian, Polish, Turkish 	No longer available!
V1.04 03.09.2007	Software optimized and new functions <ul style="list-style-type: none"> • Extended USB memory stick recognition 	No longer available!
V1.03 04.06.2007	Software optimized	No longer available!
V1.02 01.06.2007	New Function <ul style="list-style-type: none"> • Improved opening and closing movement with smooth start function 	No longer available!
V1.01 22.05.2007	New Function <ul style="list-style-type: none"> • Protocols of firing programs include SW information 	No longer available!
V1.00 14.05.2007	Basic Version	No longer available!

Version 01.04.2016