

# Programat EP 5000 & EP 5000/G2

## Software Revision Log

Version	Changes from previous versions	Status
<b>V6.1</b> 06.07.2011	<b>Software optimized</b> Optimized program selection for individual pressing programs	<b>Currently</b> delivered version
<b>V6.0</b> 04.04.2011	<b>New Software EP 5000/G2: Software optimized and new features:</b>  Compared to the predecessor model Programat EP 5000 following functions got integrated:  <ul style="list-style-type: none"> <li>- Power Saving Technology</li> <li>- TSP: Thermo Shock Protection</li> <li>- CSP: Cooling Shock Protection</li> <li>- New program-structure (130 Ivoclar Vivadent programs, 500 individual programs)</li> <li>- New Ivoclar Vivadent Program for IPS e.max Press Multi ingots</li> <li>- New Ivoclar Vivadent Programs for IPS InLine One (Programs 121 – 126)</li> <li>- New Ivoclar Vivadent Programs for IPS e.max CAD-on (Programs 83 – 84)</li> </ul>	<b>No longer available!</b>
<b>V5.0</b> 01.05.2010	<b>Software optimized</b> Firing parameters were changed	<b>No longer available!</b>
<b>V3.3</b> 19.05.2009	<b>Software optimized</b> <ul style="list-style-type: none"> <li>- Pressing parameters for IPS e.max Press HT for 100g and 200g IPS Investment ring system</li> <li>- IPF (Intelligent Press Function) parameters for IPS e.max ZirPress (for IPS Investment ring system 300g) and for IPS Inline PoM (for IPS investment ring system 300g)</li> </ul>	<b>No longer available!</b>
<b>V3.2</b> 12.12.2008	<b>Software optimized</b> <ul style="list-style-type: none"> <li>- optimized setting of the Crack Detection System (CDS) for 300 g investment ring in program "IPS InLine PoM - 300g"</li> </ul>	<b>No longer available!</b>

<b>V3.1</b> 14.10.2008	<b>New functions</b> - Selection option for 300 g. investment ring is integrated - Parameter for 300 g. investment ring is integrated	<b>No longer available!</b>
<b>V3.00</b> 15.11.2007	<b>Software optimized and new functions</b> - New STD press program "InLine PoM" - Aesthetic press programs new with stop speed of 180 rather than 250 - Updated Ivoclar standard firing parameters P1-P120 - Newly added languages: Chinese simplified, Chinese traditional, Swedish, Russian, Polish, Turkish	<b>No longer available!</b>
<b>V1.04</b> 28.08.2007	<b>Software optimized and new functions</b> Extended USB memory stick recognition	<b>No longer available!</b>
<b>V1.03</b> 04.06.2007	<b>Software optimized</b>	<b>No longer available!</b>
<b>V1.02</b> 01.06.2007	<b>New function</b> Improved opening and closing movement with smooth start function	<b>No longer available!</b>
<b>V1.01</b> 22.05.2007	<b>New function</b> Protocols of firing programs include SW information	<b>No longer available!</b>
<b>V1.00</b> 14.05.2007	<b>Basic version</b>	<b>No longer available!</b>

Status: 06.07.2011