

Programat® EP 5000 & EP 5000/G2

Software Revision Log

Version	Changes from previous version	Status
V8.0 01.04.2016	Software optimized <ul style="list-style-type: none"> New programs for IPS Style, IPS Ivocolor and IPS e.max Press MT Various software optimizations 	Currently delivered version
V7.1 30.01.2015	Software optimized <ul style="list-style-type: none"> Various software improvements 	No longer available!
V7.0 01.10.2014	Software optimized <ul style="list-style-type: none"> New programs for IPS InLine powder opaquer New programs for IPS e.max Press Multi Various software improvements 	No longer available!
V6.9 05.05.2014	Software optimized <ul style="list-style-type: none"> New service interval indicator Software improvements for service technician Program structure prepared for speed crystallization program for IPS e.max CAD Various software improvements 	No longer available!
V6.1 06.07.2011	Software optimized <ul style="list-style-type: none"> Optimized program selection for individual pressing programs 	No longer available!
V6.0 04.04.2011	New Software EP 5000/G2: Software optimized and new features Compared to the predecessor model Programat EP 5000 following functions got integrated: <ul style="list-style-type: none"> Power Saving Technology TSP: Thermo Shock Protection CSP: Cooling Shock Protection New program-structure (130 Ivoclar Vivadent programs, 500 individual programs) New Ivoclar Vivadent Program for IPS e.max Press Multi ingots New Ivoclar Vivadent Programs for IPS InLine One (Programs 121 – 126) New Ivoclar Vivadent Programs for IPS e.max CAD-on (Programs 83 – 84) 	No longer available!

Programat® EP 5000 & EP 5000/G2

V5.0 01.05.2010	Software optimized <ul style="list-style-type: none"> Firing parameters were changed 	No longer available!
V3.3 19.05.2009	Software optimized <ul style="list-style-type: none"> Pressing parameters for IPS e.max Press HT for 100 g and 200 g IPS Investment ring system IPF (Intelligent Press Function) parameters for IPS e.max ZirPress (for IPS Investment ring system 300 g) and for IPS Inline PoM (for IPS investment ring system 300 g) 	No longer available!
V3.2 12.12.2008	Software optimized <ul style="list-style-type: none"> Optimized setting of the Crack Detection System (CDS) for 300 g investment ring in program "IPS InLine PoM – 300 g" 	No longer available!
V3.1 14.10.2008	New functions <ul style="list-style-type: none"> Selection option for 300 g investment ring is integrated Parameter for 300 g investment ring is integrated 	No longer available!
V3.0 15.11.2007	Software optimized and new functions <ul style="list-style-type: none"> New STD press program "InLine PoM" Aesthetic press programs new with stop speed of 180 rather than 250 Updated Ivoclar standard firing parameters P1-P120 Newly added languages: Chinese simplified, Chinese traditional, Swedish, Russian, Polish, Turkish 	No longer available!
V1.04 28.08.2007	Software optimized and new functions <ul style="list-style-type: none"> Extended USB memory stick recognition 	No longer available!
V1.03 04.06.2007	Software optimized	No longer available!
V1.02 01.06.2007	New function <ul style="list-style-type: none"> Improved opening and closing movement with smooth start function 	No longer available!
V1.01 22.05.2007	New function <ul style="list-style-type: none"> Protocols of firing programs include SW information 	No longer available!
V1.00 14.05.2007	Basic version	No longer available!

Version 01.04.2016