

Programat EP 3000 & EP 3000/G2

Software Revision Log

Version	Changes from previous versions	Status
V6.3 07.02.2012	Software optimized Optimized temperature calibration in conjunction with the ATK2 temperature checking set	Currently delivered version
V6.1 06.07.2011	Software optimized Optimized program selection for individual pressing programs	No longer available!
V6.0 04.04.2011	New Software EP 3000/G2: Software optimized and new features: Compared to the predecessor model Programat EP 3000 following functions got integrated: <ul style="list-style-type: none"> - Power Saving Technology - TSP: Thermo Shock Protection - New program-structure (130 Ivoclar Vivadent programs, 300 individual programs) - New Ivoclar Vivadent Program for IPS e.max Press Multi ingots - New Ivoclar Vivadent Programs for IPS InLine One (Programs 121 – 126) - New Ivoclar Vivadent Programs for IPS e.max CAD-on (Programs 83 – 84) 	No longer available!
V5.0 01.05.2010	Software optimized Ivoclar Vivadent firing parameters were changed	No longer available!

V1.3 19.05.09	Software optimized - Pressing parameters for IPS e.max Press HT for 100g and 200g IPS Investment ring system	No longer available!
V1.2 12.12.08	Software optimized - optimized setting of the Crack Detection System (CDS) for 300 g investment ring in program "IPS InLine PoM - 300g"	No longer available!
V1.1 14.10.2008	New functions - Selection option for 300 g. investment ring is integrated - Parameter for 300 g. investment ring is integrated	No longer available!
V1.00 1.7.2008	Basic version	No longer available!

Version: 07.02.2012