

# Programat® EP 3000 & EP 3000/G2

## Software Revision Log

Version	Changes from previous version	Status
<b>V8.0</b> 01.04.2016	<b>Software optimized</b> <ul style="list-style-type: none"> <li>New programs for IPS Style, IPS Ivocolor and IPS e.max Press MT</li> </ul>	<b>Currently delivered version</b>
<b>V7.1</b> 30.01.2015	<b>Software optimized</b> <ul style="list-style-type: none"> <li>Various software improvements</li> </ul>	<b>No longer available!</b>
<b>V7.0</b> 01.10.2014	<b>Software optimized</b> <ul style="list-style-type: none"> <li>New programs for IPS InLine powder opaquer</li> <li>New programs for IPS e.max Press Multi</li> <li>Various software improvements</li> </ul>	<b>No longer available!</b>
<b>V6.9</b> 05.05.2014	<b>Software optimized</b> <ul style="list-style-type: none"> <li>New service interval indicator</li> <li>Software improvements for service technician</li> <li>Program structure prepared for speed crystallization program for IPS e.max CAD</li> <li>Various software improvements</li> </ul>	<b>No longer available!</b>
<b>V6.3</b> 07.02.2012	<b>Software optimized</b> <ul style="list-style-type: none"> <li>Optimized temperature calibration in conjunction with the ATK2 temperature checking set</li> </ul>	<b>No longer available!</b>
<b>V6.1</b> 06.07.2011	<b>Software optimized</b> <ul style="list-style-type: none"> <li>Optimized program selection for individual pressing programs</li> </ul>	<b>No longer available!</b>

# Programat® EP 3000 & EP 3000/G2

<b>V6.0</b> 04.04.2011	<p><b>New Software EP 3000/G2: Software optimized and new features</b></p> <p>Compared to the predecessor model Programat EP 3000 following functions got integrated:</p> <ul style="list-style-type: none"> <li>• Power Saving Technology</li> <li>• TSP: Thermo Shock Protection</li> <li>• New program-structure (130 Ivoclar Vivadent programs, 300 individual programs)</li> <li>• New Ivoclar Vivadent Program for IPS e.max Press Multi ingots</li> <li>• New Ivoclar Vivadent Programs for IPS InLine One (Programs 121 – 126)</li> <li>• New Ivoclar Vivadent Programs for IPS e.max CAD-on (Programs 83 – 84)</li> </ul>	<b>No longer available!</b>
<b>V5.0</b> 01.05.2010	<p><b>Software optimized</b></p> <ul style="list-style-type: none"> <li>• Ivoclar Vivadent firing parameters were changed</li> </ul>	<b>No longer available!</b>
<b>V1.3</b> 19.05.2009	<p><b>Software optimized</b></p> <ul style="list-style-type: none"> <li>• Pressing parameters for IPS e.max Press HT for 100g and 200g IPS Investment ring system</li> </ul>	<b>No longer available!</b>
<b>V1.2</b> 12.12.2008	<p><b>Software optimized</b></p> <ul style="list-style-type: none"> <li>• Optimized setting of the Crack Detection System (CDS) for 300g investment ring in program "IPS InLine PoM - 300g"</li> </ul>	<b>No longer available!</b>
<b>V1.1</b> 14.10.2008	<p><b>New functions</b></p> <ul style="list-style-type: none"> <li>• Selection option for 300g investment ring is integrated</li> <li>• Parameter for 300g investment ring is integrated</li> </ul>	<b>No longer available!</b>
<b>V1.0</b> 01.07.2008	<p><b>Basic version</b></p>	<b>No longer available!</b>

Version 01.04.2016