

Programat CS

Software Revision Log

Version	Changes compared to the predecessor version	Status
V2.10 14 Nov 2011	Optimized software <ul style="list-style-type: none">- Program write protection has been adapted as follows:<ul style="list-style-type: none">o Program 1, 3, 5, 7 → protection activeo Program 2, 4, 6 → protection inactive- Optimized data backup via PrograBase 2	Currently supplied version
V2.00 04 Apr 2011	New programs <ul style="list-style-type: none">- New Ivoclar Vivadent Programs for IPS e.max CAD-on (Programs 5 & 6)- New Ivoclar Vivadent Program for IPS e.max CAD Impulse (Program 7)	No longer available!
V1.06 22 June 2009	Optimized software Optimized program display for program 3.	No longer available!
V1.5 24 Apr 2009	Optimized software The software has been equipped with the new program 3: "Speed crystallization/glaze firing program for IPS e.max CAD". In addition, the parameters for P1 and P2 have been optimized.	No longer available!
V1.4 18 Dec 2007	Optimized software Smooth furnace head movement (however, this function is only active in units with serial numbers ranging from 500514 to 599999 and from 600599 to 699999)	No longer available!

V1.03 14 May 2007	Optimized software Adjustment of the parameters in the P2 program (new T 840 °C, V ₂ 839 °C)	No longer available!
V1.02 20 Mar 2007	Basic version	No longer available!

Status: 14.11.2011