

EU Safety Data Sheet

Systemp.link

Date of issue / Reference	14.12.2007	lise / Version 4		
Replaces version of	04.11.2005			
Date of printing	07.07.2009	Sheet No. 1583	Version 4	Page 1 of 5
	Edited		Deleted	

Company	Ivoclar Vivadent AG, Bendererstrasse 2, FL - 9494 Schaan Fürstentum Liechtenstein
---------	--

1 Commercial product name and supplier

1.1	Commercial product name / Designation	Systemp.link
1.2	Application / Use	Temporary composite cement
1.3	Producer	Ivoclar Vivadent AG, Bendererstrasse 2, FL - 9494 Schaan Fürstentum Liechtenstein msds@ivoclarvivadent.com
	Supplier	
1.4	TOX emergency number	
	Official	Emergency-Call: +423 / 235 35 35 Ivoclar Vivadent AG, FL-9494 Schaan, Liechtenstein

2 Hazards identification

Uncured material: Direct contact can cause eye and skin irritation. The material is contraindicated if a person is known to be allergic to any of the ingredients of the product.

3 Composition

3.1	Chemical characterization	Paste of dimethacrylates, inorganic fillers, initiators, stabilizers and pigments
3.2	Hazardous components	<p>CAS No. 94-36-0 <1 % Benzoylperoxide EINECS/ELINCS No.: 202-327-6 Xi: Irritant. E: Explosive. R2: Risk of explosion by shock, friction, fire or other sources of ignition. R36: Irritating to eyes. R43: May cause sensitisation by skin contact.</p> <p>< 56 % Dimethacrylates (CAS No. 76558-74-2 and 25852-47-5) Xi: Irritant. R36/37/38: Irritating to eyes, respiratory system and skin.</p>
3.3	Further information	None.

4 First aid measures

4.1	Eye contact	Flush with plenty of water. Consult a physician if irritation persists.
4.2	Skin contact	Wash thoroughly with water.
4.3	Ingestion	No hazards anticipated from swallowing small amounts incidentally to normal handling.

EU Safety Data Sheet

Systemp.link

Date of issue / Reference	14.12.2007	lise / Version 4		
Replaces version of	04.11.2005			
Date of printing	07.07.2009	Sheet No. 1583	Version 4	Page 2 of 5
	Edited		Deleted	

-
- | | | |
|-----|---------------------|--|
| 4.4 | Inhalation | Remove to fresh air. |
| 4.5 | Further information | If you feel unwell, seek medical advice (show this safety data sheet). |
-

5 Fire-fighting measures

- | | | |
|-----|------------------------------|---|
| 5.1 | Suitable extinguishing media | Water fog, carbon dioxide, foam, dry chemicals. |
| 5.2 | Extinguishing media to avoid | None known |
| 5.3 | Further information | None. |
-

- | | | |
|---|------------------------------------|---|
| 6 | Accidental release measures | Clean up mechanically.
Dispose of according to local and national regulations. |
|---|------------------------------------|---|
-

7 Handling and storage

- | | | |
|-----|--------------------------------|---|
| 7.1 | Handling | Only adequately trained personnel should handle this product. |
| 7.2 | Industrial hygiene | Usual hygienic measures for dental practice. |
| 7.3 | Storage | Store at 2-8 °C |
| 7.4 | Place of storage | Avoid exposure to light.
Store in a cool, dry place.
Avoid exposure to light. |
| 7.5 | Fire- and explosion-protection | Not required. |
-

8 Exposure controls / Personal protection

- | | | |
|-------|---------------------------------|--|
| 8.1 | Exposure controls | Good general ventilation should be sufficient. |
| 8.2 | Exposure limit values | None established. |
| 8.3 | Occupational exposure controls | |
| 8.3.1 | Respiratory protection | Not required. |
| 8.3.2 | Hand protection | Rubber gloves.
Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.
Check penetration time for the specific gloves used with the glove manufacturer. |
| 8.3.3 | Eye protection | Safety goggles. |
| 8.3.4 | Other | None. |
| 8.4 | Environmental exposure controls | |
-

EU Safety Data Sheet

Systemp.link

Date of issue / Reference	14.12.2007	lise / Version 4		
Replaces version of	04.11.2005			
Date of printing	07.07.2009	Sheet No. 1583	Version 4	Page 3 of 5
	Edited		Deleted	

9 Physical and chemical properties

9.1 Appearance	Paste	
9.2 Colour	colourless to light yellow	
9.3 Odour	practically odourless	
9.4 Change of physical state		Test method:

	not determined	
9.5 Density		
	Not determined.	
9.6 Vapour pressure		
	not applicable	
9.7 Viscosity		
	not determined	
9.8 Solubility		
Solubility in water	< 0.5 %	
9.9 pH		
	Not determined.	
9.10 Flash point	> 150 °C	
9.11 Ignition temperature		
	not determined	
9.12 Explosion limits		
	Lower:	
	Upper:	
	not applicable	
9.13 Further information		
Part. coeff. n-octanol/water		
Evaporat. rate		
	None.	

10 Stability and reactivity

10.1 Thermal decomposition	None, if used in accordance to instructions.
10.2 Hazardous decomposition products	None under normal conditions of storage and use.
10.3 Conditions / materials to avoid	None.
10.4 Further information	Avoid exposure of product to light.

EU Safety Data Sheet

Systemp.link

Date of issue / Reference	14.12.2007	lise / Version 4		
Replaces version of	04.11.2005			
Date of printing	07.07.2009	Sheet No. 1583	Version 4	Page 4 of 5
	Edited		Deleted	

11 Toxicological information

11.1 Acute toxicity

11.2 Subacute / Chronic toxicity Uncured material: prolonged or frequently repeated skin contact may cause allergic skin reactions.

11.3 Further information

12 Ecological information

12.1 Ecotoxicity No data available.

12.2 Mobility No data available.

12.3 Persistence and degradability No data available.

12.4 Bioaccumulative potential No data available.

12.5 Further information No ecological problems to be anticipated if properly handled and used.
nearly insoluble

13 Disposal considerations

Take to an approved landfill or a waste incineration plant, under conditions approved by the local authority.

13.1 EU waste key 20 01 39

14 Transport information

14.1 Transport at land ADR --- RID ---

UN Number --- Kemler Number ---
Packing Group ---
Proper shipping name

14.2 Transport at sea ADN --- IMDG ---

UN Number ---
EMS --- MFAG ---
Packing Group
Proper shipping name ---
Marine pollutant

EU Safety Data Sheet

Systemp.link

Date of issue / Reference 14.12.2007 lise / Version 4
Replaces version of 04.11.2005
Date of printing 07.07.2009 **Sheet No. 1583** Version 4 Page 5 of 5
Edited Deleted

14.3 Air transport ICAO / IATA-DGR ---
UN Number ---
Proper shipping name ---
Subsidiary Risk ---
Labels ---
Packing Group ---
Passenger airplane Packing Instructions ---
max. ---
Cargo Airplane Packing Instructions ---
max. ---

14.4 Further information Product is not classified as a dangerous good for transport.

15 Regulatory information The product is a medical device according to the EC-directive 93/42/EEC.

15.1 UN number ---
15.2 National regulations ---
15.3 EINECS/ELINCS number ---
15.4 Hazard symbols ---
15.5 Hazard designation Contains Benzoylperoxide, Dimethacrylates (CAS No. 76558-74-2 and 25852-47-5).
15.6 Risk phrases ---
15.7 Safety phrases ---
15.8 AGW value --- ml/m³ (ppm)
15.9 BVD classification (CH) ---
15.10 VbF (D) ---
15.11 Further information None.

16 Other information Version: 4

The above mentioned data correspond to our present state of knowledge and experience. The safety data sheet serves as description of the products in regard to necessary safety measures. The indications do not have the meaning of guarantees on properties.
