







IVOCLAR VIVADENT CAD/CAM Materials for CEREC® and inLab®

Milling Liquid Overview / Übersicht Schleifmittel

		Milling Unit / Schleifeinheit	
		CEREC® and inLab®	CEREC® and inLab® MC XL
			
Material		Milling Concentrate / Schleifzusatz Dentatec	
	IPS Empress CAD	75 ml	150 ml
	IPS e.max CAD	75 ml	150 ml
	IPS e.max ZirCAD	max. 25 ml	max. 50 ml
	IPS AcryCAD	max. 10 ml	max. 15 ml
	Telio CAD	max. 10 ml	max. 15 ml

The mentioned mixing ratios refer to a full tank.

Die aufgeführten Mischungsverhältnisse beziehen sich auf einen vollen Tank.

Les rapports de mélange indiqués s'appliquent pour un réservoir plein.

I rapporti di miscelazione riportati si riferiscono ad un serbatoio pieno.

Los porcentajes de mezcla indicados se aconsejan para un tanque lleno.

As proporções de mistura citadas referem-se a um tanque cheio.



Shake well before use
Vor Gebrauch gut schütteln
Agiter bien avant usage
Agitare il flacone
Agitar el bote antes del uso
Sacode ben o recipiente