

Programat® S1

Software Revision Log

Version	Changes from previous version	Status
V5.0 29.05.2017	Software optimized – new program structure <ul style="list-style-type: none"> • New programs for IPS e.max ZirCAD MT Multi/MT/LT • Extended function with Programat Dosto-Tray • 50 programs available (22 pre-programmed) <ul style="list-style-type: none"> ○ PLEASE NOTE: All individual programs will be deleted in case of an update! • Various software improvements 	Currently delivered version
V3.0 09.10.2014	Software optimized – new program structure <ul style="list-style-type: none"> • New programs for Wieland Dental materials (Zenotec, Zenostar) • Various software improvements 	No longer available!
V2.0 11.10.2012	Software optimized – new program structure <ul style="list-style-type: none"> • Sintering of IPS e.max ZirCAD single-tooth frameworks in 75 minutes (in the past, 90min were required for this procedure) • Sintering of IPS e.max ZirCAD bridge frameworks with up to 4 units in 90 minutes (in the past, 2h 45min were required for this procedure) • Quick sintering of full-contour ZrO2 single-tooth restorations with enhanced translucency in just 2h 55min • Quick sintering of ZrO2 bridges with enhanced translucency in just 4h 30min 	No longer available!
V1.0 07.12.2009	Basic Version of the Programat S1	No longer available!

Version 29.05.2017