

Astrobrush[®]

Lustrous finish for composite restorations

Touching
brilliance



Astrobrush®

Lustrous finish for composite restorations

Astrobrush® allows you to achieve a shiny smooth surface finish in no time at all. The brushes are ideally suited to impart a lustrous shine to pre-polished composite restorations. Three different shapes provide easy access even to difficult-to-reach areas.

Proven quality

- Special fibres offer an exceptional high-gloss finish.
- 3 shapes for flexible use in the anterior and posterior region

Smart working mechanism

- The abrasive medium silicon carbide is integrated into the special fibres of the polishers.
- As the bristles wear down, fresh abrasive medium is released on the surface.

Economical

- Autoclavable and reusable

Indication	Small Cup	Large Cup	Point
Concave surfaces	✓		✓
Convex surfaces		✓	
Interdental	✓		✓
Occlusal	✓	✓	✓



This product forms part of the "Direct Restoratives" product category. Products from this category are optimally coordinated with each other.

Delivery forms:

Assortment

Large Cup (5 x)
Small Cup (5 x)
Point (5 x)

Refills

Each shape is available separately in packages of 6.

Ivoclar Vivadent AG
Benderstr. 2
9494 Schaan
Liechtenstein
Tel. +423 235 35 35
Fax +423 235 33 60
www.ivoclarvivadent.com

Descriptions and data constitute no warranty of attributes. Printed in Liechtenstein
© Ivoclar Vivadent AG, Schaan Liechtenstein

678862/2015-06-18/EN

ivoclar
vivadent
passion vision innovation