


all digital all options

WIELAND
ZENOSTAR®

The zirconia system

Milling / Sintering.
Staining / Glazing.
Finished!

Zenostar® MT 1 – 4



Zenostar® MT 1 – 4

More efficiency and esthetics due to pre-shaded zirconium oxide

Zenostar® MT discs are available from now on in monochrome pre-shaded colour groups. This means esthetic restorations made from Zenostar MT can be produced faster and more efficiently.

The shade concept is designed so that restorations in the most popular A–D tooth colours only require one stain and one glaze firing with IPS Ivocolor after sintering. The desired tooth shade can therefore be reproduced more easily.

With Zenostar MT it is possible to produce full-contour single crowns and up to three unit bridges in the anterior and posterior region.

Highlights:

- High degree of efficiency during production
- Good colour match with the most popular A–D tooth shades
- Natural esthetics in the anterior and posterior regions

Shades:



Delivery forms:

Item No.	Product	
678842	Zenostar MT 1	98.5 - 14 mm / 1
678843	Zenostar MT 2	98.5 - 14 mm / 1
678844	Zenostar MT 3	98.5 - 14 mm / 1
678845	Zenostar MT 4	98.5 - 14 mm / 1
678850	Zenostar MT 1	98.5 - 18 mm / 1
678851	Zenostar MT 2	98.5 - 18 mm / 1
678852	Zenostar MT 3	98.5 - 18 mm / 1
678853	Zenostar MT 4	98.5 - 18 mm / 1

IPS Ivocolor for the correct tooth shade:

Tooth shade	A 1	A 2	A 3	A 3.5	B 1	B 2	C 1	C 2	C 3
Zenostar MT shade	MT 1	MT 2	MT 3	MT 3	MT 1	MT 2	MT 1	MT 4	MT 4
IPS Ivocolor Shade Dentin	Shade 1	Shade 2	Shade 2	Shade 2	Shade 4	Shade 4	Shade 6	Shade 6	Shade 6
Number of stain firings	1	1	1	3	1	1	1	1	2



LEARN MORE:
alldigitalalloptions.com

Manufacturer:
Wieland Dental+ Technik GmbH & Co. KG
Lindenstraße 2, 75175 Pforzheim, Germany
Tel. +49 7231 3705 0
www.wieland-dental.de

Marketing
www.ivoclarvivadent.com/distributors

ivoclar
vivadent
passion vision innovation