

NEW

Programat® P510 & EP 5010

With infrared technology

New software
with optimized
infrared technology

New software now available for download!

Versions: Programat® P510 Version 3.0 or Programat EP 5010 Version 2.0

- Higher reliability
- Increased process safety
- Easier operation



ivoclar
vivadent[®]
passion vision innovation

New Programat® software

With many optimized infrared technology (IRT) functions

The new Programat furnaces (Programat P510 & EP 5010) from Ivoclar Vivadent are equipped with infrared technology. This technological innovation has become well-established and gained significant importance in the dental industry since the market introduction of the Programat P510 series in September 2013. Thanks to the input and feedback of many users, we were able to further develop and optimize the Programat software as well as to add new features to it.



Enhanced, automatic object recognition – even of small restorations

The new software recognizes even the smallest restorations (facets, inlays, overlays, etc.).



New IRT plus operating mode for opaquer firing cycles and complex restorations

If the IRT plus operating mode is selected, the infrared camera regulates the predrying and closing process with a reduced drying rate. As a result, the predrying and closing processes are prolonged. The IRT plus mode is especially recommended for opaquer firing cycles and complex restorations.



Fast selection of IRT operating mode directly on the touch display

This new key function, which is positioned directly in the firing curve display, enables the user to quickly and easily select or toggle between the two available IRT operating modes. The IRT function can also be completely deactivated in this menu, if required.

Where do you find this new software?

The new software can be downloaded from the Software Download Center free of charge:

www.ivoclarvivadent.com/downloadcenter

Versions: Programat P510 Version 3.0 or Programat EP 5010 Version 2.0

Benefit from these many advantages and download the new software to your Programat furnace free of charge now. If you have any questions, please contact your Ivoclar Vivadent representative or the After Sales Service.



This product forms part of our Fixed Prosthetics category. All the products of this category are optimally coordinated with each other.