

Media Release

Schaan (FL) and Cologne (D) – 20th March 2017

Powerful LED for every use

Bluephase Style 20i: Ivoclar Vivadent is introducing a new curing light

Bluephase Style 20i is a cordless high-performance curing light that combines maximum light output with extremely short curing times. Equipped with the Polywave technology of the third, and latest, LED generation, the powerful new curing light allows all light-curing dental materials in the wavelength range from 385 to 515 nm to be polymerized rapidly and without restriction.

The high light intensity of up to 2,000 mW/cm² in the Turbo mode is especially beneficial for polymerizing indirect restorations. It enables short curing times starting from only 5 seconds and supports economic working procedures during everyday practice.

Especially suitable for indirect restorations

Depending on the material thickness, shade and opacity of an indirect restoration, the quantity of light reaching the luting composite can be significantly reduced. Especially for this indication, the high light intensity of Bluephase Style 20i is of decisive advantage. It ensures that enough energy penetrates the crown or inlay even with high-quality all-ceramic restorations made of, for example, IPS e.max or IPS Empress – for an adequate and reliable curing of the light- or dual-curing luting composite. The curing light offers a selection of four adjustable exposure times and two curing programs High Power (1,200 mW/cm²) and Turbo (2,000 mW/cm²) to match the indication. The settings can be selected using one hand only.

Now available in green and with more output

Simultaneously with the introduction of the Bluephase Style 20i, all the existing Bluephase Style curing lights are now offered with an increased light output from 1,100 mW/cm² to 1,200 mW/cm². In addition, the existing range of grey, blue and pink colours has been extended to include a green unit.

Bluephase®, IPS e.max® and IPS Empress® are registered trademarks of Ivoclar Vivadent AG.

Captions:

(Bluephase Style 20i.jpg)

Fig. 1: Bluephase Style 20i: powerful LED for every use

(Bluephase Style Line.jpg)

Fig. 2: Bluephase Style: now available in the 4 colours grey, blue, pink and green and with increased light output

For additional information:

Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan/Liechtenstein
Phone: +423 235 35 35
Fax: +423 235 33 60
E-Mail: info@ivoclarvivadent.com
www.ivoclarvivadent.com

Media Contact:

André Büssers
Public Relations Manager
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan/Liechtenstein
Phone: +423 235 36 98
Fax: +423 235 36 33
E-Mail: andre.buessers@ivoclarvivadent.com