

Patient Information



IPS Empress[®] CAD

High-tech precision makes your smile look so beautiful



IPS Empress[®] CAD


CERAMIC FOR MODERN DENTISTRY



IPS Empress CAD is a high-quality ceramic for fabricating inlays/onlays, veneers and crowns. This material meets the cutting-edge requirements of aesthetics, function and durability. In addition, the material has proven itself a

million times over in the dental office. The IPS Empress CAD all-ceramic is used without a metal substructure and represents a biocompatible alternative – not only for allergic persons.

The leucite-reinforced glass-ceramic blocks are machined using the innovative CAD/CAM technology of the CEREC® unit.





IPS Empress[®] CAD



INLAYS / ONLAYS

Today, there are attractive alternatives to the frequently used amalgam fillings. Only your dentist will recognize a restoration made from IPS Empress CAD.



Small picture:
Back teeth restored with amalgam



Large picture:
The highly aesthetic alternative:
IPS Empress CAD inlays / onlays

Dr. G. Koike, Japan

 IPS Empress[®] CAD



VENEERS

Chipped or discoloured teeth can be repaired with veneers (laminates). True-to-nature IPS Empress CAD veneers will enhance the appearance of your teeth.



Small picture:
Fractured and discoloured front teeth

Large picture:
The highly aesthetic alternative:
IPS Empress CAD veneers

Dr. K. Wiedhahn, Germany

 IPS Empress[®] CAD



CROWNS

Severely damaged teeth can be repaired with a crown. IPS Empress CAD crowns will blend in harmoniously into the natural surroundings.



Small picture:
Composite filling

Large picture:
The highly aesthetic alternative:
IPS Empress CAD posterior crown

Lee Culp, USA



IPS Empress[®] CAD



DON'T LEAVE YOUR NEW TEETH TO CHANCE!

Your dentist will be glad to inform you whether an all-ceramic restoration is the optimum functional and aesthetic solution for you. Ivoclar Vivadent offers you the suitable product depending on your situation.

**TREAT YOURSELF TO A SUCCESSFUL AND BRIGHT
SMILE AND TAKE YOUR TIME. YOU ARE GOING TO
MAKE AN IMPORTANT DECISION.**



IPS Empress[®] CAD

TO MAKE YOUR SMILE LOOK SO BEAUTIFUL

For further information on
IPS Empress CAD and other
all-ceramic materials from
Ivoclar Vivadent visit
www.ivoclarvivadent.com

Practice stamp

Ivoclar Vivadent AG
Bendererstr. 2
FL-9494 Schaan
Principality of Liechtenstein
Tel. +423 / 235 35 35
Fax +423 / 235 33 60
www.ivoclarvivadent.com


ivoclar
vivadent[®]
passion vision innovation