

SERVICE+ | PORTFOLIO MATERIAL & INDICATION

	Fixed Prosthetics				Implant Prosthetics			Removable Prosthetics	Process-Supporting Prosthetics		
	Inlays / Onlays / Veneers	Crowns / Frameworks	Bridges	Telescope crowns	Customized abutments	Implant-supported bars	Implant-supported superstructures	Superstructures ¹ on bars	Digital Denture	Splints	Wax patterns
IPS e.max® ZirCAD		✓	✓	✓			✓				
IPS e.max® CAD	✓	✓	✓				✓				
IPS Empress® CAD	✓	✓									
IPS e.max® CAD Veneering Solutions		✓	✓								
Colado® CAD CoCr 4		✓	✓	✓		✓	✓	✓			
Colado® CAD Ti2		✓	✓	✓							
Colado® CAD Ti5		✓	✓	✓		✓	✓	✓			
Custom Abutment Ti5					✓						
Telio® CAD		✓	✓				✓				
SR Vivodent® CAD									✓ ²		
Ivobase® CAD									✓ ²		
ProArt CAD Splint										✓	
ProArt CAD Wax yellow											✓
ProArt CAD Wax pink									✓ ²		

¹ Superstructures are available only with implant bar from Design & Production.

² Digital Denture will be available as of July 2017 (requirements: Dental Designer 2017 and Digital Denture add-on).