

Programat® S1 1600

Software Revision Log

Version	Changes from previous version	Status
V5.0 29.05.2017	Software optimized – new programs <ul style="list-style-type: none"> • New programs for IPS e.max ZirCAD MT Multi/MT/LT • Extended function with Programat Dosto-Tray • 50 programs available (22 pre-programmed) <ul style="list-style-type: none"> ○ PLEASE NOTE: All individual programs will be deleted in case of an update! • Various software improvements 	Currently delivered version
V4.0 15.09.2015	Basic Version of the Programat S1 1600	No longer available!

Version 29.05.2017