

Programat P300/G2

Software Revision Log

Version	Changes from previous version	Status
V6.0 04.04.2011	Software optimized New Ivoclar Vivadent Programs for IPS InLine One (Programs 121 – 126)	Currently delivered version
V5.2 28.01.2011	Software optimized Update function has been optimized	No longer available!
V5.1 10.09.2010	Software optimized Integration of the Ivoclar Vivadent Programs for the “ <i>IPS e.max CAD-on technique</i> ” (P83 and P84)	No longer available!
V5.0 02.08.2010	Basic Version of the Programat P300/G2 Compared to the predecessor model Programat P300 following functions got integrated: <ul style="list-style-type: none">- Power Saving Technology- TSP: Thermo Shock Protection- New program-structure: 120 Ivoclar Vivadent and 300 individual programs	No longer available!

Version: 08.05.2011