

Ivoclar Vivadent AG
Bendererstrasse 2
FL-9494 Schaan
Principality of Liechtenstein

Tel: +423 235 3535
Fax: +423/235 3633
Extension: +423/235 3595 or 3687
presse@ivoclarvivadent.com



Schaan, October 2007

Press Release

System solutions from Ivoclar Vivadent

Harmony[®] 3 – Compatible with SR Adoro and SR Link

The Harmony 3 casting alloy has been especially developed for use with the SR Adoro veneering composite and SR Link bonding system. This combination achieves a sound and gap-free bond for aesthetic and long-lasting restorations.

In times of a sensitive gold market, alloys, such as Harmony 3, are the right alternative. Given its gold content of 3.5% and low density, Harmony 3 is an economical choice for fabricating veneered and full-cast restorations.

The wide range of indications offers a variety of possibilities: Harmony 3 can be used for inlays, onlays, partial crowns, single crowns, telescope and conus crowns, endodontic posts, as well as for short- and long-span bridges. The good casting and flow properties enable a reliable cast and thus the accurate reproduction of details.

Harmony 3 is characterized by easy and fast finishing of the alloy surface and excellent polishing properties.

Caption: Harmony[®] 3 veneered with SR Adoro

For further information, please contact:
Ivoclar Vivadent AG
Bendererstr. 2
FL-9494 Schaan
Tel.: +423 235 3535
Fax.: +423 235 33 60
e-mail: presse@ivoclarvivadent.com