


IPS e.max[®]

ceram**POWER**

UNLEASH THE
POWER



Introducing IPS e.max Ceram Power – The next generation of **layering ceramic** – designed specifically for today's highly esthetic translucent restorations.

- High value retention – eliminates graying effect
- Dentin and Incisal formulations
- Same layering techniques as traditional Ceram

IPS e.max ceram **POWER**


Powerful esthetics for powerful results.

ivoclar
vivadent[®]
passion vision innovation

Extended range of possibilities with IPS e.max® Ceram

IPS e.max Ceram Power materials have been designed with today's all-ceramic products in mind. In the past, the goal of layering was to add translucency and depth to opaque frameworks. Now, with the high-strength and highly translucent materials on the market, adding translucent layering ceramic can often lead to a grayish appearance.

IPS e.max Ceram Power powders have been formulated for application on translucent structural ceramics such as IPS e.max lithium disilicate and zirconia to keep value levels high and to meet dentist and patient requirements. The Power materials may also be added to opaque substructures for increased brightness when the situation demands it.



Structured materials are more translucent today.



The difference in value and brightness on translucent materials is evident when comparing traditional IPS e.max Ceram to IPS e.max Ceram Power materials.



IPS e.max Ceram Power on IPS e.max Press

681629EN IPS e.max Ceram Power Dentin Intro Kit

The Power Dentin Intro Kit provides a small sampling of the Ceram Power materials along with the Ivoclar stain and glaze products to create simple but beautiful restorations. Contents:

- 1x IPS e.max Ceram Power Dentin A2 20g
- 1x IPS e.max Ceram Power Incisal 1 20g
- 1x Build-Up Liquid allround 60ml
- 1x IPS Ivocolor Shade 1 3g
- 1x IPS Ivocolor Shade Incisal 3 3g
- 1x IPS Ivocolor Essence sunset 1.8g
- 1x IPS Ivocolor Essence mahogany 1.8g
- 1x IPS Ivocolor Glaze Powder FLUO 1.8g
- 1x IPS Ivocolor Glaze Paste FLUO 3g
- 1x IPS Ivocolor Mixing Liquid allround 15ml

681627AN IPS e.max Ceram Power Dentin Kit A-D

The Power Dentin A-D kit contains popular shades of both Dentin and Incisal materials in the Power formulation for use on translucent substructure materials. Contents:

- 1x IPS e.max Ceram Power Dentin A1 20g
- 1x IPS e.max Ceram Power Dentin A2 20g
- 1x IPS e.max Ceram Power Dentin A3 20g
- 1x IPS e.max Ceram Power Dentin A3.5 20g
- 1x IPS e.max Ceram Power Dentin B1 20g
- 1x IPS e.max Ceram Power Dentin B2 20g
- 1x IPS e.max Ceram Power Dentin C1 20g
- 1x IPS e.max Ceram Power Dentin D2 20g
- 1x IPS e.max Ceram Power Incisal 1 20g
- 1x IPS e.max Ceram Power Incisal 2 20g
- 1x IPS e.max Ceram Power Incisal 3 20g

IPS e.max Ceram Power Dentin 20g Refill

- 681328 IPS e.max Ceram Power Dentin 20g A1
- 681329 IPS e.max Ceram Power Dentin 20g A2
- 681330 IPS e.max Ceram Power Dentin 20g A3
- 681331 IPS e.max Ceram Power Dentin 20g A3.5
- 681332 IPS e.max Ceram Power Dentin 20g A4
- 681333 IPS e.max Ceram Power Dentin 20g B1
- 681334 IPS e.max Ceram Power Dentin 20g B2
- 681335 IPS e.max Ceram Power Dentin 20g B3
- 681336 IPS e.max Ceram Power Dentin 20g B4
- 681337 IPS e.max Ceram Power Dentin 20g C1
- 681338 IPS e.max Ceram Power Dentin 20g C2
- 681339 IPS e.max Ceram Power Dentin 20g C3
- 681340 IPS e.max Ceram Power Dentin 20g C4
- 681341 IPS e.max Ceram Power Dentin 20g D2
- 681342 IPS e.max Ceram Power Dentin 20g D3
- 681343 IPS e.max Ceram Power Dentin 20g D4
- 681347 IPS e.max Ceram Power Dentin 20g BL1
- 681348 IPS e.max Ceram Power Dentin 20g BL2
- 681349 IPS e.max Ceram Power Dentin 20g BL3
- 681350 IPS e.max Ceram Power Dentin 20g BL4

IPS e.max Ceram Power Incisal 20g Refill

- 681344 IPS e.max Ceram Power Incisal 20g 1
- 681345 IPS e.max Ceram Power Incisal 20g 2
- 681346 IPS e.max Ceram Power Incisal 20g 3
- 681351 IPS e.max Ceram Power Incisal 20g BL



© 2016 Ivoclar Vivadent, Inc. Ivoclar Vivadent and IPS e.max are registered trademarks of Ivoclar Vivadent, Inc.

683043 Rev.1 1/2016

United States
175 Pineview Drive
Amherst, NY 14228
716.691.0010
800.533.6825
Fax: 716.691.2285
mail@ivoclarvivadent.com

Canada
1 - 6600 Dixie Road
Mississauga, Ontario L5T 2Y2
905-670-8499
800-263-8182
Fax: 905-670-3102
mail@ivoclarvivadent.com

ivoclar
vivadent
passion vision innovation