

CAD/CAM

Block Portfolio for PlanMill®





IPS Empress® CAD

Leucite-reinforced glass-ceramic

IPS Empress CAD Indications

Veneers, inlays, onlays,
partial crowns, single crowns

Multi



I 12



C 14



C 14 L


Shades

A-D



A1 A2 A3 A3.5 B1

Bleach



BL 1 BL 3

HT (High Translucency)



V 12

A-D



A1 A2 A3

Chromascop



100



I 8



I 10



I 12

A-D



A1 A2 A3 A3.5



B1 B2 B3 C2 D3

Chromascop



100 200 300

LT (Low Translucency)



V 12

A-D



A1 A2 A3

Chromascop



100



I 10



I 12



C 14

A-D



A1 A2 A3 A3.5



B1 B2 B3 C2 D3

Chromascop



100 200 300

Bleach



BL 1 BL 2



BL 3 BL 4



Telio™ CAD


Polymethyl methacrylate

Telio CAD Indications

Temporary crowns and bridges with no more than 2 pontics,
temporary restorations on implants, therapeutic restorations
to correct TMJ and occlusal problems (max. 12 month)

LT (Low Translucency)

Shades



B 40 L

A-D



A1 A2 A3 A3.5 B1

Bleach



BL 3



B 55



IPS e.max[®] CAD

Lithium disilicate (LS₂) glass-ceramic

IPS e.max CAD Indications

Thin veneers, veneers, inlays, onlays, partial crowns, crowns, Anterior bridges, multiple-unit bridges, 3-unit-bridge (up to the second premolar as the terminal abutment)

HT (High Translucency)



I 12



C 14

Shades

A-D



Bleach



BL 1 BL 2 BL 3 BL 4

MT (Medium Translucency)



C 14

A-D



A1 A2 A3 B1

Bleach



BL 2 BL 3 BL 4

LT (Low Translucency)



I 12



C 14

A-D



Bleach



BL 1 BL 2 BL 3 BL 4



C 16

B 32

A-D



A1 A2 A3 A3.5 B1 B2 C1 C2 D2

Bleach



BL 2

Impulse



C 14



Opal 1 Opal 2

Ivoclar Vivadent Block Dimensions

IPS Empress® CAD				
Abbreviation	IPS Empress CAD	Size		
		Width	Length	Height
Inlay	I 8	8.0	8.0	15.0
Inlay	I 10	8.5	10.0	15.0
Inlay	I 12	10.0	12.0	15.0
Crown	C 14	12.0	14.0	18.0
Crown	C 14 L	14.4	14.4	18.0

IPS e.max® CAD				
Abbreviation	IPS e.max CAD	Size		
		Width	Length	Height
Inlay	I 12	10.4	12.5	15.0
Crown	C 14	12.4	14.5	18.0
Crown	C16	17.8	15.8	18.0
Bridge	B 32	14.5	14.5	32.0

Telio™ CAD				
Abbreviation	Telio CAD	Size		
		Width	Length	Height
Bridge	B 40 L	15.4	19.0	39.0
Bridge	B 55	15.5	19.0	55.0

100% CUSTOMER SATISFACTION
GUARANTEED!



ivoclarvivadent.com

Call us toll free at 1-800-533-6825 in the U.S., 1-800-263-8182 in Canada.
©2016 Ivoclar Vivadent, Inc. Ivoclar Vivadent and IPS e.max, are registered trademarks of Ivoclar Vivadent, Inc.
664497 Rev. 3 02/2016

