

Issue Date: 30.11.2018 Version: 1.0.0 Revision Date: 29.11.2023 Print Date: 30.11.2018 1. IDENTIFICATION OF THE SUBSTANCE

1. IDENTIFICATION OF THE SUBSTANCE/MIXTURE/PRODUCT AND MANUFACTURER/IMPORTER 1.1 Product identifier:-

Product	name:
Product	number:

SAFETY DATA SHEET (GHS)

Heliomolar Flow

Ivoclar Vivadent AG Bendererstrasse 2 FL-9494 Schaan

Tel: +423 235 35 35 Fax: +423 235 33 60

Ivoclar Vivadent Ltd

www.ivoclarvivadent.co.nz

PRINCIPALITY OF LIECHTENSTEIN

0800 764 766 (National Poison Centre)

Poisons Hotline (24 hours / 7 days)

12 Omega St, Rosedale, Auckland, New Zealand Phone +64 9 914 9999 Fax +64 9 914 9990

556744/ 557029AN/ 557030AN/ 557031AN/ 557032AN/ 557033AN/ 557035AN/701140

- **1.2 Other means of identification:-**Not applicable.
- 1.3 Recommended use of the chemical and restrictions on use:-Not applicable.
 Identified uses:
 Light-curing resin-based dental restorative material
- 1.4 Details of the manufacturer and importer:-Manufacturer:

Importer:

1.5 Emergency phone number:

2. HAZARD(S) IDENTIFICATION

2.1 GHS Classification:-

- Skin Irrit. 2 H315 Causes skin irritation. Eye Irrit. 2 H319 Causes serious eye irritation. Skin Sens. 1 H317 May cause an allergic skin reaction.
- 2.2 GHS Label elements, including precautionary statements:-Hazard Pictogram:

GHS07 Warning Bis-GMA urethane dimethacrylate triethylene glycol dimethacrylate H315 Causes skin irritation. H319 Causes serious eye irritation. H319 Causes serious eye irritation. H317 May cause an allergic skin reaction. P280 Wear protective gloves/protective clothing/eye protection/face protection. P363 Wash contaminated clothing before reuse. P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

Signal word: Hazard-determining components of labelling:

Hazard statements:

Precautionary statements:



passion vision innovation P337+P313 If eye irritation persists: Get medical advice/attention.

2.3 Additional information:

Medical devices as defined in Directive 93/42/EEC and which are invasive or used in direct physical contact with the human body, are exempted from the provisions of Regulation (EC) No 1272/2008 (CLP/GHS) usually if they are in the finished state and intended for the final user. Other hazards:-

Results of PBT and vPvB assessment;

PBT: Not applicable. vPvB: Not applicable.

3 **COMPOSITION/INFORMATION ON INGREDIENTS** Mixture of substances listed below with nonhazardous additions

CAS No.	Classification	Concentration
1565-94-2	Skin Irrit. 2, H315; Eye Irrit. 2, H319; Skin Sens. 1, H317	10-<25%
13760-80-0	Skin Irrit. 2, H315; Eye Irrit.	10-<25%
72869-86-4	Eye Irrit. 2, H319; Skin Sens. 1, H317; Aquatic	10-<25%
109-16-0	Skin Irrit. 2, H315; Eye Irrit. 2, H319; Skin Sens. 1, H317; STOT SE 3, H335	3-<10%
	1565-94-2 13760-80-0 72869-86-4	1565-94-2 Skin Irrit. 2, H315; Eye Irrit. 2, H319; Skin Sens. 1, H317 13760-80-0 Skin Irrit. 2, H315; Eye Irrit. 2, H319; STOT SE 3, H335 72869-86-4 Eye Irrit. 2, H319; Skin Sens. 1, H317; Aquatic Chronic 3, H412 109-16-0 Skin Irrit. 2, H319; Skin Sens. 1, H317; STOT SE 3,

For the full text of the H-Statements mentioned in this Section, refer to Section 16.

FIRST AID MEASURES 4.

4.1 Description of necessary first aid measures:-

	If inhaled:	Supply fresh air and to be sure call for a doctor. In case of unconsciousness place patient stably in side position for transportation.
	In case of skin contact:	Rinse with water
	In case of eye contact:	Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
	If swallowed:	If symptoms persist consult doctor.
4.2	Symptoms caused by exposure:-	No further relevant information available
4.3	Medical attention and special treatment:-	No further relevant information available.
5. 5.1	FIRE FIGHTING MEASURES Suitable extinguishing equipment:- Suitable extinguishing media:	CO2, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
	Unsuitable extinguishing media:	No further relevant information available.
5.2	Specific hazards arising from the substance/mixture/product:-	No further relevant information available.
5.3	Special protective equipment and precautic Special personal protective equipment: Precautions:	ons for fire fighters:- No further relevant information available No further relevant information available

6. ACCIDENTAL RELEASE MEASURES

6.1 Personal precautions, protective equipment and emergency procedures:-



Not required.

6.2 Environmental precautions:-Do not allow to enter sewers/ surface or ground water.

6.3 Methods and materials for containment and cleaning up:-Pick up mechanically. See Section 7 for information on safe handling. See Section 8 for information on personal protection equipment.

See Section 8 for information on personal protection equipmen See Section 13 for disposal information

7. HANDLING AND STORAGE

7.1 Precautions for safe handling:-

Only adequately trained personnel should handle this product. For use in dentistry only.

7.2 Conditions for safe storage, including any incompatibilities:-Store only in the original receptacle Protect from exposure to the light. Protect from heat and direct sunlight. Keep container tightly sealed.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

8.1 Exposure control measures:-

Occupational exposure limits:

Component	CAS No.	Value	Parameters	Basis
Void			Void	

Ingredients with biological limit values:

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

8.2 Biological monitoring:-

The lists valid during the making were used as basis.

Exposure controls / **Personal protective equipment** / **General protective and hygienic measures:** Usual hygienic measures for dental practice.

Wash hands before breaks and at the end of work.

Avoid contact with the eyes and skin.

Keep away from foodstuffs, beverages and feed.

8.3 Control banding:-

Use good industrial hygiene practice and general ventilation.

8.4 Engineering controls:-

Void

8.5 Individual protection measures include PPE:-Eye/face protection:



Skin protection:

Safety glasses

Required Use tightly fitting safety glasses as per Australian Standard AS 1336 and AS/NZS 1337.

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/the chemical mixture. Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

· Material of gloves



passion vision innovation

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· Penetration time of glove material

Not Required

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

Respiratory protection:



9. PHYSICAL/CHEMICAL PROPERTIES

- 9.1 Information on physical/chemical properties:-
- a) Appearance/Form:
- b) Colour:
- c) Odour:
- d) Odour threshold:
- e) pH value:
- f) Melting point/melting range:
- g) Boiling point/boiling range:
- h) Flash point:
- i) Ignition temperature:
- j) Self-igniting:
- k) Danger of explosion:
- I) Upper/lower flammability or explosive limits:
- m) Vapour pressure:
- n) Density at 20°C:
- o) Relative density:
- p) Vapour density:
- q) Evaporation rate:
- r) Solubility in/miscibility with water:
- s) Partition coefficient: n- octanol/water:
- t) Viscosity:
- u) Solids content:

10. STABILITY AND REACTIVITY

- 10.1 Reactivity:-
 - No further relevant information available.

10.2 Chemical stability:-Stable under normal handling and storage conditions. Photoreactive (Danger of polymerization). Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.

- Pastv According to product specification Characteristic Not determined. Not applicable. Not determined Undetermined. Not applicable. Not determined. Product is not self-igniting. Not determined. Lower Not determined. Upper Not determined. Not applicable. 1.94 g/cm³ Not determined. Not applicable.
- Not applicable. Nearly insoluble. Not determined. Dynamic | Not applicable. Kinematic | Not applicable. Not determined.

556744/ 557029AN/ 557030AN/ 557031AN/ 557032AN/ 557033AN/ 557035AN/701140



		passion vision innova
10.3	Possibility of hazardous reactions:-	
	No dangerous reactions known.	
10.4	Conditions to avoid:-	
	Protect from exposure to the light.	
	Keep away from heat and direct sunlight.	
10.5	Incompatible materials:-	
	No further relevant information available.	
10.6	Hazardous decomposition products:-	
	None under normal conditions of storage and	use.
11.	TOXICOLOGICAL INFORMATION	
11.1	Information on toxicological effects:-	
	Acute toxicity / Values relevant for classifi	cation:
	No further relevant information available.	I
	Skin corrosion/irritation:	Irritant to skin and mucous membranes.
	Serious eye damage/eye irritation:	Irritating effect.
	Respiratory or skin sensitization:	Sensitisation possible through skin contact.
	Germ cell mutagenicity:	No further relevant information available.
	Carcinogenicity:	No further relevant information available.
	Reproductive toxicity:	No further relevant information available.
	Specific target organ toxicity - single	No further relevant information available.
	exposure:	No further relevant information available.
	Specific target organ toxicity - repeated	
	exposure: Aspiration hazard:	No further relevant information available.
	Additional information:	No further relevant information available
11.2	Information on possible routes of exposur	
11.2	Short Term (Acute) Exposure:	No further relevant information available.
	Swallowed:	No further relevant information available.
	Eyes:	No further relevant information available.
	Skin:	No further relevant information available.
	Inhaled:	No further relevant information available.
	Long Term (Chronic) Exposure:	No further relevant information available.
	Swallowed:	No further relevant information available.
	Eyes:	No further relevant information available.
	Skin:	No further relevant information available.
	Inhaled:	No further relevant information available.
11.3	Early onset symptoms related to	No further relevant information available.
	exposure:-	
11.4	Delayed health effects from exposure:-	No further relevant information available.
11.5	Exposure levels and health effects:-	No further relevant information available.
11.6	Interactive effects:-	No further relevant information available.
11.7	Other:-	No further relevant information available.
12.	ECOLOGICAL INFORMATION	
12.1	Ecotoxicity:-	
	No further relevant information available.	
12.2	Persistence/degradability:-	
	No further relevant information available.	
12.3	Bioaccumulative potential:-	
40.4	No further relevant information available.	
12.4	Mobility in soil:-	
10 E	No further relevant information available. Other adverse effects:-	
12.5		ition of it to roach ground water, water course or course
	system.	ities of it to reach ground water, water course or sewage
	oyotom.	



Additional ecological information / General notes:-

No further relevant information available.

12.6 Other adverse effects:-

No further relevant information available.

13. DISPOSAL CONSIDERATIONS

13.1 Disposal methods:-

Take to an approved landfill or a waste incineration plant, under conditions approved by the local authority. Disposal must be made according to official regulations.

14.	TRANSPORTINFORMATION		
	UN number ADR / IMDG / IATA:-	Void.	
	UN proper shipping name or technical name:-		
	ADR:	Void.	
	IMDG, IATA:	Void.	
	Transport hazard class(es):	Void.	
	Label:	Void.	
	Packaging group:	Void.	
	Environmental hazards:	Not applicable.	
	Special precautions for user:	Not applicable.	
	Danger code:	Void.	
	EMS Number:	Void.	
	Transport in bulk according to Annex II of	Not applicable.	
	MARPOL73/78 and the IBC Code:		
	Transport/Additional information:	Product is not classified as a dangerous good for transport (ADR, IMDG, IATA).	

15. REGULATORY INFORMATION

15.1 Safety, health and environmental regulations/legislation specific for the substance/mixture/product:-

The product is a medical device according to the Directive 93/42/EEC.

15.2 Chemical safety assessment:

A Chemical Safety Assessment has not been carried out.

16. OTHER INFORMATION

Key to abbreviations/acronyms used in SDS:-

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H319 Causes serious eye irritation.

H335 May cause respiratory irritation.

H412 Harmful to aquatic life with long lasting effects.

Key literature references/data sources used to compile SDS:-Void

Copyright statement:-

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

Abbreviations and acronyms:-

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road). IMDG: International Maritime Code for Dangerous Goods.



IATA: International Air Transport Association. GHS: Globally Harmonised System of Classification and Labelling of Chemicals. EINECS: European Inventory of Existing Commercial Chemical Substances. ELINCS: European List of Notified Chemical Substances. CAS: Chemical Abstracts Service (division of the American Chemical Society). LC50: Lethal concentration, 50 percent. LD50: Lethal dose, 50 percent. Flam. Liq. 2: Flammable liquids, Hazard Category 2. Flam. Lig. 3: Flammable liquids, Hazard Category 3. Acute Tox, 4: Acute toxicity, Hazard Category 4. Skin Irrit. 2: Skin corrosion/irritation. Hazard Category 2. Eye Irrit. 2: Serious eye damage/eye irritation, Hazard Category 2. Repr. 2: Reproductive toxicity, Hazard Category 2. STOT SE 3: Specific target organ toxicity - Single exposure, Hazard Category 3. STOT RE 2: Specific target organ toxicity - Repeated exposure, Hazard Category 2. Asp. Tox. 1: Aspiration hazard, Hazard Category 1.

Disclaimer:-

Ivoclar Vivadent Ltd believes the information in this document to be accurate as at the date of preparation, but, to the maximum extent permitted by law, Ivoclar Vivadent Ltd accepts no responsibility for any loss or damage caused by any person acting or refraining from action because of this information. The provision of this information should not be construed by anyone as a recommendation to use this product. In particular, no one should use any product in violation of any patent or other intellectual proprietary rights or in breach of any statute or regulation.

Users should rely on their own knowledge and inquiries and make their own determination as to the applicability of this information in relation to their particular purposes and specific circumstances. Each user should read this SDS and consider the information in the context of how the product will be handled and used in the workplace and in conjunction with other substances or products.

These data are based on our present knowledge. However, they shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.