

Safety Data Sheet according to WHS Regulations

Printing date 14.01.2021

Version number 7

Revision: 14.01.2021

Not classified as hazardous according to criteria of Australian Safety and Compensation Council.

1 Identification · Product identifier Trade name: IPS e.max Ceram Glaze and Stain Liquid allround · CAS Number: 107-88-0 · EC number: 203-529-7 · Relevant identified uses of the substance or mixture and uses advised against No further relevant information available. · Application of the substance / the mixture Auxiliary for manufacture of dental prothesis • Details of the supplier of the safety data sheet · Manufacturer/Supplier: Ivoclar Vivadent AG Bendererstrasse 2 9494 Schaan PRINCIPALITY OF LIECHTENSTEIN *Tel:* +423 235 35 35 / *Fax:* +423 235 33 60 Importer: Ivoclar Vivadent Pty. Ltd. 1-5 Overseas Drive Noble Park North VIC 3174 Tel: + 61 3 9795 9599 / Fax: + 61 3 9795 9645 • Further information obtainable from: **Regulatory** Affairs sds@ivoclarvivadent.com • Emergency telephone number: 131 126 (Poisons Information Centre - 24 hours / 7 days) **2** Hazard(s) Identification · Classification of the substance or mixture The substance is not classified, according to the Globally Harmonised System (GHS). · Label elements · GHS label elements Void · Hazard pictograms Void · Signal word Void · Hazard statements Void · Other hazards · Results of PBT and vPvB assessment · PBT: Not applicable. · vPvB: Not applicable. **3** Composition and Information on Ingredients · Chemical characterisation: Substances · CAS No. Description

107-88-0 1,3-Butanediol

(Contd. on page 2)

AU

Safety Data Sheet according to WHS Regulations

Printing date 14.01.2021

Version number 7

Revision: 14.01.2021

Trade name: IPS e.max Ceram Glaze and Stain Liquid allround

(Contd. of page 1)

- Identification number(s)
- EC number: 203-529-7

4 First Aid Measures

- · Description of first aid measures
- · General information: No special measures required.
- · After inhalation: Supply fresh air; consult doctor in case of complaints.
- After skin contact:
- Rinse with water.

Generally the product does not irritate the skin.

- After eye contact:
- Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
- After swallowing:

Rinse out mouth and then drink plenty of water.

If symptoms persist consult doctor.

- · Information for doctor:
- · Most important symptoms and effects, both acute and delayed No further relevant information available.
- · Indication of any immediate medical attention and special treatment needed
- No further relevant information available.

5 Fire Fighting Measures

- · Extinguishing media
- · Suitable extinguishing agents:

CO2, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

- Special hazards arising from the substance or mixture No further relevant information available.
- Advice for firefighters
- · Protective equipment: No special measures required.

6 Accidental Release Measures

· Personal precautions, protective equipment and emergency procedures Not required.

- · Environmental precautions: Do not allow to enter sewers/ surface or ground water.
- Methods and material for containment and cleaning up:
- Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- **Reference to other sections** See Section 7 for information on safe handling. See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

7 Handling and Storage

· Handling:

- · Precautions for safe handling
- Only adequately trained personnel should handle this product.
- For use in dentistry only.
- Information about fire and explosion protection: No special measures required.
- · Conditions for safe storage, including any incompatibilities

· Storage:

· Requirements to be met by storerooms and receptacles: Store only in the original receptacle.

(Contd. on page 3)

Safety Data Sheet

according to WHS Regulations

Printing date 14.01.2021

Version number 7

Revision: 14.01.2021

Trade name: IPS e.max Ceram Glaze and Stain Liquid allround

- · Information about storage in one common storage facility: Not required.
- (Contd. of page 2)

- · Further information about storage conditions:
- *This product is hygroscopic. Keep container tightly sealed.*
- Specific end use(s) No further relevant information available.

8 Exposure controls and personal protection

• Additional information about design of technical facilities: No further data; see item 7.

- *Ingredients with limit values that require monitoring at the workplace:* The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- Additional information: The lists valid during the making were used as basis.
- · Exposure controls
- Personal protective equipment:
- General protective and hygienic measures:
- Usual hygienic measures for dental practice and dental laboratories.
- Wash hands before breaks and at the end of work.
- Keep away from foodstuffs, beverages and feed.
- **Respiratory protection:** Not required.
- Protection of hands:

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

• Material of gloves

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer.

Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

· Eye protection: Safety glasses

Information on basic physical and c General Information	chemical properties	
Appearance:		
Form:	Fluid	
Colour:	Colourless	
Odour:	Nearly odourless	
Odour threshold:	Not determined.	
pH-value:	Not determined.	
Change in condition		
Melting point/freezing point:	-50 °C	
Initial boiling point and boiling ra	inge: Undetermined.	
Flash point:	109 °С	

[·] Control parameters

Safety Data Sheet

according to WHS Regulations

Printing date 14.01.2021

Version number 7

Revision: 14.01.2021

Trade name: IPS e.max Ceram Glaze and Stain Liquid allround

	(Co	ntd. of page 3)
· Flammability (solid, gas):	Not applicable.	
· Ignition temperature:	375 °C	
• Auto-ignition temperature:	Not determined.	
· Explosive properties:	Product does not present an explosion hazard.	
· Explosion limits:		
Lower:	Not determined.	
Upper:	Not determined.	
· Vapour pressure:	Not determined.	
· Density at 20 °C:	l g/cm ³	
· Relative density	Not determined.	
· Vapour density	Not determined.	
Evaporation rate	Not determined.	
· Solubility in / Miscibility with		
water:	Fully miscible.	
Partition coefficient: n-octanol/water:	Not determined.	
· Viscosity:		
Dynamic:	Not determined.	
Kinematic:	Not determined.	
• Other information	No further relevant information available.	

10 Stability and Reactivity

· Reactivity No further relevant information available.

- · Chemical stability Stable under normal handling and storage conditions.
- Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
- · Possibility of hazardous reactions No dangerous reactions known.
- · Conditions to avoid No further relevant information available.
- *Incompatible materials:* No further relevant information available.
- · Hazardous decomposition products: None under normal conditions of storage and use.

11 Toxicological Information

- · Information on toxicological effects
- Acute toxicity

· LD/LC50 values relevant for classification:

CAS: 107-88-0 1,3-Butanediol

Oral LD50 18610 mg/kg (rat)

- · Skin corrosion/irritation No irritant effect.
- · Serious eye damage/irritation No irritating effect.
- · Respiratory or skin sensitisation No sensitising effects known.
- Additional toxicological information:
- No further relevant information available.

The substance is not subject to classification according to the latest version of the EU lists.

(Contd. on page 5)

Safety Data Sheet according to WHS Regulations

Printing date 14.01.2021

Version number 7

Trade name: IPS e.max Ceram Glaze and Stain Liquid allround

(Contd. of page 4)

12 Ecological Information

- · Toxicity
- · Aquatic toxicity: No further relevant information available.
- · Persistence and degradability No further relevant information available.
- · Behaviour in environmental systems:
- · Bioaccumulative potential No further relevant information available.
- · Mobility in soil No further relevant information available.
- Additional ecological information:
- · General notes:

Water hazard class 1 (German Regulation) (Assessment by list): slightly hazardous for water Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

- · Results of PBT and vPvB assessment
- · **PBT:** Not applicable.
- · vPvB: Not applicable.
- Other adverse effects No further relevant information available.

13 Disposal considerations

- · Waste treatment methods
- · Recommendation

Take to an approved landfill or a waste incineration plant, under conditions approved by the local authority.

- · Uncleaned packaging:
- *Recommendation: Disposal must be made according to official regulations.*

UN-Number	
ADG, ADN, IMDG, IATA	Void
UN proper shipping name	
ADG, ADN, IMDG, IATA	Void
Transport hazard class(es)	
ADG, ADN, IMDG, IATA	
Class	Void
Packing group	
ADG, IMDG, IATA	Void
Environmental hazards:	
Marine pollutant:	No
Special precautions for user	Not applicable.
Transport in bulk according to Annex II o	f Marpol
and the IBC Code	Not applicable.
Transport/Additional information:	Product is not classified as a dangerous good for
-	transport (ADR, IMDG, IATA).
UN ''Model Regulation'':	Void

(Contd. on page 6)

Safety Data Sheet

according to WHS Regulations

Printing date 14.01.2021

Version number 7

Revision: 14.01.2021

Trade name: IPS e.max Ceram Glaze and Stain Liquid allround

(Contd. of page 5)

15 Regulatory information

· Safety, health and environmental regulations/legislation specific for the substance or mixture

·Australian Inventory of Industrial Chemicals

Substance is listed.

· Standard for the Uniform Scheduling of Medicines and Poisons

Substance is not listed.

· Australia: Priority Existing Chemicals

Substance is not listed.

• GHS label elements Void

· Hazard pictograms Void

· Signal word Void

· Hazard statements Void

· Directive 2012/18/EU

· Named dangerous substances - ANNEX I Substance is not listed.

· National regulations:

• Other regulations, limitations and prohibitive regulations

The product is a medical device according to the Directive 93/42/EEC.

• Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

• Abbreviations and acronyms:

RID: Règlement international concernant le transport des marchandises dangereuses par chemin de fer (Regulations Concerning the International Transport of Dangerous Goods by Rail) ICAO: International Civil Aviation Organisation ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road) IMDG: International Maritime Code for Dangerous Goods IATA: International Maritime Code for Dangerous Goods IATA: International Air Transport Association EINECS: European Inventory of Existing Commercial Chemical Substances CAS: Chemical Abstracts Service (division of the American Chemical Society) LC50: Lethal concentration, 50 percent ID50: Lethal dose, 50 percent PBT: Persistent, Bioaccumulative and Toxic vPvB: very Persistent and very Bioaccumulative