ivoclar vivadeni:

Safety Data Sheet according to WHS Regulations

Printing date 12.12.2019

Version number 6

Revision: 12.12.2019

Not classified as hazardous according to criteria of Australian Safety and Compensation Council.

1 Identification · Product identifier · Trade name: IPS e.max CAD · CAS Number: 66402-68-4 · EC number: 266-340-9 · Relevant identified uses of the substance or mixture and uses advised against No further relevant information available. • Application of the substance / the mixture Manufacture of dental prothesis Lithium disilicate glass-ceramic ingot for the CAD/CAM technology · Details of the supplier of the safety data sheet · Manufacturer/Supplier: Ivoclar Vivadent AG Bendererstrasse 2 9494 Schaan PRINCIPALITY OF LIECHTENSTEIN Tel: +423 235 35 35 / Fax: +423 235 33 60 Importer: Ivoclar Vivadent Pty. Ltd. 1-5 Overseas Drive Noble Park North VIC 3174 Tel: + 61 3 9795 9599 / Fax: + 61 3 9795 9645 • Further information obtainable from: Regulatory Affairs sds@ivoclarvivadent.com • Emergency telephone number: 131 126 (Poisons Information Centre - 24 hours / 7 days) 2 Hazard(s) Identification

• *Classification of the substance or mixture* The substance is not classified, according to the Globally Harmonised System (GHS).

- · Label elements
- GHS label elements Void
- · Hazard pictograms Void
- · Signal word Void
- · Hazard statements Void
- · Other hazards Avoid breathing sanding dust /swarf.
- · Results of PBT and vPvB assessment
- **PBT:** Not applicable.
- **vPvB:** Not applicable.

3 Composition and Information on Ingredients

- · Chemical characterisation: Substances
- · CAS No. Description

66402-68-4 ceramic materials

(Contd. on page 2)

AU

Printing date 12.12.2019

Version number 6

Revision: 12.12.2019

Trade name: IPS e.max CAD

(Contd. of page 1)

- · Identification number(s) ---
- **EC number:** 266-340-9
- Description: Lithium disilicate glass-ceramic ingot

4 First Aid Measures

- · Description of first aid measures
- · General information: No special measures required.
- · After inhalation: Grinding dust: Supply fresh air; consult doctor in case of complaints.
- *After skin contact: Generally the product does not irritate the skin.*
- After eye contact:
- Grinding dust: Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
- Mechanical effects only.
- After swallowing:
- Grinding dust: Rinse out mouth and then drink plenty of water.
- If symptoms persist consult doctor.
- Information for doctor:
- · Most important symptoms and effects, both acute and delayed No further relevant information available.
- · Indication of any immediate medical attention and special treatment needed
- No further relevant information available.

5 Fire Fighting Measures

- · Extinguishing media
- Suitable extinguishing agents:
- The product is not flammable.
- Use fire extinguishing methods suitable to surrounding conditions.
- Special hazards arising from the substance or mixture No further relevant information available.
- Advice for firefighters
- · Protective equipment: No special measures required.

6 Accidental Release Measures

- · Personal precautions, protective equipment and emergency procedures Not required.
- · Environmental precautions: No special measures required.
- Methods and material for containment and cleaning up: Pick up mechanically.
- · Reference to other sections
- See Section 7 for information on safe handling.
- See Section 8 for information on personal protection equipment.
- See Section 13 for disposal information.

7 Handling and Storage

- · Handling:
- · Precautions for safe handling
- *Only adequately trained personnel should handle this product. For use in dentistry only.*
- Information about fire and explosion protection: No special measures required.
- · Conditions for safe storage, including any incompatibilities

• Storage:

• *Requirements to be met by storerooms and receptacles:* Store only in the original receptacle.

(Contd. on page 3)

AU

Printing date 12.12.2019

Version number 6

Revision: 12.12.2019

(Contd. of page 2)

Trade name: IPS e.max CAD

- · Information about storage in one common storage facility: Not required.
- Further information about storage conditions: None.
- Specific end use(s) No further relevant information available.

8 Exposure controls and personal protection

• Additional information about design of technical facilities: No further data; see item 7.

· Control parameters

- · Ingredients with limit values that require monitoring at the workplace: Not required.
- Additional information: The lists valid during the making were used as basis.
- · Exposure controls
- · Personal protective equipment:
- General protective and hygienic measures:
- Usual hygienic measures for dental practice and dental laboratories.
- Respiratory protection: Use respiratory protective device against the effects of dust.
- Protection of hands: Not required.
- · Material of gloves -
- · Penetration time of glove material -
- Eye protection:

Always wear safety goggles during mechanical processing (grinding, sawing /cutting, drilling, milling).

Information on basic physical and chemic General Information	ical properties	
Appearance:		
Form:	Solid	
Colour:	According to product specification	
Odour:	Odourless	
Odour threshold:	Not determined.	
pH-value:	Not applicable.	
Change in condition		
Melting point/freezing point:	Undetermined.	
Initial boiling point and boiling range:	Undetermined.	
Flash point:	Not applicable.	
Flammability (solid, gas):	Product is not flammable.	
Auto-ignition temperature:	Not determined.	
Explosive properties:	Product does not present an explosion hazard.	
Explosion limits:		
Lower:	Not determined.	
Upper:	Not determined.	
Vapour pressure:	Not applicable.	
Density:	Not determined.	
Relative density	Not determined.	
Vapour density	Not applicable.	
Evaporation rate	Not applicable.	

Printing date 12.12.2019

Version number 6

Revision: 12.12.2019

Trade name: IPS e.max CAD

		(Contd. of page
· Solubility in / Miscibility with		
water:	Insoluble.	
· Partition coefficient: n-octanol/water:	Not determined.	
· Viscosity:		
Dynamic:	Not applicable.	
Kinematic:	Not applicable.	
• Other information	No further relevant information available.	

10 Stability and Reactivity

· Reactivity No further relevant information available.

- · Chemical stability Stable under normal handling and storage conditions.
- Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
- · Possibility of hazardous reactions No dangerous reactions known.

 \cdot Conditions to avoid No further relevant information available.

- · Incompatible materials: No further relevant information available.
- · Hazardous decomposition products: None under normal conditions of storage and use.

11 Toxicological Information

- · Information on toxicological effects
- · Acute toxicity
- · Skin corrosion/irritation No irritant effect.
- · Serious eye damage/irritation No irritating effect.
- · Respiratory or skin sensitisation No sensitising effects known.
- Additional toxicological information: No further relevant information available.

12 Ecological Information

- · Toxicity
- · Aquatic toxicity: No further relevant information available.
- · Persistence and degradability No further relevant information available.
- Behaviour in environmental systems:
- · *Bioaccumulative potential* No further relevant information available.
- · Mobility in soil No further relevant information available.
- · Additional ecological information:
- · General notes: Not hazardous for water.
- · Results of PBT and vPvB assessment
- *PBT:* Not applicable.
- · vPvB: Not applicable.
- · Other adverse effects No further relevant information available.

13 Disposal considerations

- · Waste treatment methods
- Recommendation

Take to an approved landfill or a waste incineration plant, under conditions approved by the local authority.

(Contd. on page 5)

Printing date 12.12.2019

Version number 6

Revision: 12.12.2019

Trade name: IPS e.max CAD

(Contd. of page 4)

· Uncleaned packaging:

• *Recommendation: Disposal must be made according to official regulations.*

UN-Number	
ADG, ADN, IMDG, IATA	Void
UN proper shipping name	
ADG, ADN, IMDG, IATA	Void
Transport hazard class(es)	
ADG, ADN, IMDG, IATA	
Class	Void
Packing group	
ADG, IMDG, IATA	Void
Environmental hazards:	
Marine pollutant:	No
Special precautions for user	Not applicable.
Transport in bulk according to Annex II	of Marpol
and the IBC Code	Not applicable.
Transport/Additional information:	Product is not classified as a dangerous good for
-	transport (ADR, IMDG, IATA).
UN "Model Regulation":	Void

15 Regulatory information

· Safety, health and environmental regulations/legislation specific for the substance or mixture

· Australian Inventory of Chemical Substances

Substance is not listed.

· Standard for the Uniform Scheduling of Medicines and Poisons

Substance is not listed.

- GHS label elements Void
- Hazard pictograms Void
- Signal word Void
- Hazard statements Void
- · Directive 2012/18/EU
- · Named dangerous substances ANNEX I Substance is not listed.
- · National regulations:
- Other regulations, limitations and prohibitive regulations
- *The product is a medical device according to the Directive 93/42/EEC.*
- Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

(Contd. on page 6)

Printing date 12.12.2019

Version number 6

Revision: 12.12.2019

Trade name: IPS e.max CAD

(Contd. of page 5)

AT

• Abbreviations and acronyms: RID: Règlement international concernant le transport des marchandises dangereuses par chemin de fer (Regulations Concerning the International Transport of Dangerous Goods by Rail) ICAO: International Civil Aviation Organisation ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

EINECS: European Inventory of Existing Commercial Chemical Substances CAS: Chemical Abstracts Service (division of the American Chemical Society)

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative